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



-80°C

8P8801.02

TECHNOPATH
Ballina, Co. Tipperary, Ireland

Maria Solana Heredia
Apoderada

Jorge Luis Marun
Farmacéutico
Co-Director Técnico

PK-SPEC-TP-048 Label Specification-Multichem S Plus assayed ICQ 08P88-10/11/12 Rev 02

Multichem S Plus Assayed ICQ (08P88-10) Kit Box Level 1

0000 **NS**

Rev: 08P88-10.02

An assayed quality control serum (human).

INFORMATION FOR USA ONLY

T E C H N O P A T H

TECHNOPATH

REF 08P88-10

Multichem S Plus (Assayed)

IVD

CE

LOT 01XXMMYY01

CONTROL **1** 12 x 5 ml

-80°C **-20°C** **YYYY-MM-DD**

www.technopathclinicaldiagnostics.com

GTIN

PRODUCT OF IRELAND

TECHNO-PATH MANUFACTURING LTD.
Fort Henry Business Park, Ballina,
Co. Tipperary, Ireland. + 353 61 525700

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

0000 **NS**

Rev: 08P88-11.02

An assayed quality control serum (human).

INFORMATION FOR USA ONLY

TECHNOPATH

REF 08P88-11

Multichem S Plus (Assayed)

IVD

CE

CONTROL **2** 12 x 5 ml


LOT 01XXMMYY02

-80°C -20°C

YYYY-MM-DD

www.technopathclinicaldiagnostics.com

GTIN

PRODUCT OF IRELAND

TECHNO-PATH MANUFACTURING LTD.
Fort Henry Business Park, Ballina,
Co. Tipperary, Ireland. + 353 61 525700

DISTRIBUTED BY

Abbott Laboratories, Inc.
Abbott Park, IL 60064 USA &
ABBOTT 65205 Wiesbaden, Germany

(01)05391523440783(240)08P88-11(17)000000(10)01XXMMYY02


Maria Solana Heredia
Apodaca
Abbott Laboratories Argentina
División Diagnósticos

Jorge Luis Marun
Farmacéutico
Co-Director Técnico
Abbott Laboratories Argentina
División Diagnósticos

PK-SPEC-TP-048 Label Specification-Multichem S Plus assayed ICQ 08P88-10/11/12 Rev 02

Multichem S Plus Assayed ICQ (8P8803) Vial label Level 3

5mL

REF

8P8803

Multichem S Plus (Assayed)

IVD

CONTROL

3

LOT


????????????

CE

YYYY-MM-DD-20°C

SN

?????-80°C

8P8803.02

TECHNOPATH

Ballina, Co. Tipperary, Ireland

Maria Solana Heredia
 Apoderada
 Abbott Laboratories Argentina
 División Diagnósticos

Jorge Luis Marun
 Farmacéutico
 Co-Director Técnico
 Abbott Laboratories Argentina
 División Diagnósticos

PK-SPEC-TP-048 Label Specification-Multichem S Plus assayed ICQ 08P88-10/11/12 Rev 02

Multichem S Plus Assayed ICQ (08P88-12) Kit Box Level 3

0000 **NS**

Rev: 08P88-12.02

An assayed quality control serum (human).

INFORMATION FOR USA ONLY

TECHNOPATH

REF 08P88-12

Multichem S Plus (Assayed)

IVD

CE


CONTROL **3** 12 x 5 ml

LOT 01XXMMYY03

-80°C -20°C **YYYY-MM-DD**

www.technopathclinicaldiagnostics.com

GTIN

PRODUCT OF IRELAND

TECHNO-PATH MANUFACTURING LTD.
Fort Henry Business Park, Ballina,
Co. Tipperary, Ireland. + 353 61 525700

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





0000 NS

Rev: 08P87-11.02

An unassayed quality control serum (human).

INFORMATION FOR USA ONLY

TECHNOPATH

TECHNOPATH

REF 08P87-11

Multichem S Plus (Unassayed)

IVD

CE

LOT 01XXMMYY02

CONTROL 2 12 x 5 ml

-80°C -20°C

YYYY-MM-DD

www.technopathclinicaldiagnostics.com

GTIN


PRODUCT OF IRELAND

TECHNO-PATH MANUFACTURING LTD.
Fort Henry Business Park, Ballina,
Co. Tipperary, Ireland. + 353 61 525700

DISTRIBUTED BY

Abbott Laboratories, Inc.
Abbott Park, IL 60064 USA &
ABBOTT 65205 Wiesbaden, Germany

(01)05391523440813(240)08P87-11(17)000000(10)01XXMMYY02

Jorge Luis Marun
Farmacéutico
Co-Director Técnico

María Solana Heredia

Page 4 of 7

Apoderada
Abbott Laboratories Argentina
División Diagnósticos

Valid only on the day of printing, 23 October 2017

Abbott Laboratories Argentina
División Diagnósticos

PK-SPEC-TP-049 Label Specification-Multichem S Plus Unassayed ICQ 08P87-

10/11/12 Rev 02

Multichem S Plus Unassayed ICQ (8P8703) Vial label Level 3

5mL

REF 8P8703

Multichem S Plus (Unassayed) IVD

CONTROL 3 

LOT ???????????? **CE**

EXP YYYY-MM-DD  -20°C


SN ?????? -80°C 


8P8703.02


TECHNOPATH

Ballina, Co. Tipperary, Ireland


 Maria Solana Heredia
 Responsable
 Abbott Laboratories Argentina
 División Diagnósticos

 Jorge Luis Marun
 Farmacéutico
 Co-Director Técnico
 Abbott Laboratories Argentina
 División Diagnósticos

TECHNOPATH

PK-SPEC-TP-049 Label Specification-Multichem S Plus Unassayed ICQ 08P87-10/11/12 Rev 02

Multichem S Plus Unassayed ICQ (08P87-12) Kit Box label Level 3

0000 **NS**

Rev: 08P87-12.02

An unassayed quality control serum (human).

INFORMATION FOR USA ONLY

TECHNOPATH

TECHNOPATH

REF 08P87-12

Multichem S Plus (Unassayed)

IVD

CE

LOT 01XXMMYY03

CONTROL **3** 12 x 5 ml

-20°C
 -80°C

YYYY-MM-DD

www.technopathclinicaldiagnostics.com

GTIN

PRODUCT OF IRELAND

TECHNO-PATH MANUFACTURING LTD.
Fort Henry Business Park, Ballina,
Co. Tipperary, Ireland. + 353 61 525700

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
1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organisation. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organisation; 2004.

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Multichem S Plus (Assayed)

INTENDED USE

Multichem S Plus (Assayed) is intended for use as an assayed quality control serum to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.

SUMMARY AND PRINCIPLE

The use of quality control material is indicated as an objective assessment of the precision of methods and techniques and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the analytical range.

COMPOSITION

This product is prepared from human serum to which purified biochemical material (extracts of human and animal origin), chemicals, drugs, preservatives and stabilizers have been added. The control is used in liquid form for convenience.

PRECAUTIONS

Warning: Biological source material. Treat as potentially infectious.

- For In Vitro Diagnostic Use.
- Each donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV, and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases.
- This product may also contain other human source material for which there is no approved test. It is recommended that this material be handled in accordance with the OSHA Standard on Bloodborne pathogens¹, Biosafety Level 2² or other appropriate biosafety practices³.
- Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.
- Material Safety Data Sheet (MSDS) is available for professional user on request.

MATERIALS INDICATED BUT NOT PROVIDED

Mixing device, e.g. roller mixer

PROCEDURE FOR USE

Multichem S Plus Control should be treated the same as patient specimens. It can be used in one of two ways with the Alinity c analyser and the Alinity c reagents:

- On board storage:** Multichem S Plus is stored on board the Alinity c analyser for the duration of the predetermined on board stability.
 - Allow the control vial to stand at room temperature (18 to 25 °C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.
 - Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogeneous solution.
 - After mixing, remove the cap and discard.
 - Place the control vial into an appropriate rack and place on the analyser within the reagent carousel. Refer to the Alinity ci-series Operations Manual for onboard use instructions.
 - The control vial remains on the analyser until the on board stability duration has expired.
 - See "Stability" section below for on board stability duration and exceptions.

- Off board storage:** Multichem S Plus can be stored in a 2 to 8 °C storage location (not on the analyser) for the duration of the predetermined off board stability.
 - Allow the control vial to stand at room temperature (18 to 25 °C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.
 - Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogeneous solution.
 - After mixing, remove the cap and transfer the appropriate volume of control material into a sample cup. Refer to the Alinity ci-series Operations Manual for instructions on loading sample cups onto the analyser.
 - Promptly replace the cap and return to the 2 to 8 °C storage location (not on the analyser) in an upright position.
 - Prior to subsequent use, remove the control vial from 2 to 8 °C storage and mix gently before sampling.
 - Continue to use the off board control vial until the off board stability duration has expired or the contents of the vial has been emptied.
 - See "Stability" section below for off board stability duration and exceptions.

- Storage**
For optimum performance, avoid storage in a frost-free freezer or in a freezer which is opened regularly for other applications such as storage of clinical samples. Store vials away from the light. Bilirubin values may decrease during the product shelf life. The Control material can only be used in the On board storage or the Off board storage application. Not a combination of both.

- Stability**
The control material is stable at -20 to -80 °C until its expiry date.

Storage Condition	Temperature	Maximum Storage Time	Exceptions	
			Analyte	Duration
On Board Storage	Alinity c System Temperature	5 days	Calcium, UIBC	3 days
Analytes not suitable for storage on board the analyzer: Carbon Dioxide, Ethanol, Carbamazepine				
Storage Condition	Temperature	Maximum Storage Time	Exceptions	
			Analyte	Duration
Off Board Storage	2 to 8 °C (Off the analyser)	10 days (after opening)	XX	X Day X Hour XX Days X Hours

- Do not use the product past the expiration date.
- If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- This product is not intended for use as a standard or calibrator.
- The Control material can only be used in the On board storage or in the Off board storage application. Not a combination of both.
- Aliquoting of this product is not recommended.

ASSIGNMENT OF VALUES

The values provided in the data sheet were derived from replicate analyses and are specific for a particular lot of product. These values have been generated using third party manufacturers' instrument systems and are specific to one measurement procedure. Technopath make no accuracy claims regarding these values. Tests were performed by the control manufacturer and/or by independent laboratories, for various methods and instrument systems. As a tool to assist in establishing their own mean, laboratories can import the values into their Alinity c system. For more details and to register for access to this file please visit www.technopathcd.com. Values are provided only as guidelines, each laboratory should

Storage Condition	Temperature	Maximum Storage Time	Exceptions	
			Analyte	Duration
On Board Storage	Alinity c System Temperature	5 days	Calcium, UIBC	3 days
Analytes not suitable for storage on board the analyzer: Carbon Dioxide, Ethanol, Carbamazepine				
Storage Condition	Temperature	Maximum Storage Time	Exceptions	
			Analyte	Duration
Off Board Storage	2 to 8 °C (Off the analyser)	10 days (after opening)	XX	X Day X Hour XX Days X Hours

LIMITATIONS

- Do not use the product past the expiration date.
- If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- This product is not intended for use as a standard or calibrator.
- The Control material can only be used in the On board storage or in the Off board storage application. Not a combination of both.
- Aliquoting of this product is not recommended.

The values provided in the data sheet were derived from replicate analyses and are specific for a particular lot of product. These values have been generated using third party manufacturers' instrument systems and are specific to one measurement procedure. Technopath make no accuracy claims regarding these values. Tests were performed by the control manufacturer and/or by independent laboratories, for various methods and instrument systems. As a tool to assist in establishing their own mean, laboratories can import the values into their Alinity c system. For more details and to register for access to this file please visit www.technopathcd.com. Values are provided only as guidelines, each laboratory should

REF

REF	Level / Niveau / Nivel / Seviye / Múc:	12 x 5mL	Level 1: Enter Lot
08P88-10 (8P8801)	Level / Niveau / Nivel / Seviye / Múc: 1	12 x 5mL	Level 1: Enter Lot
08P88-11 (8P8802)	Level / Niveau / Nivel / Seviye / Múc: 2	12 x 5mL	Level 2: Enter Lot
08P88-12 (8P8803)	Level / Niveau / Nivel / Seviye / Múc: 3	12 x 5mL	Level 3: Enter Lot

establish its own statistical limits. Laboratory means may vary from the values listed during the shelf life of the control. Technopath monitors the values over the shelf life of the control and provides update(s) at www.technopathcd.com or contact your local Abbott customer service representative.

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organisation. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organisation; 2004.

Customer Service: For questions, including suspected shipping damage, contact your local area Abbott Customer Service representative or find country specific contact information at www.abbottdiagnostics.com.

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CE

INTENDED USE

Multichem S Plus (Unassayed) is intended for use as an unassayed quality control serum to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.

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1. For In Vitro Diagnostic Use.

2. Each donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV, and nonreactive for HBsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases.

3. This product may also contain other human source material for which there is no approved test. It is recommended that this material be handled in accordance with the OSHA Standard on Bloodborne pathogens¹, Biosafety Level 2 or other appropriate biosafety practices.

4. Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

5. Material Safety Data Sheet (MSDS) is available for professional user on request.

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1. **On board storage:** Multichem S Plus is stored on board the Alinity c analyser for the duration of the predetermined on board stability.

• Allow the control vial to stand at room temperature (18 to 25 °C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.

• Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogeneous solution.

• After mixing, remove the cap and discard.

• Place the control vial into an appropriate rack and place on the analyser within the reagent carousel. Refer to the Alinity ci-series Operations Manual for onboard use instructions.

• The control vial remains on the analyser until the on board stability duration has expired.

• See "Stability" section below for on board stability duration and exceptions.

2. **Off board storage:** Multichem S Plus can be stored in a 2 to 8 °C storage location (not on the analyser) for the duration of the predetermined off board stability.

• Allow the control vial to stand at room temperature (18 to 25 °C) for 30 minutes or until completely thawed. The control must be completely solubilized before continuing to the next step.

• Place the thawed control vial on a suitable mixing device, e.g. roller mixer for 15 minutes, to achieve a homogeneous solution.

• After mixing, remove the cap and transfer the appropriate volume of control material into a sample cup. Refer to the Alinity ci-series Operations Manual for instructions on loading sample cups onto the analyser.

• Promptly replace the cap and return to the 2 to 8 °C storage location (not on the analyser) in an upright position.

• Prior to subsequent use, remove the control vial from 2 to 8 °C storage and mix gently before sampling.

• Continue to use the off board control vial until the off board stability duration has expired or the contents of the vial has been emptied.

• See "Stability" section below for off board stability duration and exceptions.

STORAGE

For optimum performance, avoid storage in a frost-free freezer or in a freezer which is opened regularly for other applications such as storage of clinical samples. Store vials away from the light. Bilirubin values may decrease during the product shelf life. The Control material can only be used in the On board storage or the Off board storage application. Not a combination of both.

STABILITY

The control material is stable at -20 to -80 °C until its expiry date.

Storage	Temperature	Maximum Storage Time	Exceptions	
			Analyte	Duration
On Board Storage	Alinity c System Temperature	5 days	Calcium, UIBC	3 days
			Analytes not suitable for storage on board the analyzer: Carbon Dioxide, Ethanol, Carbamazepine	
Storage	Temperature	Maximum Storage Time	Exceptions	
Off Board Storage	2 to 8 °C (Off the analyser)	10 days (after opening)	XX	X Day X Hour XX Days X Hours

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ASSIGNMENT OF VALUES

The values provided in the data sheet were derived from replicate analyses and are specific for a particular lot of product. These values have been generated using third party manufacturers' instrument systems and are specific to one measurement procedure.

Technopath make no accuracy claims regarding these values. Tests were performed by the control manufacturer and/or by independent laboratories, for various methods and instrument systems. As a tool to assist in establishing their own mean, laboratories can import the values into their Alinity c system. For more details and to register for access to this file please visit www.technopathcd.com.

Values are provided only as guidelines, each laboratory should establish its own statistical limits. Laboratory means may vary from the values listed during the shelf life of the control. Technopath monitors the values over the shelf life of the control and provides update(s) at www.technopathcd.com or contact your local Abbott customer service representative.

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