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Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Mindray building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China

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# Anti-HBc

## Antibody to Hepatitis B Core Antigen (CLIA)

### Order Information

Catalog No.	Package Size
Anti-HBc121	2x50 tests
Anti-HBc122	2x100 tests

### Intended Use

The CLIA series Anti-HBc assay is a Chemiluminescent Immunoassay (CLIA) for the qualitative determination of antibody to hepatitis B core antigen (Anti-HBc) in human serum or plasma.

The CLIA series Anti-HBc assay is intended to be used as an aid in the diagnosis of HBV infection and as a screening test for donated blood and blood components.

### Summary

Hepatitis B core antigen (HBcAg) is a structural protein of hepatitis B virus (HBV) and comprises 183-185 amino acids with strong immunogenicity. After HBV infection, antibodies to Hepatitis B core antigen (anti-HBc) can be detected in peripheral blood soon after the appearance of HBsAg, and often persist for life. Accordingly, Anti-HBc is an indicator for acute or past infection. Determination of Anti-HBc in association with other HBV serological markers permits the diagnosis and monitoring of HBV infection. In the absence of other HBV markers, Anti-HBc may be the only indicator of an ongoing or past HBV infection<sup>1-5</sup>.

### Assay Principles

The CLIA series Anti-HBc assay is a competitive binding assay to determine the level of antibody to hepatitis B core antigen.

In the first step, sample, sample treatment solution, neutralizing reagent are added into a reaction vessel. After incubation, anti-HBc in the sample will bind to recombinant HBcAg (rHBcAg) from neutralizing reagent.

In the second step, paramagnetic microparticles coated with specific anti-HBc monoclonal antibody (mouse IgG) and another alkaline phosphatase (ALP) labeled Anti-HBc monoclonal antibody (mouse IgG) conjugate are added the reaction vessel. After incubation, the free rHBcAg binds to Anti-HBc monoclonal antibody coated on microparticles and another anti-HBc monoclonal antibody labeled with ALP to form a complex (Microparticle Anti-HBc MAb)-rHBcAg-(Anti-HBcMAb -ALP). Afterwards, microparticles are magnetically captured while other unbound substances are removed by washing.

In the third step, the substrate solution is added to the reaction vessel. The substrate is catalyzed

by ALP in the immune complex captured on microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. A direct relationship exists between the amount of Anti-HBc in the sample and the RLUs generated during the reaction.


The presence or absence of Anti-HBc in the sample is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from calibration. If the chemiluminescent signal of the reaction is less than the cutoff signal ( $COI \leq 1.0$ ), then the sample is considered as Reactive for Anti-HBc. If the chemiluminescent signal of the reaction is higher than the cutoff signal ( $COI > 1.0$ ), then the sample is considered as Non-reactive for Anti-HBc.

### Reagents

The reagent kit is composed of four components: Ra, Rb, Rc, and Rd. The component cannot be exchanged, and the detailed information of each component is listed below:

Ra	3.8mL(50 Tests) or 6.6mL(100 Tests) paramagnetic microparticles coated with anti-HBc monoclonal antibody (mouse IgG) in 50mM TRIS buffer. Minimum concentration: 0.01% solids. Preservatives: 0.09% sodium azide.
Rb	3.5mL(50 Tests) or 6.3mL(100 Tests) alkaline Phosphatase labeled anti-HBc monoclonal antibody (mouse IgG) in MES buffer with protein (bovine serum albumin, mouse IgG) stabilizers. Minimal concentration: 0.1µg/mL. Preservative: 0.05%ProClin 300.
Rc	2.8mL(50 Tests) or 4.8mL(100 Tests) recombinant HBcAg (rHBcAg) in 50mM TRIS buffer with protein (bovine serum albumin and mouse IgG) stabilizers. Minimal concentration: 0.1µg/mL. Preservative: 0.09% sodium azide.
Rd	2.0mL(50 Tests) or 3.3mL(100 Tests) 50mM citrate buffer.

The position of each component is shown in the figure below (front view on the left and top view on the right) :

### Storage Stability

The unopened Anti-HBc reagent kit is stable up to the stated expiration date when stored at 2-8°C. The actual expiration date is stated on the label. The Anti-HBc reagent can be stored onboard at 2-8°C and used for a maximum of 28 days after opening for use.

### Reagent Preparation

Ra: Ready to use  
Rb: Ready to use  
Rc: Ready to use  
Rd: Ready to use

### Materials required but not provided

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000i.

Cat. No. Anti-HBc221: Anti-HBc Calibrators, 1x2.0mL for each level of calibrator C0 and C1;  
Cat. No. Anti-HBcN321 or Anti-HBcN322: Anti-HBc Negative Control 3x2.0mL or 6x2.0mL;  
Cat. No. Anti-HBcP321 or Anti-HBcP322: Anti-HBc Positive Control, 3x2.0mL or 6x2.0mL;  
Cat. No. CS511: Substrate Solution, 4x115mL;  
Cat. No. WB411: Wash Buffer; 1x10L;  
Cat.No.115-025456-00: Reaction Vessel.

### Instrument System

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 and CLIA 2000.

### Specimen Collection and Preparation

Human serum, heparin plasma or citrate plasma is suitable for the test.

Human serum is recommended for this assay. Centrifuge the specimens at 3500 rpm for a minimum of 10 minutes after clot formation is complete. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment. Test the specimens within two hours after centrifugation. If testing is not completed within 24 hours, transfer the supernatant into tubes for storage. Specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 7 days, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles, which may cause sample deterioration. Specimen can be used after a maximum of five cycles of freeze and thaw. Previously frozen specimens must be mixed gently and thoroughly after thawing and centrifuged at 10000g for 10 minutes.

Do not use specimen with the following conditions:

- heat inactivated
- grossly hemolyzed
- obvious microbial contamination
- visible fibrin or other cell debris

### Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully to get sufficient information, such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well. Before loading the Anti-HBc reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to

resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Do not invert opened reagent bottle.

Unsealed reagent kit can be stored onboard (in the reagent carousel) or at 2~8 °C after unloading from the reagent carousel.

Before re-loading a reagent, manually rotate the gear below the magnetic bead bottle by hand for mixing until the magnetic beads deposited at bottom have been well mixed.

Refer to the "On-line load of Reagents", "Off-line load of Reagents" and "Unloading Reagents" sections of system operation manual for details of the procedures.

Reagent pack opened for more than 28 days is not recommended for use.

This assay requires 50 µL of sample volume for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

### Calibration

The calibration information is stored in the barcode attached in the reagent and calibrator pack. When performing the calibration, scan the information from the barcodes into the system first, and then test the calibrators at two levels. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified ranges. For detailed instruction of calibration, refer to the system operation manual.

### Quality Control

For verifying the reliability of assay system, quality controls should be run at least once every 24 hours, after a new reagent pack is loaded, or after every calibration. The recommended quality controls for this assay are Wiener lab. Anti-HBc Negative Control and Anti-HBc Positive Control.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the system operation manual to check up the system. If the quality control results are still out of the specified ranges, please contact Wiener lab. Customer Service for help.

### Calculation

The Wiener lab. CLIA series Chemiluminescence analyzer calculates the mean RLU signals from

replicate tests of both calibrator C1 and C0. The Wiener lab. CLIA series Chemiluminescence analyzer calculates Anti-HBc result based on the ratio of the sample RLU to the cutoff RLU (Cutoff Index, COI) for each specimen and control.

Cutoff RLU = Mean RLU of C0 - [(Mean RLU of C0 - Mean RLU of C1) x Calibration Coefficient]

Calibration Coefficient is specific for each reagent and calibrator lot.

COI = Sample RLU/Cutoff RLU

### Interpretation of Results

1. Samples with Anti-HBc COI > 1.00 are Non-reactive. These samples are considered as negative for Anti-HBc, and do not need further testing.
2. Samples with Anti-HBc COI ≤ 1.00 are Initially Reactive.

Initially Reactive samples need to be retested. Samples should be transferred to a centrifuge tube and centrifuged at ≥ 10,000 RCF (Relative Centrifugal Force) for 10 minutes, and then retested in duplicate with Wiener lab. Anti-HBc kit.

Initially reactive samples with a COI of > 1.00 in both retests are considered as Non-reactive, and is negative for Anti-HBc. Initially reactive samples with a COI of ≤ 1.00 in either of the retests are considered as Repeatedly Reactive.

### Limitation of Measurement

The Anti-HBc result of any given specimen can vary with assays from different manufacturers, due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, and clinical history, etc. to make clinical decisions.

### Performance Characteristics

#### Sensitivity

The **Antibody to Hepatitis B Core Antigen (CLIA)** assay sensitivity at the cutoff is ≤ 0.6 PEI U/mL.

A total of 522 specimens including 119 serum and 403 plasma samples from individuals infected HBV were tested. Of the 522 specimens, 521 (99.81%) were reactive

Category	Number of Tested	Number of Reactive (% of total)
Serum	119	118
Plasma	403	403
Total	522	521 (99.81%)

#### Precision

The CLIA series Anti-HBc assay is designed to have a precision of ≤ 10% (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2.<sup>6</sup> Two levels of quality controls were tested in duplicate in two separate

runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean Anti-HBc (COI)	Within-Run CV	Between-Run CV	Within-Device CV
Negative Control	2.79	2.92%	3.62%	6.45%
Positive Control	0.42	2.68%	3.85%	5.96%

#### Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 20 mg/dL, triglycerides up to 1000 mg/dL, total protein up to 10 g/dL, rheumatoid factor (RF) up to 1500 IU/mL, anti-DNA antibody (ANA) up to 2000 U/L, and mouse anti-human antibody (HAMA) will not interfere with the Anti-HBc assay. These substances show less than 10% interference at indicated concentrations. 5096 specimens from blood donors, which are negative on a commercially available reference Anti-HBc assay, were tested. Of the 5096 specimens, 5077 (99.63%) were nonreactive. 201 specimens from hospitalized patients and 100 with potentially interfering substances which are negative on a commercially available reference Anti-HBc assay, were tested. All 300 specimens were identified as nonreactive.

Category	Number of Tested	Number of Nonreactive (% of total)
Negative Blood Donors	5096	5077 (99.63%)
Negative Hospitalized Patients	201	200 (100%)
Potentially Interfering Substances <sup>a</sup>	100	100 (100%)

<sup>a</sup>: Category included the following: Anti-HCV(20), pregnant women (11), Pregnant women multipara (9), Anti-HTLV 1/2(10), Rheumatoid factors(15) Anti-HIV(20) and Anti-HEV(15).

#### Warnings and Precautions

1. For in vitro diagnostic use only.
2. The operators must be trained professionals and the operation procedure should follow the system operations manual.
3. Follow all the rules in handling laboratory reagents and take necessary safety precautions.

 **Warning:** May cause an allergic skin reaction.


Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

4. Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Wiener

lab. system, or using Wiener lab. reagent kits on other systems.

5. Do not use reagent kits beyond the expiration date.
6. Do not use reagents mixed from different reagent lots.
7. Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
8. Reagent pack opened for more than 28 days is not recommended for use.
9. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
10. All the specimen and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
11. The Material Safety Data Sheet (MSDS) is available upon request.

### Graphical Symbols

### Bibliography

1. Koziol DE, Holland PV, Alling DW, et al. Antibody to Hepatitis B Core Antigen as a Paradoxical Marker for Non-A, Non-B Hepatitis Agents in Donated Blood. *Ann Intern Med* 1986; 104:488-495.
2. Boscato, LM and Sturat, MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988; 34:27.
3. Szmuness W, Hoofnagle JH, Stevens CE, et al. Antibody against the hepatitis type B core antigen: A new tool for epidemiologic studies. *Am J Epidemiol* 1976; 104: 256-262.
4. Hoofnagle JH, Seeff LB, Buskell-Bales Z, et al. Type B hepatitis after transfusion with blood containing antibody to hepatitis B core antigen. *N Engl J Med* 1978; 198: 1379-1383.
5. Kumar S, Pound DC, Serologic diagnosis of viral hepatitis. Postgraduate *Medicine* 1992; 92:55-65
6. CLSI. EP5-A2: Vol. 24, No. 25, *Evaluation of Precision Performance of Quantitative Measurement Method*; Approved Guideline - Second Edition

© 2012 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
All rights Reserved  
**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

**E-mail Address:** service@mindray.com.cn

**Tel:** +86-755-26582888

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH (Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

**Tel:** 0049-40-2513175

**Fax:** 0049-40-255726

Distributed by:  
Wiener Laboratorios S.A.I.C.  
Riobamba 2944  
2000 Rosario - Argentina  
<http://www.wiener-lab.com>  
A.N.M.A.T. Registered product  
PM-1102-170  
Tech. Dir.: Viviana E. Cétola  
Biochemist

# Anti-HBc

## Antibody to Hepatitis B Core Antigen (CLIA)

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
**Símbolos gráficos**

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1. Kozioł DE, Holland PV, Alling DW, et al. Antibody to Hepatitis B Core Antigen as a Paradoxical Marker for Non-A, Non-B Hepatitis Agents in Donated Blood. *Ann Internal Med* 1986; 104:488-495.
2. Boscato, LM y Sturat, MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27.
3. Szmuness W, Hoofnagle JH, Stevens CE, et al. Antibody against the hepatitis type B core antigen: A new tool for epidemiologic studies. *Am J Epidemiol* 1976; 104: 256-262.
4. Hoofnagle JH, Seeff LB, Buskell-Bales Z, et al. Type B hepatitis after transfusion with blood

WIENER Laboratorios S.A.I.C.

*Mey*  
Dra. VIVIANA E. CETOLA  
DIRECTORA TECNICA

Anti-HBc

Anti-HBc

2 x 50 tests  
P/N: XXX-XXXXXX-XX

### Antibody to Hepatitis B Core Antigen (CLIA)

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

### Intended Use

Anti-HBc Positive Control is used for monitoring the accuracy and precision in the qualitative measurement of Anti-HBc by Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000. It is also intended for clinical laboratories to assess their capability in quality control of lab work.

### Summary

Anti-HBc Positive Control has similar components as the test samples. Therefore, the control can be tested by following the same procedure as the specimen, measured values of the control is compared to the target values and range to see if the control meets its specifications.

### Components

Anti-HBc Positive Control contains Anti-HBc Positive human Plasma in the range of 0.10~0.90 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative.

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- The control values are lot-specific with each reagent pack as listed in the value sheet.
- The controls should be run daily in

parallel with the patient samples, after each calibration, or after the reagent lot has been changed.

- Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.
- Raw material for control production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.
- As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard<sup>1</sup>.
- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material in compliance with local guidelines.
- The target value and range of Wiener lab. Anti-HBc Positive Control are for reference only. Each laboratory should establish its own internal quality control procedure and set up its own target value and range.
- Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

### Preparation

- Take the vial out of 2~8 °C refrigerator and allow standing on bench for no more than 30 min to reach room temperature.
- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Dispense the required volume into a sample cup, and then run QC test and analyze it in the same way as patient samples. Control can be set up manually or imported and the system automatically analyzes information from the bar code.

### Storage and Stability

Anti-HBc Positive Control is stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 7 days at 2~8°C if capped tightly to avoid microbial contamination.

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Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBc121 or Anti-HBc122 : Antibody to Hepatitis B core Antigen (CLIA), 2x50 tests or 2x100 tests;

Cat. No. Anti-HBc221: Anti-HBc Calibrators, 1x2.0mL for each level of calibrator CO<sub>1</sub>, C1;

Cat. No. Anti-HBcN321 or Anti-HBcN322: Anti-HBc Negative Control, 3x2.0mL or 6x2.0mL;

Cat. No. CS511: Substrate Solution, 4x115mL;

Cat. No. WB411: Wash Buffer, 1x10L;

Cat.No.115-025456-00: Reaction Vessel.


### Assay Values

The control values (the target value and range) determined by Wiener lab. standard measurement procedure are listed in the target value sheet. The target value was obtained from Wiener lab. measurement system and the range was calculated as the target value  $\pm$  3 standard deviations. The control value is lot-specific, Please check the lot number of each control before use.

### Quality Control

Use Wiener lab. Anti-HBc Positive Control to monitor the status of each calibration. The value of the calibrator is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. The control result should be within the defined range as shown in the control target value sheet. If the control is out of the range, the measurement system and reagents should be checked. For instance, expiration date or storage condition of the calibrators, reagent, and controls; possible contamination of reagent; position of reagent; sample placement on the analyzer; the parameter setting in software; and operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and measurement procedure.

### Graphical Symbols

### References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard.

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**Tel:** +86-755-81888998

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH(Europe)

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
## Símbolos gráficos


In vitro diagnostic medical device Batch code Catalogue Number

European Conformity Authorized representative in the European Community

Control Consult Instructions for Use Biological Risks

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


## Intended Use

Anti-HBc Negative Control is used for monitoring the accuracy and precision in the qualitative measurement of Anti-HBc by Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000. It is also intended for clinical laboratories to assess their capability in quality control of lab work.

## Summary

Anti-HBc Negative Control has similar components as the test samples. Therefore, the control can be tested by following the same procedure as the specimen; measured value of the control is compared to the target values and range to see if the control meets its specifications.

## Components

Anti-HBc Negative Control is prepared from Anti-HBc negative human Plasma with 0.09% sodium azide and 0.05% ProClin 300 as the preservative.

## Warnings and Precautions

- For in vitro diagnostic use only.
- The operators must be professionally trained and qualified with certain certificates, and the operation procedure should follow the system operations manual.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.

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- The control value is lot-specific with each reagent pack as listed in the value sheet.
- The controls should be run daily in parallel with the patient samples,

after each calibration, or after the reagent lot has been changed.

- Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.
- Raw material for control production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.
- As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard<sup>1</sup>.
- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material in compliance with local guidelines.
- The target value and range of Wiener lab. Anti-HBc Negative Control are for reference only. Each laboratory should establish its own internal quality control procedure and set up its own target value and range.
- Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

## Preparation

- Take the vial out of 2~8 °C refrigerator and standing on bench for no more than 30 min to reach room temperature.
- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Dispense the required volume into a sample cup, and then run QC test and analyze it in the same way as patient samples. Control can be set up manually or imported and the system automatically analyzes information from the bar code.

## Storage and Stability

Anti-HBc Negative Control is stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 7 days at 2~8°C if capped tightly to avoid microbial contamination.

## Materials required but not provided

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBc121 or Anti-HBc122: Antibody to Hepatitis B core Antigen (CLIA), 2x50 tests or 2x100 tests;

Cat. No. Anti-HBc221: Anti-HBc Calibrators, 1x2.0mL for each level of calibrator C0 , C1;

Cat. No. Anti-HBcP321 or Anti-HBcP322: Anti-HBc Positive Control, 3x2.0mL or 6x2.0mL;

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
## Quality Control

Use Wiener lab. Anti-HBc Negative Control to monitor the status of each calibration. The value of the calibrator is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. The control result should be within the defined range as shown in the control target value sheet. If the control is out of the range, the measurement system and reagents should be checked. For instance, expiration date or storage condition of the calibrators, reagent, and controls; possible contamination of reagent; position of reagent; sample placement on the analyzer; the parameter setting in software; and operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and measurement procedure.

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**Address:** Eiffeustraße 80, Hamburg 20537, Germany

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## Anti-HBc Calibrators

### Order Information

Catalog No.	Package Size
Anti-HBc221	C0:1×2.0 mL/vial C1:1×2.0 mL/vial

### Intended Use

The Wiener lab. Anti-HBc Calibrators (Anti-HBc CAL) are intended to calibrate the Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000 for the qualitative determination of human antibodies to hepatitis B core antigen (Anti-HBc) in human serum and plasma.

### Summary

Anti-HBc Calibrators are prepared by adding human anti-HBc into HBV negative human plasma. The Wiener lab. CLIA measurement system calculates an Anti-HBc result based on the ratio of the sample RLU to the cutoff RLU (S/CO) for each specimen and control.

Cutoff RLU=Mean RLU of C0 - (Mean RLU of C0 - Mean RLU of C1) × calibration factor.

The calibration factor is stored in calibration card for each lot. Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence analyzers, Wiener lab. reagent kits, calibrators, and controls.

### Components

C0	Anti-HBc negative human plasma in the range of 1.50~10.00 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative
C1	Anti-HBc Positive human plasma in the range of 0.00~0.50 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative
Calibration Card	1

### Warnings and Precautions

- For in vitro diagnostic use only.
- The operators must be professionals trained and qualified with certain certificates, and the operation procedure should follow the system operations manual.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.

**Warning:** May cause an allergic skin reaction.

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- The calibration factor is lot-specific with each reagent pack as listed in the Calibration Card.
- Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure. Recalibrate the system when controls are out of specified ranges.
- Do not mix calibrators from different lots.
- Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.
- Raw material for calibrator production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.
- As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazards<sup>1</sup>.
- If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material in compliance with local guidelines.
- Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

### Testing Method

- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Refer to calibration card for accurate information about lot-specific Anti-HBc calibration factor.

### Preparation

- Take the vial out of 2-8°C refrigerator and allow standing on bench for no more than 30 min to reach room temperature.
- Dispense required volume of each calibrator into a sample cup, follow the calibration procedure described in the instrument manual. Use the Calibrators' lot-specific calibration factor to set calibration parameters.
- Once the calibrator information is input, you are allowed to assign calibrator position only when the system status is Standby or Stopped.

### Storage and Stability

The calibrators are stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 30 days at 2~8°C when capped tightly to avoid microbial contamination.

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Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBc121 or Anti-HBc122: Antibody to Hepatitis B core Antigen (CLIA), 2×50 tests or 2×100 tests ;

Cat. No. Anti-HBcN321 or Anti-HBcN322: Anti-HBc Negative Control 3×2.0mL or 6×2.0mL;

Cat. No. Anti-HBcP321 or Anti-HBcP322: Anti-HBc Positive Control, 3×2.0mL or 6×2.0mL;

Cat. No. CS511: Substrate Solution, 4×115mL;

Cat. No. WB411: Wash Buffer, 1×10L;

Cat.No.115-025456-00: Reaction Vessel.

### Traceability of Calibrator


The traceability process is based on EN ISO 17511:2003<sup>3</sup>, the analyte in the calibrators is traceable to the Anti-HBc Reference Material of the Paul Ehrlich Institute, Langen (Germany). The calibration factor of calibrators is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. If other system or method is used, values of the calibrators must be re-assigned, as bias may exist among different systems and methods. Details about traceability are available upon request.

### Quality Control

Use Wiener lab. controls as the samples to monitor the status of each calibration. The control results should be within the defined

ranges as shown in the control value sheet. If the controls are out of specified ranges, the measurement system should be checked. For instance, expiration date or storage condition of the calibrators, reagent and control; contamination of reagent; position of reagent or sample placed on the analyzer; the parameter setting in software; or performance of the analyzer.

### Graphical Symbols

### References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard.
- ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.

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**E-mail** **Address:** service@mindray.com.cn

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2016 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved