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Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Mindray building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China

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Digitally signed by Gestion Documental  
Electronica  
Date: 2023.04.17 20:01:25 -03:00

# Anti-HBe

## Antibody to Hepatitis B e

### Antigen (CLIA)

#### Order Information

Catalog No.	Package Size
Anti-HBe121	2×50 tests
Anti-HBe122	2×100 tests

#### Intended Use

The CLIA series Anti-HBe assay is a Chemiluminescent Immunoassay (CLIA) for the qualitative determination of antibody to hepatitis B virus e antigen (Anti-HBe) in human serum or plasma.

The CLIA series Anti-HBe assay is intended to be used as an aid in the diagnosis and monitoring of hepatitis B viral infection.

#### Summary

Hepatitis B e antigen (HBeAg) is a cleavage product of pre-core protein encoded by pre-C/C gene of hepatitis B virus (HBV)<sup>1</sup>. HBeAg and its antibody (anti-HBe) are commonly found in HBV infection. After HBV infection, HBeAg rises sharply following the appearance of HBsAg, which indicates the active reproduction of HBV. Serconversion from HBeAg to anti-HBe is usually indicative of recovery from acute infection with a reduced level of infectivity. A negative HBeAg result may occur in early acute infection, before the peak of viral replication or during early convalescence when HBeAg has declined below detectable levels. Measurement of anti-HBe can distinguish between these two phases. A subset of hepatitis B patients is serologically HBeAg negative, but anti-HBe positive. These patients may also be positive for serum hepatitis B virus DNA<sup>2</sup>.

#### Assay Principles

The CLIA series Anti-HBe assay is a competitive binding assay to determine the level of antibody to hepatitis B e antigen.

In the first step, sample, sample treatment solution and neutralizing reagent are added into a reaction vessel. After incubation, anti-HBe in the sample will bind to recombinant HBeAg (rHBeAg) from neutralizing reagent.

In the second step, paramagnetic microparticles coated with specific anti-HBe monoclonal antibody (mouse) and another alkaline phosphatase (ALP) labeled anti-HBe monoclonal antibody (mouse) conjugate are added to the reaction vessel. After incubation, the free rHBeAg binds to anti-HBe monoclonal antibody coated on microparticles and another anti-HBe monoclonal antibody labeled with ALP to form a complex: (Microparticle- Anti-HBe mAb)-rHBeAg-(Anti-HBe mAb-

ALP). Microparticles are magnetically captured while other unbound substances are removed by washing. In the third step, the substrate solution is added to the reaction vessel. The substrate is catalyzed by ALP in the immune complex captured on the microparticles. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. A direct relationship exists between the amount of Anti-HBe in the sample and the RLUs generated during the reaction.


The presence or absence of Anti-HBe in the sample is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from system calibration. If the chemiluminescent signal of the reaction is less than the cutoff signal (COI≤1.0), then the sample is considered as reactive for Anti-HBe. If the chemiluminescent signal of the reaction is higher than the cutoff signal (COI>1.0), then the sample is considered as Non-reactive for Anti-HBe.

#### Reagents Components

The reagent kit is composed of three components: Ra, Rb, and Rc. The component cannot be exchanged, and the detailed information of each component is listed below:

Ra	3.8 mL (50 Tests) or 6.6 mL (100 Tests) anti-HBe monoclonal antibody (mouse IgG) coated paramagnetic microparticles in 50 mM TRIS buffer. Minimum concentration: 0.01% solids. Preservatives: 0.05% ProClin 300 and 0.09% NaN <sub>3</sub> .
Rb	3.5 mL (50 Tests) or 6.3 mL (100 Tests) Alkaline Phosphatase labeled anti-HBe monoclonal antibody (mouse IgG) in MES buffer with protein (bovine serum albumin, mouse IgG) stabilizers. Minimal concentration: 0.1 µg/mL. Preservative: 0.05% ProClin 300.
Rc	3.5 mL (50 Tests) or 6.3 mL (100 Tests) Recombinant HBeAg (rHBeAg) in 50 mM Tris buffer with protein (bovine serum albumin) stabilizers. Minimal concentration: 0.1µg/mL. Preservative: 0.05%ProClin 300 and 0.09%NaN <sub>3</sub> .
Rd	3.5mL (50 Tests) or 6.3mL (100 Tests) Sample treatment solution (50mM citrate buffer). Preservative: 0.05%ProClin 300.

The position of each component is shown in the figure below (front view on the left and top view on the right):

#### Storage Stability

The unopened Anti-HBe (CLIA) reagent is stable up

to the stated expiration date when stored at 2-8°C. The actual expiration date is stated on the label.

The Anti-HBe (CLIA) reagent can be stored onboard at 2-8°C and used for a maximum of 28 days after opening for use.

#### Reagent Preparation

Ra: Ready to use

Rb: Ready to use

Rc: Ready to use

Rd: Ready to use

#### Materials required but not provided

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBe221: Anti-HBe Calibrators, 1×2.0mL for each level of calibrator C0 and C1:

Cat. No. Anti-HBeN321 or Anti-HBeN322: Anti-HBe Negative Control 3×2.0mL or 6×2.0mL;

Cat. No. Anti-HBeP321 or Anti-HBeP322: Anti-HBe Positive Control, 3×2.0mL or 6×2.0mL;

Cat. No. CS511: Substrate Solution, 4×115mL;

Cat. No. WB411: Wash Buffer;

Cat.No.115-025456-00: Reaction Vessel.

#### Instrument System

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 and CLIA 2000.

#### Specimen Collection and Preparation

Human serum or heparin plasma is suitable for the test.

Human serum is recommended for this assay. Centrifuge the specimens at 3500 rpm for a minimum of 10 minutes after clot formation is complete. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment. Test the specimens within two hours after centrifugation. If testing is not completed within 24 hours, transfer the supernatant into tubes for storage. Specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 7 days, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles, which may cause sample deterioration. Specimen can be used after a maximum of five cycles of freeze and thaw. Previously frozen specimens must be mixed gently and thoroughly after thawing and centrifuged at 10000g for 10 minutes.

Do not use specimen with the following conditions:

- heat inactivated
- grossly hemolyzed
- obvious microbial contamination
- visible fibrin or other cell debris

#### Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully to get sufficient information, such as operation instructions, sample preservation and management, safety precautions, and maintenance. Prepare all required materials for the assay as well.

Before loading the Anti-HBe reagent kit on the instrument for the first time, invert unopened reagent bottle gently for at least 30 times to resuspend the microparticles, which have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be homogenized, it is recommended not to use this bottle of reagent. Contact Wiener lab. Customer Service for help. Do not invert opened reagent bottle.

Unsealed reagent kit can be stored onboard (in the reagent carousel) or at 2~8°C after unloading from the reagent carousel.

Before re-loading a reagent, manually rotate the gear below the magnetic bead bottle by hand for mixing until the magnetic beads deposited at bottom have been well mixed.

Refer to the "On-line load of Reagents", "Off-line load of Reagents" and "Unloading Reagents" sections of system operation manual for details of the procedures.

Reagent pack opened for more than 28 days is not recommended for use.

This assay requires 50 µL of sample volume for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

#### Calibration

The calibration information is stored in the barcode attached in the reagent and calibrator pack. When performing the calibration, scan the information from the barcodes into the system first, and then test the calibrators at two levels. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

#### Quality Control

For verifying the reliability of assay system, quality controls should be run at least once every 24 hours, after a new reagent pack is loaded, or after every calibration. The recommended quality controls for this assay are Wiener lab. Anti-HBe Negative Control and Anti-HBe Positive Control.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the system operation manual to check up the system. If the quality control results are still out of the specified ranges, please contact

## Quality Control

For verifying the reliability of assay system, quality controls should be run at least once every 24 hours, after a new reagent pack is loaded, or after every calibration. The recommended quality controls for this assay are Wiener lab. Anti-HBe Negative Control and Anti-HBe Positive Control.

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Refer to the system operation manual to check up the system. If the quality control results are still out of the specified ranges, please contact Wiener lab. Customer Service for help.

## Calculation

The Wiener lab. CLIA series Chemiluminescence analyzers calculates the mean RLU signal from the duplicate tests of Calibrator C0 and C1, Anti-HBe result is based on the ratio of the mean RLU to the cutoff RLU of each specimen (Cutoff Index, COI)

$$\text{Cutoff RLU} = \text{Mean RLU of C0} - [(\text{Mean RLU of C0} - \text{Mean RLU of C1}) \times \text{Calibration Coefficient}]$$

Calibration Coefficient is specific for each reagent and calibrator lot.

$$\text{COI} = \text{Sample RLU} / \text{Cutoff RLU}$$

## Interpretation of Results

Samples with anti-HBe COI > 1.00 are Non-reactive, and is negative for anti-HBe assay. No further testing is needed for these samples.

Samples with anti-HBe COI ≤ 1.00 are Initially Reactive.

Initially Reactive samples need to be retested. Samples should be transferred to centrifuge tubes and centrifuged at ≥ 10,000 RCF (Relative Centrifugal Force) for 10 minutes, and retested in duplicate.

Initially reactive samples with COI of > 1.00 in both retests are considered as Non-reactive, and is negative for Anti-HBe. No further testing is needed.

Initially reactive samples with a COI of ≤ 1.00 in either of the retests are considered as Repeatedly Reactive.

## Limitation of Measurement

The Anti-HBe result of a given specimen can vary, in assays from different manufacturers, due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, and clinical history, etc. to make clinical decisions.

## Performance Characteristics

### Sensitivity

The CLIA series Anti-HBe assay sensitivity at the cutoff is ≤ 0.45 PEI U/mL.

### Diagnostic Sensitivity

A total of 200 specimens from individuals of different stages of HBV infection were tested. Of the 200 specimens, 200 (100%) were reactive

### Precision

The CLIA series Anti-HBe assay is designed to have a precision of ≤ 10% (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2<sup>3</sup>. Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean Anti-HBe (COI)	Within-Run CV	Between-Run CV	Within-Device CV
Negative Control	1.30	3.34%	2.58%	3.40%
Positive Control	0.33	3.45%	2.74%	3.51%

### Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 10 mg/dL, triglycerides up to 1500 mg/dL, total protein up to 10 g/dL, rheumatoid factor (RF) up to 375 IU/mL, anti-DNA antibody (ANA) up to 1000 U/L, and mouse anti-human antibody (HAMA) will not interfere with the Anti-HBe assay. These substances show less than 10% interference at indicated concentrations.

200 specimens from blood donors, which were negative on a commercially available reference Anti-HBe assay, were tested. All 200 specimens were nonreactive.

200 specimens from hospitalized patients and 50 with potentially interfering substances which were negative on a commercially available reference Anti-HBe assay, were tested. All 250 specimens were identified as nonreactive.

Category	Number of Tested	Number of Nonreactive (% of total)
Negative Blood Donors	200	200 (100%)
Negative Hospitalized Patients	200	200 (100%)
Potentially Interfering Substances	50	50 (100%)

a: Category included the following: Anti-HCV(10), pregnant women (5), Pregnant women multipara (5), Anti-HTLV 1/2(5), Rheumatoid factors(10), Anti-HIV(10) and Anti-HEV(5).

## Warnings and Precautions

- For in vitro diagnostic use only.
- The operators must be trained professionals and the operation procedure should follow the system operations manual.


- Follow all the rules in handling laboratory reagents and take necessary safety precautions.

 **Warning:** May cause an allergic skin reaction.

Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/ attention. Take off contaminated clothing and wash it before reuse.

- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Wiener lab. system, or using Wiener lab. reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. Specimens and reaction wastes should be handled in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

## Graphical Symbols

## Bibliography

- Bruss V, et al. Formation of transmembranous hepatitis B e-antigen by cotranslation in vitro processing of the viral precore protein. *Virology* 1988;163:268-275
- Brunetto MR. et al. A new hepatitis B virus strain in patients with severe anti-HBe positive chronic hepatitis B. *J Hepatol* 1990;10:258-261
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline –Second Edition.

© 2012 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

**E-mail Address:** service@mindray.com.cn

**Tel:** +86-755-26582888

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH (Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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
Anti-HBe

418101302040100101


Antibody to Hepatitis B e Antigen (CLIA) Anti-HBe

50 test

Ra: 3.8 mL; Rb: 3.5 mL;  
Rc: 3.5 mL; Rd: 3.5 mL

Anti-HBe

418101302040100101


Antibody to Hepatitis B e Antigen (CLIA) Anti-HBe

100 test

Ra: 6.6 mL; Rb: 6.3 mL;  
Rc: 6.3 mL; Rd: 6.3 mL

# Anti-HBe Calibrators

## Order Information

Catalog No.	Package Size
Anti-HBe221	C0:1×2.0 mL/vial C1:1×2.0 mL/vial

## Intended Use

The Wiener lab. Anti-HBe Calibrators (Anti-HBe CAL) are intended to calibrate the Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000 for the qualitative determination of human antibodies to hepatitis B e antigen (Anti-HBe) in human serum and plasma.

## Summary

Anti-HBe Calibrators are prepared by adding human antibodies to hepatitis B e antigen into HBV negative human plasma. The Wiener lab. CLIA measurement System calculates an Anti-HBe result based on the ratio of the sample RLU to the cutoff RLU (S/CO) for each specimen and control.

Cutoff RLU=Mean RLU of C0 - (Mean RLU of C0 - Mean RLU of C1) × calibration factor.

The calibration factor is stored in calibration card for each lot. Wiener lab. CLIA measurement system is composed of Wiener lab. CLIA series Chemiluminescence analyzers, Wiener lab. reagent kits, calibrators, and controls.

## Components

C0	Anti-HBe negative human plasma in the range of 1.10~4.00 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative
C1	Anti-HBe Positive human plasma in the range of 0.00~0.50 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative
Calibration Card	1

## Warnings and Precautions

- For in vitro diagnostic use only.
- The operators must be professionals trained and qualified with certain certificates, and the operation procedure should follow the system operations manual.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.

**Warning:** May cause an allergic skin reaction.

Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

4. The calibration factor is lot-specific with each reagent pack as listed in the Calibration Card.

5. Calibration is recommended after each reagent lot change, special maintenance, or troubleshooting procedure. Recalibrate the system when controls are out of specified range.

6. Do not mix calibrators from different lots.

7. Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.

8. Raw material for calibrator production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.

9. As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazards<sup>1</sup>.

10. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.

11. Material safety data sheet is available upon request.

12. Disposal of all waste material in compliance with local guidelines.

13. Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

## Testing Method

- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Refer to calibration card for accurate information about lot-specific Anti-HBe calibration factor.

## Preparation

- Take the vial out of 2-8°C refrigerator and allow standing on bench for no more than 30 min to reach room temperature.
- Dispense required volume of each calibrator into a sample cup, follow the calibration procedure described in the instrument manual. Use the Calibrators' lot-specific calibration factor to set calibration parameters.
- Once the calibrator information is input, you are allowed to assign calibrator position only when the system status is Standby or Stopped.

## Storage and Stability

The calibrators are stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 30 days at 2~8°C when capped tightly to avoid microbial contamination.

## Materials required but not provided

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBe121 or Anti-HBe122 : Antibody to Hepatitis B e Antigen (CLIA), 2×50 tests or 2×100 tests:

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Cat.No.115-025456-00: Reaction Vessel.

## Traceability of Calibrator


The traceability process is based on EN ISO 17511:2003<sup>3</sup>, the analyte in the calibrators is traceable to the Anti-HBe Reference Material of the Paul Ehrlich Institute, Langen (Germany). The calibration factor of calibrators is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. If other system or method is used, values of the calibrators must be re-assigned, as bias may exist among different systems and methods. Details about traceability are available upon request.

## Quality Control

Use Wiener lab. controls as the samples to monitor the status of each calibration. The control results should be within the defined ranges as shown in the control value sheet.

If the controls are out of specified ranges, the measurement system should be checked. For instance, expiration date or storage condition of the calibrators, reagent and control; contamination of reagent; position of reagent or sample placed on the analyzer; the parameter setting in software; or performance of the analyzer.

## Graphical Symbols

## References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard.
- ISO 17511:2003. In vitro diagnostic medical devices—Measurement of quantities in biological samples—Metrological traceability of values assigned to calibrators and control materials.

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**Tel:** +86-755-81888998

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH(Europe)

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


## Anti-HBe Positive Control

### Order Information

Catalog No.	Package Size
Anti-HBeP321	3×2.0 mL/vial
Anti-HBeP322	6×2.0 mL/vial

### Intended Use

Anti-HBe Positive Control is used for monitoring the accuracy and precision in the qualitative measurement of Anti-HBe by Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000. It is also intended for clinical laboratories to assess their capability in quality control of lab work.

### Summary

Anti-HBe Positive Control has similar components as the test samples. Therefore, the control can be tested by following the same procedure as the specimen; measured value of the control is compared to the target value and range to see if the control meets its specification.

### Components

Anti-HBe Positive Control contains Anti-HBe Positive human Plasma in the range of 0.10~0.90 COI, with 0.09% sodium azide and 0.05% ProClin 300 as the preservative.

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- The control value is lot-specific with each reagent pack as listed in the value sheet.

- The control should be run daily in parallel with the patient samples, after each calibration, or after the reagent lot has been changed.
- Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.
- Raw material for control production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.
- As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard<sup>1</sup>.
- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material in compliance with local guidelines.
- The target value and range of Wiener lab. Anti-HBe Positive Control are for reference only. Each laboratory can establish its own internal quality control procedure and set up its own target value and range.
- Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

### Preparation

- Take the vial out of 2~8 °C refrigerator and allow standing on bench for no more than 30 min to reach room temperature.
- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Dispense the required volume into a sample cup, and then run QC test and analyze it in the same way as patient samples. Control can be set up manually or imported and the system automatically analyzes information from the bar code.

### Storage and Stability

Anti-HBe Positive Control is stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 7 days at 2~8°C if capped tightly to avoid microbial contamination.

### Materials required but not provided

Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000.

Cat. No. Anti-HBe121 or Anti-HBe122: Antibody to Hepatitis B e Antigen (CLIA), 2×50 tests or 2×100 tests;

Cat. No. Anti-HBe221: Anti-HBe Calibrators, 1×2.0mL for each level of calibrator C0, C1;

Cat. No. Anti-HBeN321 or Anti-HBeN322: Anti-HBe Negative Control, 3×2.0mL or 6×2.0mL;

Cat. No. CS511: Substrate Solution, 4×115mL;

Cat. No. WB411: Wash Buffer, 1×10L;

Cat.No.115-025456-00: Reaction Vessel.

### Assay Values

The control values (the target value and ranges) determined by Wiener lab. standard procedure and routine method are listed in the target value sheet. The target value was obtained from Wiener lab. measurement system and the range was calculated as the target values ± 3 standard deviations. The control value is lot-specific. Please check the lot number of each control before use.

### Quality Control

Use Wiener lab. Anti-HBe Positive Control to monitor the status of each calibration. The value of the calibrator is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. The control result should be within the defined range as shown in the control target value sheet. If the control is out of the range, the measurement system and reagents should be checked. For instance, expiration date or storage condition of the calibrators, reagent, and controls; possible contamination of reagent; position of reagent; sample placement on the analyzer; the parameter setting in software; and operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and measurement procedure.

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<b>Anti-HBe Positive Control</b>		<b>Wiener lab.</b>
LOT 	2-8°C 2 mL	

WIENER Laboratorios S.A.I.C.

  
 Dra. VIVIANA E. CETOLA  
 DIRECTORA TÉCNICA

# Anti-HBe Negative Control

## Order Information

Catalog No.	Package Size
Anti-HBeN321	3x2.0 mL/vial
Anti-HBeN322	6x2.0 mL/vial

## Intended Use

Anti-HBe Negative Control is used for monitoring the accuracy and precision in the qualitative measurement of Anti-HBe by Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000. It is also intended for clinical laboratories to assess their capability in quality control of lab work.

## Summary

Anti-HBe Negative Control has similar components as the test samples. Therefore, the control can be tested by following the same procedure as the specimen, measured value of the control is compared to the target value and range to see if the control meets its specifications.

## Components

Anti-HBe Negative Control is prepared from Anti-HBe negative human Plasma with 0.09% sodium azide and 0.05% ProClin 300 as the preservative.

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- The control values are lot-specific with each reagent pack as listed in the value sheet.

- The control should be run daily in parallel with the patient samples, after each calibration, or after the reagent lot has been changed.
- Please perform maintenance and standard operation as scheduled, including calibration and system check-up to ensure the performance of the measurement system.
- Raw material for control production has been tested with Chinese Food and Drug Administration-approved kits for antibodies to HIV-1, HIV-2, HCV, Syphilis TP, and for HBsAg. The results are all negative.
- As no testing method can rule out the potential risk of infection with absolute certainty, this material should be handled as potential biohazard<sup>1</sup>.
- If there is any evidence of microbial contamination or excessive turbidity in the product, discard the vial.
- Material safety data sheet is available upon request.
- Disposal of all waste material in compliance with local guidelines.
- The target value and range of Wiener lab. Anti-HBe Negative Control are for reference only. Each laboratory should establish its own internal quality control procedure and set up its own target value and range.
- Sodium azide may react with tube fittings made of lead or copper, generating highly explosive metal azide. Please wash thoroughly with a large volume of water after disposal of these liquid to prevent azide accumulation<sup>2</sup>.

## Preparation

- Take the vial out of 2~8 °C refrigerator and standing on bench for no more than 30 min to reach room temperature.
- Mix contents by gently inverting the vial multiple times before use. Avoid bubble formation.
- Dispense the required volume into a sample cup, and then run QC test and analyze it in the same way as patient samples. Control can be set up manually or imported and the system automatically analyzes information from the bar code.

## Storage and Stability

Anti-HBe Negative Control is stable up to the expiration date as indicated on the label, when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 7 days at 2~8°C if capped tightly to avoid microbial

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Wiener lab. CLIA series Chemiluminescence analyzers CLIA 1000 or CLIA 2000

Cat. No. Anti-HBe121 or Anti-HBe122 : Antibody to Hepatitis B e Antigen (CLIA), 2x50 tests or 2x100 tests;

Cat. No. Anti-HBe221: Anti-HBe Calibrators, 1x2.0mL for each level of calibrator C0, C1; Cat. No. Anti-HBeP321 or Anti-HBeP322: Anti-HBe Positive Control, 3x2.0mL or 6x2.0mL;

Cat. No. CS511: Substrate Solution, 4x115mL;

Cat. No. WB411: Wash Buffer, 1x10L;

Cat.No.115-025456-00: Reaction Vessel.

## Assay Values

The control values (the target value and range) determined by Wiener lab. standard measurement procedure are listed in the target value sheet. The target value was obtained from Wiener lab. measurement system and the range was calculated as the target value  $\pm$  3 standard deviations. The control value is lot-specific, Please check the lot number of each control before use.

## Quality Control

Use Wiener lab. Anti-HBe Negative Control to monitor the status of each calibration. The value of the calibrator is lot-specific on Wiener lab. CLIA series Chemiluminescence analyzers. The control result should be within the defined range as shown in the control target value sheet. If the control is out of the range, the measurement system and reagents should be checked. For instance, expiration date or storage condition of the calibrators, reagent, and controls; possible contamination of reagent; position of reagent; sample placement on the analyzer; the parameter setting in software; and operation status of the analyzer. Each laboratory should establish its own internal quality control scheme and measurement procedure.


## Graphical Symbols

In vitro diagnostic medical device Batch code Catalogue Number

European Conformity Authorized representative in the European Community

## References

- HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories.
- DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard.

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2016 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved