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1. Kingsnorth A, LeBlanc K: Hernias – inguinal and incisional. Seminar. Lancet 2003; 352: 1561–71
2. Arnd P K: Lichtenstein tension-free hernioplasty: Its inception, evolution, and principles. Hernia (2004); 8: 1–7
3. Bittner R, Leibl BJ, Jager C, Kraft B, Ulrich M, Schwarz J: TAPP – Stultjart technique and result of a large single center series. J Min Access Surg 2006; 2: 155–9
4. M. Korenkov, A. Paul, S. Sauerland et al., Classification and surgical treatment of incisional hernia: results of an experts' meeting. Langenbecks Arch Surg 386 (2001), pp. 65–73
5. Licheri S, Erdas E, Pisano G, et al. Chvrel technique for midline incisional hernia: still an effective procedure. Hernia. Published online: 31 October 2007
6. Schumpelick V, Klinge U, Junge K, Stumpf M. Incisional abdominal hernia – the open mesh repair. Langenbecks Arch Surg 389: 1–5

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LITERATURA

1. Kingsnorth A, LeBlanc K, Herrias. Inguinal and incisional. Seminar Lancet 2003; 362: 1561-71
2. Amid P K: Lichtenstein tension-free hernioplasty: Its inception, evolution, and principles. Hernia (2004); 8: 1-7
3. Bittner R, Leibl B, Jager C, Kraft G, Ulrich M, Schwarz J. TAPP - Stuttgart technique and result of a large single center series. J Min Access Surg 2006;2:155-9
4. M. Kerenkov, A. Paul, S. Sauerland et al. Classification and surgical treatment of incisional hernia: results of an experts' meeting. Langenbecks Arch Surg 386 (2001); pp. 65-73
5. Licheri S, Erdas E, Pisano G, et al. Chevrel technique for midline incisional hernia: still an effective procedure. Hernia. Published online: 31 October 2007
6. Schumpelick V, Klinge U, Junge K, Stumpf M. Incisional abdominal hernia: the open mesh repair. Langenbecks Arch Surg 389:1-5

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Typ	Rozměr (mm)	Katalogové číslo
PP MESH light	150 x 100	092-1510
PP MESH light	200 x 150	092-2015
PP MESH light	300 x 150	092-3015
PP MESH standard	100 x 100	092-3030
PP MESH standard	60 x 110	091-0611
PP MESH standard	150 x 100	091-1510
PP MESH standard	150 x 150	091-1515
PP MESH standard	200 x 150	091-2015
PP MESH standard	300 x 150	091-3015
PP MESH standard	300 x 300	091-3030
PP MESH large pore	150 x 100	093-1510
PP MESH large pore	200 x 150	093-2015
PP MESH large pore	300 x 150	093-3015
PP MESH large pore	300 x 300	093-3030
PP MESH extra large pore	150 x 100	094-1510
PP MESH extra large pore	200 x 150	094-2015
PP MESH extra large pore	300 x 150	094-3015
PP MESH extra large pore	300 x 300	094-3030

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INSTRUCTIONS FOR USE

Polypropylene Mesh for Surgical Applications PP MESH

PRODUCT DESCRIPTION

PP MESH is a mesh made of a polypropylene monofilament yarn i.e. of a non-absorbable synthetic material considered as a physically and chemically inert and stable material; it is non-toxic and does not show an immunogenic activity. Thanks to its successful use for many years, most authors consider polypropylene as a golden standard in the field of herniology.

MESH TYPES

PP MESH Light

Its weight per unit area is 40 g/m², it corresponds to the lighter form of polypropylene meshes. The outstanding features of the mesh are its good porosity (pores are greater than 0.8 mm) and elasticity that is close to the elasticity of the abdominal wall.

PP MESH Standard

Its weight per unit area is 50 g/m², it corresponds to the standard weight of polypropylene meshes. The outstanding features of the mesh are a higher porosity (pores are greater than 1.5 mm) and elasticity that is close to the elasticity of the abdominal wall.

PP MESH Large Pore

Its weight per unit area is 67 g/m², it corresponds to the standard weight of polypropylene meshes. The outstanding features of the mesh are its large pores (2.6x1.5) mm and an elasticity comparable to that of the abdominal wall.

PP MESH Extra Large Pore

Its weight per unit area is 47 g/m², it corresponds to the lighter form of polypropylene meshes. The outstanding features of the mesh are its extra large pores (3.0x2.8) mm and an elasticity comparable to that of the abdominal wall.

INTENDED USE

The structure of the PP MESH is optimally designed for permanent support, augmentation or complete replacement of aponeurotic structures of the abdominal wall (e.g. in the case of large ventral hernias) and it safely resists the pressures to which the abdominal wall is subjected. Thus the PP MESH is indicated for use in surgeries of all types of abdominal herniae.

INDICATIONS

Inguinal hernia surgeries

The producer recommends the PP MESH to be implanted in conformity with the current recommendations of leading experts of the European Hernia Society and the American Hernia Society. (1) Thus the described material is suitable for application in "open" inguinal hernioplasty with

SK

SK

GB

Medicenter
Valeria Martinez

LEANDRO OLSINA

Director Técnico
MEDICENTER

De Valeria Soledad Martínez
IF-2018-01169887-APN-DNPM#ANMAT

a mesh (Lichtenstein's operation). The producer recommends observing the principles of Lichtenstein's operative technique – an appropriate extent of preparation, sufficient mesh size ensuring an overlap of defect edges, suitable type of mesh fixation. It is recommended that the fixation sutures penetrate the mesh at least 10 mm from the edge of the mesh. (2) The PP MESH is further indicated also for laparoscopic TAPP operations. The producer recommends a sufficient overlap at the edges of the fascial defect mesh placement without unnecessary warping and wrinkles, adequate attachment of the mesh. (3)

Ventral hernia surgeries (incisional, umbilical, epigastric and the like)

The producer recommends the PP MESH to be implanted in conformity with the recommendations of the European Hernia Society (EHS). (4) The mesh is suitable for all types of extraperitoneal ventral hernioplasty, therefore the material should be in no case left in contact with the intraperitoneal organs of the abdominal cavity (small intestine, colon, stomach, etc.). The choice of the type of plastic surgery – onlay or sublay – depends on the operating surgeon's preference; the PP MESH is suitable for both modifications of hernioplasty of ventral hernia. To reduce the risk of hernia recurrence it is advisable to use a mesh of sufficient size, the overlap of the edges of the fascial defect should be at least 5 cm to each side, the surgeon should choose an adequate fixation of the material preventing warping or dislocation of the mesh, fixation sutures should be made of a non-absorbable material and should penetrate the mesh at least 10 mm from the edge of the mesh. (5, 6)

Recommendations for fixation

Fixation of the mesh is an important part of the surgery; a mesh left without fixation is exposed to the danger of dislocation or inadequate placement (mesh warping, formation of undesirable wrinkles, bent edges of the mesh, etc.). To fix the mesh it is necessary to use a suture made of a non-absorbable material of sufficient strength (2-0 to 0 or 1), a monofilament material is preferred. The distance of the suture position from the edge of the mesh: the suture shall penetrate the mesh at least 10 mm from its edge.

CONTRAINDICATIONS

The PP MESH should not be used if the operative field is contaminated, e.g. in case of an unplanned opening of the intestinal lumen in the course of preparation. The PP MESH must always be separated from the organs of the abdominal cavity by the peritoneum, therefore it is contraindicated for intraperitoneal placement.

ADJUSTMENT OF MESH SHAPE AND SIZE

According to the respective operative situation the PP MESH can be cut to the required size, the mechanical characteristics of the mesh apart from its original edge do not change and so they remain maintained.

STERILIZATION

The product is sterilized by ethylene oxide.
Do not use the mesh if the sterile package has been damaged.
THE PRODUCT CANNOT BE RE-STERILIZED BY NO MEANS!
The product is intended for one use only.

Note:

COUNCIL DIRECTIVE 93/42/EEC, Annex No. 1, article 8.3, medical devices delivered in a sterile state must be designed, manufactured and packed in a non-reusable pack and/or according to appropriate procedures to ensure that they are sterile when placed on the market and remain sterile, under the storage and transport conditions laid down, until the protective packaging is opened or damaged. For this reason the medical device must not be reused because opening of the original pack (nonreusable pack) means damage to the manufacturer's original pack. Once the original pack has been opened and the product taken out of this pack the sterility of the product is broken. Once the original pack has been opened the product must not be repackaged and resterilized. The manufacturer cannot guarantee that a contamination of the product by microorganisms will not occur after removal of the product from its original pack and during handling.

EXPIRATION PERIOD

Three years from the date of sterilization.
The expiration date is given on the package.

STORAGE AND HANDLING CONDITIONS

The PP MESH polypropylene mesh shall be stored in its original packaging in clean rooms at ambient-air temperature 15 °C to 25 °C.

DISPOSAL OF NON-USED MEDICAL DEVICE

Incineration: together with municipal waste.

LITERATURE

1. Kingsnorth A, LeBlanc K: Hernias - inguinal and incisional. Seminar. Lancet 2003; 362: 1561-71
2. Amid P. K: Lichtenstein tension-free hernioplasty: its inception, evolution, and principles. Hernia (2004) 8: 1-7
3. Bittner R, Leibl B, Jager C, Kraff B, Ulrich M, Schwarz J: TAPP - Stuttgart technique and result of a large single center series. J Min Access Surg 2006; 2: 155-9
4. M. Korenkov, A. Paul, S. Sauerland et al., Classification and surgical treatment of incisional hernia: results of an experts' meeting. Langenbecks Arch Surg 386 (2001), pp. 65-73
5. Uichen S, Erdas E, Pizarro G, et. al. Chevrel technique for midline incisional hernia: still an effective procedure. Hernia. Published online: 31 October 2007
6. Schumpelick V, Klinge U, Junge K, Stumpf M. Incisional abdominal hernia - the open mesh repair. Langenbecks Arch Surg 389:1-5

SIZE AND CATALOGUE CODE

Type	Size (mm)	Catalogue code
PP MESH light	150 x 100	091-1510
PP MESH light	200 x 150	092-2015
PP MESH light	300 x 150	092-3015
PP MESH light	300 x 300	092-3030
PP MESH standard	60 x 110	091-0911
PP MESH standard	150 x 100	091-1510
PP MESH standard	150 x 150	091-1515
PP MESH standard	200 x 150	091-2015
PP MESH standard	300 x 150	091-3015
PP MESH standard	300 x 300	091-3030
PP MESH large pore	150 x 100	091-1510
PP MESH large pore	200 x 150	091-2015
PP MESH large pore	300 x 150	091-3015
PP MESH large pore	300 x 300	091-3030
PP MESH extra large pore	150 x 100	094-1510
PP MESH extra large pore	200 x 150	094-2015
PP MESH extra large pore	300 x 150	094-3015
PP MESH extra large pore	300 x 300	094-3030

Other dimensions are produced on the basis of customer specification

EXPLANATORY NOTES ON THE SYMBOLS ON LABELS:	
Use by	Keep away from sunlight
REF Catalogue code	Temperature limitation
Sterilized by ethylenoxid	Do not use if packages damaged
Lot number (batch)	Waste bin
Attention, see instructions for use!	Plastics
Date of manufacture	Polyethylene terephthalate / polyethylene
Do not reuse	CE 1023 Conformity marking

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