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SD Biosensor, Inc.

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74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbukdo, 28161, REPUBLIC OF KOREA

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Digitally signed by Gestion Documental  
Electronica  
Date: 2023.07.04 18:09:58 -03:00

EN	REF	F-TNI-02
----	-----	----------

## STANDARD F TnI Pro FIA

STANDARD™ F TnI Pro FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

### KIT CONTENTS

### MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer

### SPECIMEN COLLECTION AND PREPARATION

#### Whole blood

- Collect the venous whole blood into the commercially available EDTA tube by venipuncture.
- It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored at room temperature and refrigerator conditions (2-8°C/36-46°F), the specimen can be used for testing within 8 hours after collection.
- Do not use hemolyzed blood specimens.


#### Serum

- Collect the whole blood into the commercially available plain tube. NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- Serum specimen may be stored at room temperature and refrigerator conditions (2-8°C/36-46°F) for 1 day prior to testing.

### TEST PROCEDURE

#### Preparation

- Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for the STANDARD F TnI Pro FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch and check the test device in the foil pouch.

#### Analysis of control

STANDARD TEST™ mode	STANDARD F200 and F2400 Analyzer
1. Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.	


STANDARD F2400 Analyzer	Workplace → "Run Test" → Insert patient ID and / or operator ID on the analyzer
STANDARD F200 Analyzer	Standard Test™ mode → Insert patient ID and / or operator ID

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

- Using a disposable dropper, collect the 100µl of the specimen (whole blood or serum) to the black line of the disposable dropper. Add the collected specimen into the extraction buffer and mix thoroughly.

- Apply the 100µl of specimen mixture to the specimen well of the test device.

- After applying the control mixture, immediately press the **TEST START** button.

### PERFORMANCE CHARACTERISTICS

#### Analytical performance

- Accuracy (Method comparison)

Results of the serum / whole blood measured by STANDARD™ F TnI Pro FIA meet the acceptance criteria as below:

Regression Analysis		
	SERUM	WHOLE BLOOD
Slope	0.9704	1.0022
Y-intercept	47.971	7.425
R	0.9955	0.9945
R <sup>2</sup>	0.9910	0.9890
n	210	210

- The analyzer will automatically display the test result after 10 minutes.

### INTERPRETATION OF TEST RESULTS

Troponin I concentrations were determined using serum specimens obtained from 732 apparently healthy individuals. The upper reference limit (99th percentile) for TnI was 70 ng/L, and the lowest TnI concentration with CV 10% was found to be 70 ng/L.

TnI Concentration (ng/L)			
Population	N	99 <sup>th</sup> PERCENTILE URL	90% CONFIDENCE INTERVAL
	732	70.0	54.0 - 81.0

Elevated troponin I specific to myocardial injury and may be due to conditions other than MI. Acute MI can be established if there is a rise and/or fall in cardiac troponin with at least one value above the 99th percentile URL, together with a high clinical likelihood or ECG evidence of myocardial ischemia.

The default results are measured in ng/L. Should an alternate result unit be desired, the user can change unit through the "Setting" menu of the analyzer or manually calculate using the following conversion factors:

(Concentration in ng/L) × (0.001) = ng/mL = µg/L  
(Concentration in ng/L) × (1) = pg/mL

The reportable range is 10 ng/L to 20,000 ng/L. If result is below 10 ng/L, it will be reported as < 10 ng/L. If result is above 20,000 ng/L, it will be reported as > 20,000 ng/L.

### EXPLANATION AND SUMMARY

#### Introduction

The Troponin (Tn), one of the subunits of the troponin regulatory complex, binds to actin and inhibits interactions between actin and myosin. Troponin complex plays an important role in the regulation of skeletal and cardiac muscle contraction. Troponin complex consists of three subunits: Troponin T (TnT), Troponin I (TnI) and Troponin C (TnC). The subunits are held together by non-covalent interactions. TnI is the troponin-binding subunit that regulates interaction between the troponin complex and the thin filament. The TnI subunit is responsible for inhibiting actomyosin formation at low intracellular Ca<sup>2+</sup> concentrations. The third subunit, TnC, binds Ca<sup>2+</sup> ions during the excitation of the muscle and changes the conformation of troponin complex, thus enabling the formation of actomyosin complex and the subsequent muscle contraction. Cardiac Troponin I (cTnI) has been known as reliable markers of cardiac cell death. It is now widely used for the diagnosis of acute myocardial infarction (AMI), postsurgery myocardium trauma, chemotherapy cardiotoxicity, as well as some other diseases related to cardiac muscle injury. Being expressed cTnI only in cardiac tissue, cardiac TnI have been the preferred biomarkers for myocardial infarction for a long time. STANDARD F TnI Pro FIA, employing immunofluorescent detection system with STANDARD F analyzer, provides significantly fast, easy and accurate system to identify cardiac Troponin I (cTnI) in serum and whole blood as an aid in the diagnosis of myocardial infarction and enables supportive treatment decisions.

#### Intended use

STANDARD F TnI Pro FIA is a fluorescence immunoassay for the quantitative determination of cardiac Troponin I (cTnI) levels in human serum and whole blood using STANDARD F Analyzers, manufactured by SD BIOSENSOR. This test is an *in vitro* diagnostic use and intended for use as an aid in the screening and monitoring of acute myocardial infarction (MI).

#### Test principle

STANDARD F TnI Pro FIA uses antibodies that specifically bind to the TnI protein. When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If Troponin I (TnI) is present, it will be bound to the detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F Analyzer. If it is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of the TnI in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

### KIT STORAGE AND STABILITY

STANDARD F TnI Pro FIA is only for *in vitro* diagnostics use only. Kit materials are stable until expiration date printed on the outer box. Do not freeze the kit.

### WARNINGS AND PRECAUTION

- STANDARD F TnI Pro FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this instructions for use before testing.
- STANDARD F TnI Pro FIA should be used with STANDARD F Analyzer.
- STANDARD F TnI Pro FIA should remain in its original sealed pouch until ready to use. Do not use the test device if pouch is damaged or the seal is broken.
- STANDARD F TnI Pro FIA test device is single use only. Do not re-use it.
- Do not use hemolyzed specimens or frozen specimens.
- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating Desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Check the expiration date printed at the pouch or package.
- Check the volume (100µl) of extraction buffer.
- Use the STANDARD F TnI Pro FIA at room temperature (15-32°C/59-89.6°F).
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

### LIMITATION OF THE TEST

- The test should be used for the detection of Troponin I levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/ or produce invalid results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

### QUALITY CONTROL

#### Calibration

The calibration set test of STANDARD F Analyzer should be conducted according to the analyzer's manual.

#### When to use calibration set

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of an analyzer and test device.

#### How to use calibration set

Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.

- Select the 'Calibration' set in main menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

The STANDARD F Analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in 'STANDARD TEST™' mode. If 'EE' message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the 'EE' message still appears.

#### Internal quality control

- The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.
- The invalid result denotes that the fluorescence signal is not within the pre-set range. If the screen of STANDARD F Analyzer shows 'Invalid', turn off and turn on the analyzer again and re-test with a new test device.

#### External quality control

- Quality control testing should be run to check the performance of STANDARD F TnI Pro FIA and STANDARD F Analyzer. STANDARD F TnI Control manufactured by SD BIOSENSOR can be used for the external quality control. Control test should be conducted in accordance with the instructions of STANDARD F TnI Control.
- External quality control test should be run:
- once for new lot
  - once for each untrained operator
  - as required by test procedures in the Instructions for use of STANDARD F TnI Control and in accordance with local, state and federal regulations or accreditation requirements

ES	REF	F-TNI-02
----	-----	----------

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