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1105 Autoroute Chomedey
Laval, Canada
www.pegamedical.com

Authorized European Representative:
Syrimed - Z.A. de l'Angle 19370
Chamberet, France www.syrimed.com

NON STERILE

CE0413

The Fassier-Duval Telescopic IM System
El Sistema Telescópico Intramedular De Fassier-Duval
Le Systeme Intramedullaire Telescopique Fassier- Duval

4.8 SPS

REF : FD-048(SPS)-SS
Size: Ø 4.8mm X 404mm Lg.
LOT F: 281006-01 M: 280229-11 QTY: 1
INOX 316L/ACERO INOXIDABLE 316L/ss316L

1105 Autoroute Chomedey
Laval, Canada
www.pegamedical.com

Authorized European Representative:
Syrimed - Z.A. de l'Angle 19370
Chamberet, France www.syrimed.com

NON STERILE

CE0413


NOVAX DMA S.A.
DANIEL FIZ
PRESIDENTE

Novax DMA S.A.

[Handwritten signature]
NOVAX DMA S.A.
CALLE FRAGA 923
C1427BTS
BUENOS AIRES
MAY 2008

The Fassier-Duval Telescopic IM System
El Sistema Telescópico Intramedular De Fassier-Duval
Le Systeme Intramedullaire Telescopique Fassier- Duval

REF : FD-040(S)-SS
Size: Ø 4.0mm X 333mm Lg.
Lot: F: 280529-01 M: 281113-01 QTY: 1
INOX 316L/ACERO INOXIOABLE 316L/ss316L

Pega Medical

1105 Autoroute Chomedey
Laval, Canada
www.pegamedical.com

Authorized European Representative:
Synimed - Z. A. de l'Angle 19370
Chamberet, France www.synimed.com


NON STERILE

CE0413

The Fassier-Duval Telescopic IM System
El Sistema Telescópico Intramedular De Fassier-Duval
Le Systeme Intramedullaire Telescopique Fassier- Duval

REF : FD-032(L)-SS
Size: Ø 3.2 mm X 202mm Lg.
Lot: F: 281205-01 M: 280529-05 QTY: 1
INOX 316L/ACERO INOXIOABLE 316L/ss316L

Pega Medical

1105 Autoroute Chomedey
Laval, Canada
www.pegamedical.com

Authorized European Representative:
Synimed - Z. A. de l'Angle 19370
Chamberet, France www.synimed.com

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


NOVAX DMA S.A.
DANIEL FIZ
PRESIDENTE

Novax DMA S.A.

Pega Medical Inc., 1105 Autoroute Chomedey, Laval, Canada

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Device Description

The Fassier-Duval Telescopic IM System is a self ext- bone fractures. The design of the nail includes a fe to the proximal cortex) and a male component (whi Anchorage of the components is achieved through built-in feature that allows for extension of its len- normal patient growth occurs.

The nail is manufactured in both Ti alloy (Ti6Al4V AS Steel (316L, ASTM 138). It is available in five diameters. 420 mm in length.

Intended Use

This implant is indicated as a temporary implant to aid in the healing of long diaphysis fractures, osteotomies, malunions and nonunions and to prevent further fractures in femur, tibia and humerus in pediatric patients suffering from Osteogenesis Imperfecta without disrupting the bone growth plate. It can be used in procedures such as bone lengthening/shortening concomitantly with external fixators in pediatric or small stature patients with limb length discrepancy.

Contraindications

- Active or latent infection
- Osteoporosis, insufficient quality or quantity of bone/soft tissue
- Material sensitivity
- Sepsis
- Patients with abnormal neurological or mental conditions

Warnings

- Device breakage or damage can occur when implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing
- Improper insertion of the device during implantation can increase the possibility of loosening or migration


ement of the nails is recommended. The surgical allowed.

ization Instructions For Implant Com 115

used for cleaning and sterilizing items supplied must be cleaned and sterilized before initial use and after each subsequent use. After use, brush thoroughly using a broad spectrum cleaning and disinfecting agent; rinse with distilled water. Steam autoclave sterilization following the HIMA/AORN recommended practices for in-hospital sterilization is the only processing method suggested for the Fassier-Duval Telescopic IM System. All metallic implants and instruments can be steam sterilized following the parameters listed below:

Method	Steam	Steam
Sterilization type	Gravity	Prevacuum
Preconditioning pulse	—	3
Minimal temperature	250°F (121°C)	273°F (134°C)
Minimal cycle time	60 minutes ± 5 min.	20 minutes
Minimal drying time	10 minutes	30 minutes

It is recommended to dry and/or cool the parts to prevent condensation or perform a drain cycle after the steam cycle. Other sterilization methods and cycles may also be suitable. However, validation of any alternative method using appropriate laboratory techniques is advised. Never reuse surgical implants. No metal implants should ever be reimplanted. Although appearing undamaged, the device may have small defects or internal stresses which may eventually lead to implant failure. Correct implant handling is extremely important. Avoid contouring of metallic implants. Notching or scratching should be avoided when handling the device.

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© 2000-2009 Pega Medical Inc.
NCVAX D M A S.A.