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
Sirtex Medical (Europe) Limited
Walter-Flex Straße 2
53113 Bonn
Germany

Sirtex Medical Singapore Pte Ltd
50 Science Park Road
#01-01 The Kendall
Singapore Science Park II
Singapore 117406

EU Authorised Representative
Priory Analysts Ltd
The Atrium, Park St West, Luton, LU1 3BE, UK
Tel: 44 (0) 1908 319295
Fax: 44 (0) 1908 320776

SIRTEX

SIR-Spheres®
(Yttrium-90 Microspheres)

0086

1. DESCRIPTION

SIR-Spheres® consists of biocompatible microspheres containing yttrium-90 with a size between 20 and 60 microns in diameter. Yttrium-90 is a high-energy beta-emitting isotope with no primary gamma emission. The maximum energy of the beta particles is 2.27MeV with a mean of 0.93MeV. The maximum range of emissions in tissue is 11mm with a mean of 2.5mm. The half-life is 64.1 hours. In therapeutic use, requiring the isotope to decay to infinity, 94% of the radiation is delivered in 11 days. SIR-Spheres® is a permanent implant. Each device is moist heat sterilized and for single use only.

SIR-Spheres® is provided in a vial with water for injection. Each vial contains 3GBq of yttrium-90 (at the time of calibration) to a total of 5ml water for injection. The vial is shipped in a 6.4mm thick, lead pot. The package consists of a crimp-sealed SIR-Spheres® glass vial within the lead pot, and a package insert within Type A packing bucket.

The vial and its contents should be stored inside its transportation container at room temperature. (15-25°C, 59-77°F).

The calibration date (for radioactive contents) and the expiry information are on the vial label. SIR-Spheres® is usable up to 24 hours from the time of calibration.

SIR-Spheres® is implanted into hepatic tumours by delivery via either the common hepatic artery or the right or left hepatic artery using a catheter or implanted port. SIR-Spheres® distributes non-uniformly in the liver due to the physiology of hepatic arterial flow, the tumour to normal liver ratio of the tissue vascularity, and tumour size. The rim of the tumour gets higher density per unit distribution of SIR-Spheres® than the normal liver. Once SIR-Spheres® is implanted into the liver, it is not metabolised or excreted and it stays permanently in the liver.

2. INTENDED USE

SIR-Spheres® is intended for implantation into hepatic tumours via the hepatic artery.

3. INDICATIONS FOR USE

SIR-Spheres® is indicated for the treatment of patients with advanced non-operable liver cancer.

4. ACCESSORIES

Accessories that may be used for the implant procedure include:

- SIR-Spheres® Delivery Set and Sirtex V-vial
- Delivery Box and V-vial Holder
- Syringe Shield

5. INSTRUCTIONS FOR USE

5.1 PATIENT SELECTION AND PRE-TREATMENT TESTING

Patients with non-resectable tumours may be considered for treatment with SIR-Spheres Patient selection for treatment with SIR-Spheres® requires a medical opinion that control of tumour within the liver will result in patient benefit.

Patient Tests Before Treatment with SIR-Spheres®

- The following tests are recommended before treatment:
- A hepatic angiogram should be performed to establish arterial anatomy of the liver.
 - A nuclear medicine break-through scan (Intrahepatic Technetium MAA Scan) to determine the percent lung shunting. This is performed through the hepatic artery catheter or implanted port.
 - Serologic tests of liver function should be performed to determine the extent of liver function damage.

Other imaging studies, such as chest x-ray, CT scan of chest and abdomen, abdominal ultrasound and a bone scan are recommended to determine the extent of disease.

5.2 RADIATION DOSIMETRY

The radiation dosimetry of the SIR-Spheres® can be a complex and difficult task due to the non-uniform distribution of the microspheres in the normal liver and the tumours. In general, 1 GBq (27 mCi) of Yttrium-90/g of tissue provides the equivalent of 50 Gy of radiation dose.¹ However, because of the non-uniform distribution of the dose between the tumour and the normal liver tissue, a proportionally larger amount of radiation will be delivered to the tumour tissue.

Example, a patient has a 1500g liver with a 4cm size tumour in the right lobe, and a 3cm size nodule in the left lobe. The technetium-99 scan shows a 5:1 density ratio per unit volume between the tumour and the liver. The patient receives 2 GBq of SIR-Spheres®. The radiation dose to the tumour is 294 Gy and to the liver tissue is 58.5 Gy.

The radiation dose for other organs would be minimal or negligible, except for the organs adjacent to the liver, such as the stomach, large intestine and gall bladder, and the lung. The radiation dose may increase significantly when there is shunting of the arterial blood to the lung, or if SIR-Spheres® are inadvertently delivered to other organs such as the stomach or pancreas.

5.3 TECHNIQUE FOR PERFORMING THE INTRA-HEPATIC TECHNETIUM MAA SCAN

Technetium-99 MAA (Macro-Aggregated Albumin) 150MBq (4mCi) is injected via the hepatic artery catheter or implanted port to assess the fraction that passes through the liver to the lungs and the relative distribution of MAA (and hence SIR-Spheres®) between tumour and normal liver. The patient is placed supine. Anterior and posterior images of the abdomen and thorax, and right lateral images of the abdomen are taken. Regions of interest are drawn around the whole of lung fields and the whole of liver field.

$$\text{The percent lung shunting} = \frac{\text{counts of total lung}}{\text{counts of total lung} + \text{counts of liver}} \times 100$$

This same technique can be used to calculate the relative distribution of MAA (and hence SIR-Spheres®) between tumour and normal liver. This is only possible when clearly defined tumour and normal liver 'areas-of-interest' can be determined on the MAA scan.

If percent lung shunting is >10% then there is need for dose reduction of SIR-Spheres® (see Table 1 below)

Table 1 – Dose Reduction Recommendations

Percent Lung Shunting	Activity of SIR-Spheres®
< 10%	Deliver full amount of SIR-Spheres®
10% to 15%	Reduce amount of SIR-Spheres® by 20%
15% to 20%	Reduce amount of SIR-Spheres® by 40%
> 20 %	Do not give SIR-Spheres®

5.4 CALCULATION OF INDIVIDUAL DOSE

There are three accepted methods for calculating the patient radiation dose, these being the empirical model, the BSA model and the patient model.

Empirical

The empirical model accepts the safety margins of the doses known from the previously published clinical data and can be determined according to Table 2.

Table 2 – The Recommended Patient Dose

The % Involvement by the Tumour in the Liver	Recommended Y-90 Dose*
> 50 %	3.0 GBq
25 % - 50 %	2.5 GBq
< 25 %	2.0 GBq

*When there is 10 % or more lung shunting, the patient dose should be reduced according to Table 1.

BSA

The BSA method varies yttrium-90 activity according to the size of the patient and the size of the tumour within the liver and is the most widely used method.

The BSA method uses the patient's Body Surface Area (BSA) (calculated from the patient's weight and height) and the percentage of the liver (by volume) that is replaced with tumour (calculated from the CT scan). Most patients will receive between 1.3-2.5GBq of yttrium-90 if the whole liver is to be treated.

¹ Russell, Carden, Herron: 'Dosimetry Calculations of Yttrium-90 used in the treatment of liver cancer.' Endocurietherapy/Mayotherm Oncol 1988;4:171-186

Date of Issue: November 2010 (CDP 179)

Juan J. D.
TECNONUCLEAR S.A.
ERNESTO S. FURNARI
PRESIDENTE

Activity of SIR-Spheres® in GBq

$$= (\text{BSA} - 0.2) + \left(\frac{\% \text{ tumour involvement}}{100} \right)$$

Body Surface Area (BSA) is calculated from a weight/height chart
 $\% \text{ tumour involvement} = \frac{\text{volume of tumour} \times 100}{\text{volume of tumour} + \text{liver}}$

Partition Model

This method involves selecting safe radiation doses to the normal liver and lung and implanting the maximum activity that will not exceed these limits. The radiation dose to the normal liver parenchyma should not exceed 80Gy in patients with normal liver and 70Gy in patients with cirrhosis. The dose to the lung should not exceed 25Gy and preferably be less than 20Gy. The dose received by the tumour has no upper limit.

The technique requires two measurements to be made:

1. Measurement of the volume of tumour and normal liver determined from a CAT scan
2. Measurement of the proportions of technetium-99 labelled MAA activity that lodges in the tumour, normal liver and lung as determined from a gamma scan.

As the lung is largely filled with air, the CAT scan cannot be used to measure the volume of the lung parenchyma, and hence an estimation of 1000cc is made. For the purpose of calculating tissue mass, all tissue densities are all estimated at 1gm/cc.

Equation 1 is used to calculate the radiation dose received by an organ after SIR-Spheres® has been delivered to that organ.

$$\text{Equation 1: } \text{Tissue Radiation Dose (Gy)} = \frac{49670 \times \text{Total yttrium-90 activity in the organ or tissue (in GBq)}}{\text{Mass of the organ or tissue (in grams)}}$$

Therefore, to calculate the activity to be implanted, it is necessary to

1. from the liver CT scan, calculate the volumes of the normal liver and tumour
2. convert each volume to mass on the basis of 1g/cc
3. from a lung CT scan, determine the volume of the lung and convert to mass (or estimate it as 1000g)
4. from the nuclear medicine break-through scan, determine the activity in the lung, tumour and liver
5. determine the T/N activity ratio (calculated as activity per unit mass of organ or tissue) using Equation 2
6. determine the percentage shunted to the lungs using Equation 3.

To determine the T/N ratio, the following equation should be used.

$$\text{Equation 2: } T/N = \frac{A_{\text{tumour}}/M_{\text{tumour}}}{(A_{\text{liver}}/M_{\text{liver}})}$$

Where
 A_{tumour} is the activity in tumour
 M_{tumour} is the mass of tumour
 A_{liver} is the activity in the normal liver
 M_{liver} is the mass of the normal liver

$$\text{Equation 3: } \text{Lung Activity (} A_{\text{lung}} \text{)} = \frac{A_{\text{total}} \times L}{100}$$

Where A_{total} is the total activity implanted for the nuclear medicine break-through scan and L is the percentage lung shunting.

$$\text{Percent lung shunting} = 100 \times A_{\text{lung}} / (A_{\text{lung}} + A_{\text{liver}} + A_{\text{tumour}})$$

To calculate the total activity to be implanted, use the following equations. The activity required should be calculated using the lung dose as the limiting factor, and then again using the normal liver dose as the limiting factor. The lower of the two activities calculated should be used.

To determine the activity implanted to accommodate a limiting lung dose:

$$\text{Equation 4: } A_{\text{total}} = \frac{D_{\text{lung}} M_{\text{lung}} / 100}{49670}$$

Where:
 D_{lung} is the dose to the lung
 M_{lung} is the mass of the lung
 A_{lung} is the activity to the lung
 A_{total} is the total activity
 L = the percentage lung shunting

To determine the activity implanted to accommodate a limiting normal liver dose:

$$\text{Equation 5: } A_{\text{total}} = \frac{(D_{\text{liver}} M_{\text{liver}}) / 100}{(49670 (1-L/100))}$$

Page 1 of 9

Farm. Vilma R. Ceraso
MAT. 10,050
DIRECTOR TECNICO
TECNONUCLEAR S.A.

Sirtex Medical Limited
Level 33, 101 Miller Street
North Sydney, NSW 2060
Australia
Tel: +61 2 9964 8400

Sirtex Medical (Europe) Limited
Walter-Flex Straße 2
53113 Bonn
Germany

Sirtex Medical Singapore Pte Ltd
50 Science Park Road
#01-01 The Kendall
Singapore Science Park II
Singapore 117406

EU Authorised Representative
Priority Analysts Ltd
The Atrium, Park St West, Luton, LU1 3BE, UK
Tel: 44 (0) 1908 319295
Fax: 44 (0) 1908 320776

The partition model can only be used where the tumour mass is a discrete area within the liver. This is more likely patients with Primary Hepatocellular Carcinoma (HCC) where there is often a large single tumour mass. Patients with metastatic disease usually have multiple areas of metastatic spread that precludes defining the tumour and normal parenchymal compartments.

5.5 DOSE PREPARATION PROCEDURE

- Unpack SIR-Spheres®, leaving shipping vial in lead pot.
- Place the bench top in a lead or acrylic shielded box if available.
- Remove the centre of aluminium seal from sterile v-vial with forceps, and clean the septum with an alcohol swab.
- Place the v-vial in an empty lead pot (10 cm x 6 cm) for stability and shielding.
- Insert a short 25gauge needle through the septum of the v-vial until it just pierces the septum to create a vent.
- Remove the SIR-Spheres® shipping vial from the lead pot and shake vigorously to disperse the SIR-Spheres®.
- Using a dose calibrator, determine the activity in the shipping vial and return it to the lead pot.
- Determine the volume to be withdrawn to provide the required patient radiation dose.
- Partially remove the aluminium seal of the SIR-Spheres® shipping vial, clean with alcohol swab.
- Insert a 25 gauge needle through the septum of the shipping vial to create a vent, ensuring the needle is well clear of the contents in the shipping vial.
- Use a shielded 5ml syringe with a 20-22 gauge spinal needle at least 70mm long to puncture the septum of the SIR-Spheres® shipping vial, and quickly draw back and forth several times in order to mix the SIR-Spheres® thoroughly.
- Quickly withdraw the pre-calculated patient radiation dose, and transfer into the vented v-vial in the other lead pot. Withdraw the required amount quickly before the contents of the shipping vial start to settle.
- Verify the patient dose in the v-vial by re-measuring the activity in the shipping vial with the dose calibrator, and correct, if necessary.
- Put the v-vial, containing the confirmed patient dose into the dedicated acrylic shield.

The patient dose is now ready for transport to the SIR-Spheres® implantation room.

5.6 IMPLANT PROCEDURE

[Doctors must refer to the Sirtex Medical User's Manual for delivering SIR-Spheres® before attempting to implant this device].

SIR-Spheres® can be implanted via the hepatic artery using an implanted catheter with port or trans femorally.

Hepatic Artery Port Implantation

This method is generally used if the port is being used for other treatment, such as regional hepatic perfusion chemotherapy. A surgeon who is totally familiar with this technique must undertake insertion of the hepatic artery port. Attention to small surgical details can have a dramatic effect on the success or complications of the procedure.

Several additional factors are to be noted if SIR-Spheres® is to be implanted through the port. These include:

- The hepatic artery catheter should be placed into the arterial supply of the liver so that all the liver is perfused by the catheter.
- There are frequently small arteries that pass from the common hepatic artery (and sometimes even from the right or left hepatic arteries) to the stomach and duodenum that must be ligated at the time of inserting the port/pump. Failure to ligate these vessels may result in SIR-Spheres® lodging in the stomach and duodenum at the time of implant and this may result in severe complications.
- The catheter is usually placed into the hepatic artery by inserting it through the gastro-duodenal artery, but may need to be placed into another artery.
- The catheter should be at least 0.8mm internal diameter. If smaller diameter catheters are used, they may block during the delivery of SIR-Spheres®.
- The gallbladder should always be removed to prevent SIR-Spheres® causing radiation necrosis of the gallbladder in conjunction with HAC.
- The patient must recover from any surgical operations before being treated with SIR-Spheres®.
- It is important to deliver the SIR-Spheres® slowly into the hepatic artery to prevent the microspheres refluxing back down the hepatic artery and lodging in the pancreas, stomach or other organs. The catheter should be flushed at regular intervals during the delivery procedure to ensure the microspheres do not block the catheter.
- If a pump has been inserted, SIR-Spheres® is delivered through the side port of the pump. In some types of pumps the side port can only be accessed with a gauge 24 or smaller needle. While SIR-Spheres® can be delivered through this small needle, there is an increased risk of the spheres clogging the needle. Therefore, the operator

should deliver a very dilute suspension of SIR-Spheres® to prevent clumping of the needle.

- If the pump does not have a separate side port, then it cannot be used to deliver SIR-Spheres®.

Transfemoral Implantation

The hepatic artery catheter is inserted via the femoral artery under X-ray guidance. If this is the preferred method of implantation, a trained interventional radiologist must perform the procedure.

This method allows complete control of exactly where the catheter is placed and allows routine checking of the catheter position throughout the implant procedure. A transfemoral catheter can also be inserted further into the liver and helps avoid the possibility of reflux of SIR-Spheres® into small arteries supplying the gut. This is not possible with implanted catheters attached to port/pumps.

The procedure for delivering the SIR-Spheres® is similar to using a port/pump once the catheter has been correctly sited and the end of the catheter is connected to the SIR-Spheres® Delivery Set. The radiologist must repeatedly check the position of the catheter during the procedure to ensure it remains correctly sited and that reflux of the SIR-Spheres® into other organs does not occur. This is performed by injecting contrast medium through the left hand port of the Delivery Set during the delivery of SIR-Spheres®.

SIR-Spheres® must be delivered slowly at a rate of no more than 5ml per minute as rapid delivery may cause reflux back down the hepatic artery and into other organs. At the conclusion of the procedure, the catheter is removed.

Radiological Placement of Catheter

The radiologist must be familiar with the frequent arterial abnormalities in the blood supply to the liver and from the liver to the gut. Every attempt should be made to deliver the SIR-Spheres® into the main hepatic artery so that radiation is distributed to both lobes of the liver. If the tumours are limited to one lobe, the catheter can be selectively inserted into the lobar artery supplying only that lobe, thus sparing the normal lobe.

It is essential that SIR-Spheres® not be delivered to other organs, in particular the pancreas, stomach or duodenum. The catheter must be placed well distal to the gastro-duodenal artery (GDA) and any other artery that is supplying blood to the gut in order to prevent SIR-Spheres® going to the duodenum and stomach. If there is any possibility of SIR-Spheres® passing down the GDA then the implantation must not proceed. It may be preferable to block the GDA with an intraluminal coil or other agent to prevent SIR-Spheres® from flowing to the duodenum. No harm will occur if the gastro-duodenal artery is blocked.

Note: Virtually all complications from SIR-Spheres® arise from the inadvertent delivery of SIR-Spheres® into small blood vessels that go to the pancreas, stomach or duodenum.

The apparatus can all be assembled on a steel tray and placed at the side of the patient. Directions for the use of the Delivery Set are included with the device. These directions should be read in their entirety prior to use.

6. RADIATION SAFETY

Regulatory and local radiation usage guidelines should be followed concerning implantation and post-implantation care.

Personnel Exposure (thermo luminescent dosimetry - TLD)

Table 3 – Typical Exposure Dose Per Patient for Implant Preparation (Technologist) for a 3GBq Device (30 minutes)

	Trunk mSv (rem)	Lens of the Eye mSv (rem)	Hands mSv (rem)
Shallow Dose (0.07mm)	0.027 (2.7)	0.026 (2.6)	0.35 (35)
Deep Dose (10 mm)	0.003 (0.3)	0.004 (0.4)	

Table 4 – Typical Exposure Dose Per Patient for Implant Procedure (Physician) for a 2GBq implant (20 minutes)

	Trunk mSv (rem)	Lens of the Eye mSv (rem)	Hands mSv (rem)
Shallow Dose (0.07mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

Post-Implant Exposure from patients implanted with an average of 2.1 GBq approximately 5-6 hours post implantation at various distances from the patient's abdomen: (1 mSv = 100 rem)

0.25m	18.8 µSv/hr
0.5m	9.2 µSv/hr
1m	1.5 µSv/hr
2m	0.4 µSv/hr
4m	<0.1 µSv/hr

7. ADVERSE EVENTS

The common adverse events after receiving the SIR-Spheres® result from a mild post-embolisation syndrome and include fever, mild to moderate abnormality of liver function tests (mild increase in SGOT, alkaline phosphatase, bilirubin), abdominal pain, nausea, vomiting, and diarrhoea.

Potential Serious Adverse Events Due to High Radiation

- **Acute pancreatitis** — causes immediate severe abdominal pain. Verify by gamma camera imaging of the abdomen and test for serum amylase.
- **Radiation Pneumonitis** — causes excessive nonproductive cough. Verify by X-ray evidence of pneumonitis.
- **Acute Gastritis** — causes abdominal pain. Verify by standard methods to diagnose gastritis/ulceration.
- **Radiation Hepatitis** — causes unexplained progressive deterioration of liver function. Verify by exclusion of other causes and core biopsy of the liver.
- **Acute cholecystitis** — causes significant pain and may require cholecystectomy for resolution.

8. WARNINGS

- Inadvertent delivery of SIR-Spheres® to the gastrointestinal tract or pancreas will cause acute abdominal pain, acute pancreatitis or peptic ulceration. This may occur with greater frequency if SIR-Spheres® is delivered via an implanted hepatic artery port where there is less control over catheter placement.
- High levels of implanted radiation and/or excessive shunting to the lung may lead to radiation pneumonitis. Excessive radiation to the normal liver parenchyma may result in radiation hepatitis.
- Inadvertent delivery of SIR-Spheres® to the gallbladder may result in abdominal pain and cholecystitis that may require a cholecystectomy for resolution.

9. CONTRAINDICATIONS

- SIR-Spheres® is contraindicated in patients who have:
- had previous external beam radiation therapy to the liver;
 - ascites or are in clinical liver failure;
 - markedly abnormal synthetic and excretory liver function tests (LFTs);
 - greater than 20% lung shunting of the hepatic artery blood flow determined by Technetium MAA scan;
 - pre-assessment angiogram that demonstrates abnormal vascular anatomy that would result in significant reflux of hepatic arterial blood to the stomach, pancreas or bowel;
 - been treated with capecitabine within the two previous months, or who will be treated with capecitabine at any time following treatment with SIR-Spheres®;
 - complete portal vein thrombosis.

10. PRECAUTIONS

- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- A SPECT scan of the upper abdomen may be performed immediately after implantation of SIR-Spheres®. The SPECT scan will detect the Bremsstrahlung radiation from the yttrium-90 to confirm placement of the microspheres in the liver.
- This product is radioactive. Local regulations must be followed when handling this device.
- Some patients may develop gastritis following treatment. Gastric acid blocking drugs may be used the day before implantation of SIR-Spheres® and continued as needed to reduce gastric complications.
- Many patients may experience abdominal pain immediately after administration of SIR-Spheres® and pain relief may be required.
- SIR-Spheres® demonstrated a mild sensitisation potential when tested dermally in an animal model.

Sirtex Medical Limited
Level 33, 101 Miller Street
North Sydney, NSW 2060
Australia
Tel: +61 2 9964 8400

Sirtex Medical (Europe) Limited
Walter-Flex Straße 2
53113 Bonn
Germany

Sirtex Medical Singapore Pte Ltd
50 Science Park Road
#01-01 The Kendall
Singapore Science Park II
Singapore 117406

EU Authorised Representative
Priority Analysts Ltd
The Atrium, Park St West, Luton, LU1 3BE, UK
Tel: 44 (0) 1908 319295
Fax: 44 (0) 1908 320776

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Nombre del fabricante: Sirtex Wilmington LLC
Se extiende a Tecnonuclear S.A. el Certificado PM-2058-1, en la Ciudad de