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GenoType MTBDRsl

VER 2.0

for identification of the *M. tuberculosis* complex and its resistance to fluoroquinolones and aminoglycosides/cyclic peptides from sputum specimens or cultivated samples

0078/1

Store Kit Component 2
(included in box)
immediately
upon receipt at -20°C

CE

REF

31796A

1/2

2°C - 8°C

2/2

-20°C - -18°C

GenoType MTBDRsl

VER 2.0

1/2

2°C - 8°C

2/2

-20°C - -18°C

Store Kit Component 2
(included in box)
immediately
upon receipt at -20°C

 Hain Lifescience GmbH | Hardwiesenstr. 1 | 72147 Nehren | Germany
www.hain-lifescience.de | E-Mail: info@hain-lifescience.de | Tel.: +49 (0) 74 73- 94 51- 0

GenoType MTBDRsl VER 2.0		IVD	HAIN LIFESCIENCE
LOT AAXxxxxxx	YYYY-MM-DD	YYYY-MM-DD	
1/2	96	STRIPS	membrana strips, Membranstreifen, bandelettes, strisce di membrana, tiras de membranas, tiras de membrana, paski membranowe
	1.92 ml	DEN	denaturation solution, Denaturierungsreagent, solution de dénaturation, soluzione di denaturazione, solució de desnaturalizació, solució de desnaturalaçó, roztwór denaturacyjny
	96 ml	HYB	hybridization buffer, Hybridisierungspuffer, solution d'hybridation, tampone di ibridazione, tampón de hibridación, tampão de hibridação, bufor hybridyzacyjny
	96 ml	STR	stringent wash solution, Stringent-Waschlösung, solution de lavage stringent, soluzione di lavaggio stringente, solució de lavado asstringent, solució de lavagem asstringente, roztwór płuczacy stringent
	288 ml	RIN	rinse solution, Rinse-Lösung, solution de rinçage, soluzione di risciacquo, solució de aclarado, solució rinse, roztwór płuczacy
	960 µl	CON-C	conjugate concentrate, Konjugat-Konzentrat, conjugué concentré, conjugato concentrato, conjugado concentrado, stężony koniugat
	96 ml	CON-D	conjugate buffer, Konjugat-Puffer, tampón conjugué, tampone di diluizione del conjugato, tampón del conjugado, tampão do conjugado, bufor koniugatu
	960 µl	SUB-C	substrate concentrate, Substrat-Konzentrat, substrat concentré, substrato concentrato, substrato concentrado, stężony substrat
	96 ml	SUB-D	substrate buffer, Substrat-Puffer, tampón substrat, tampone di diluizione del substrato, tampón del substrato, tampão do substrato, bufor substratu
2/2	960 µl	AM-A	amplification mix A, Amplifikations-Mix A, mélange d'amplification A, mix di amplificazione A, mezcla de amplificación A, mistura de amplificação A, meszanina de amplifikacji A
	3.36 ml	AM-B	amplification mix B, Amplifikations-Mix B, mélange d'amplification B, mix di amplificazione B, mezcla de amplificación B, mistura de amplificação B, meszanina de amplifikacji B

001A/1

DEN: Warning/Achtung/Attention/Attenzione/Atención/Atenção/Uwaga
contains/enthält/contient/contiene/contém/zawiera: <2% NaOH
H315-319, P280-305+351+338-313

HYB, SUB-C: EUH210

Hain Lifescience GmbH | Hardwiesenstr. 1 | 72147 Nehren | Germany | Tel.: +49 (0) 74 73- 94 51- 0

002B/3

MARISOL MASINO
BIOQUIMICA M.N. 9483
DT - TECNOLAB S.A.

Director Técnico
Firma y Sello

IFU-317A-04

Valid edition of the instructions for use (IFU): see above
If you do not have this edition, you can obtain it from our website or by calling (see below).

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Performance Characteristics

The performance evaluation of the **GenoType MTBDRsl** VER 2.0 was carried out according to the instructions for use on hand.

Diagnostic performance

1. Clinical specimens

Diagnostic performance characteristics of the **GenoType MTBDRsl** VER 2.0 were determined in a study with 352 sputum specimens. The study specimens were collected in a high MDR-TB burden country. Both untreated patients as well as patients with previous or current anti-TB-treatment were included in the study.

The **GenoType MTBDRsl** VER 2.0 was compared to culture (successful cultivation on Loewenstein-Jensen solid medium or in MGIT [BD Diagnostics, Franklin Lakes, USA] and subsequent *M. tuberculosis* complex (MTBC) identification using the **GenoType Mycobacterium CM** VER 1.0). For discrepant culture-negative specimens, the result of the **GenoType MTBDRplus** VER 2.0 was used as reference method.

Furthermore, the **GenoType MTBDRsl** VER 2.0 was compared to conventional drug susceptibility testing (DST) using BACTEC MGIT 960 [BD Diagnostics, Franklin Lakes, USA] and Loewenstein-Jensen proportion method. Specimens with discrepant results were examined by sequencing the **GenoType MTBDRsl** VER 2.0 amplification region.

Additionally, all samples were examined by microscopy.

DNA extraction from NALC-NaOH-decontaminated sputum specimens was performed with the **GenoLyse**® kit according to the instructions for use.

21 specimens were excluded due to ambiguous DST results or contaminated cultures.

For the detection of MTBC, test results were rated true-positive if the result of **GenoType MTBDRsl** VER 2.0 was consistent with an MTBC-positive result of culture/**GenoType Mycobacterium CM** VER 1.0 or, in case of a culture-negative specimen, if a TB infection was indicated by an MTBC-positive **GenoType MTBDRplus** VER 2.0 result from the respective clinical specimen.

Table 1: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of MTBC from sputum specimens compared to culture/**GenoType Mycobacterium CM** VER 1.0 (GT Myco CM) or to culture/GT Myco CM and **GenoType MTBDRplus** VER 2.0 (GT MTBDRplus V2) result from the respective clinical specimens

	GenoType MTBDRsl VER 2.0	Culture/GT Myco CM		Sens: 98.8% Spec: 89.6% PPV: 96.9% NPV: 95.8%	Culture/GT Myco CM + GT MTBDRplus V2		Sens: 98.9% Spec: 100% PPV: 100% NPV: 95.8%	
		positive	negative		positive	negative		
total		positive	251	8		positive	259	0
		negative	3	69		negative	3	69
smear-positive		positive	232	8	Sens: 99.6% Spec: /* PPV: 96.7% NPV: /*	positive	240	0
		negative	1	2		negative	1	2
smear-negative		positive	19	0	Sens: 90.5% Spec: 100% PPV: 100% NPV: 97.1%	positive	19	0
		negative	2	67		negative	2	67

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value

* no value due to low sample number

For evaluation of resistance detection, the 251 MTBC-positive samples (positive both in culture and with **GenoType MTBDRsl** VER 2.0) were used. Test results were rated true-positive if the result of **GenoType MTBDRsl** VER 2.0 was consistent with the DST result or, in case of divergent results, if the result of **GenoType MTBDRsl** VER 2.0 was confirmed by sequencing from culture material of the respective sample.

Table 2: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of fluoroquinolone (FLQ) resistance from sputum specimens compared to culture/DST (tested with ofloxacin) or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 93.1% Spec: 100% PPV: 100% NPV: 98.0%	Culture/DST + sequencing		Sens: 96.4% Spec: 100% PPV: 100% NPV: 99.0%
	FLQ-R	FLQ-S		FLQ-R	FLQ-S	
	FLQ-R	54	0	FLQ-R	54	0
	FLQ-S	4	193	FLQ-S	2	195

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;

FLQ-R, fluoroquinolone-resistant; FLQ-S, fluoroquinolone-sensitive

Table 3: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of amikacin (AMK) resistance from sputum specimens compared to culture/DST or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 94.7% Spec: 98.1% PPV: 90.0% NPV: 99.1%	Culture/DST + sequencing		Sens: 97.3% Spec: 98.1% PPV: 90.0% NPV: 99.5%
	AMK-R	AMK-S		AMK-R	AMK-S	
	AMK-R	36	4	AMK-R	36	4
	AMK-S	2	209	AMK-S	1	210

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;

AMK-R, amikacin-resistant; AM-S, amikacin-sensitive

Table 4: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of capreomycin (CAP) resistance from sputum specimens compared to culture/DST or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 90.0% Spec: 98.1% PPV: 90.0% NPV: 98.1%	Culture/DST + sequencing		Sens: 97.3% Spec: 98.1% PPV: 90.0% NPV: 99.5%
	CAP-R	CAP-S		CAP-R	CAP-S	
	CAP-R	36		4	CAP-R	
CAP-S	4	207	CAP-S	1	210	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
CAP-R, capreomycin-resistant; CAP-S, capreomycin-sensitive

Table 5: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of kanamycin (KAN) resistance from sputum specimens compared to culture/DST or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 94.8% Spec: 96.6% PPV: 92.4% NPV: 97.7%	Culture/DST + sequencing		Sens: 98.6% Spec: 96.6% PPV: 92.4% NPV: 99.4%
	KAN-R	KAN-S		KAN-R	KAN-S	
	KAN-R	73		6	KAN-R	
KAN-S	4	168	KAN-S	1	171	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
KAN-R, kanamycin-resistant; KAN-S, kanamycin-sensitive

2. Culture samples

The diagnostic performance characteristics of the **GenoType MTBDRsl** VER 2.0 were determined in a study with 100 MTBC-positive culture samples.

The study samples obtained from a culture collection comprised MTBC strains originating from high MDR-TB burden countries as well as from low MDR-TB burden countries.

The **GenoType MTBDRsl** VER 2.0 was compared to conventional drug susceptibility testing (DST) using BACTEC MGIT 960 (BD Diagnostics, Franklin Lakes, USA) and Loewenstein-Jensen proportion method. Specimens with discrepant results were examined by sequencing the **GenoType MTBDRsl** VER 2.0 amplification region.

DNA extraction was performed using the **GenoLyse**[®] kit according to the instructions for use.

For kanamycin (KAN), DST was performed retrospectively using frozen culture aliquots that were recultivated. 89/100 samples yielded results; 11 samples could not be recultivated and were therefore excluded from the evaluation of KAN resistance.

Table 6: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of fluoroquinolone (FLQ) resistance from culture samples compared to culture/DST (tested with ofloxacin) or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 93.9% Spec: 98.5% PPV: 96.9% NPV: 97.1%	Culture/DST + sequencing		Sens: 100% Spec: 100% PPV: 100% NPV: 100%
	FLQ-R	FLQ-S		FLQ-R	FLQ-S	
	FLQ-R	31		1	FLQ-R	
FLQ-S	2	66	FLQ-S	0	68	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
FLQ-R, fluoroquinolone-resistant; FLQ-S, fluoroquinolone-sensitive

Table 7: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of amikacin (AMK) resistance from culture samples compared to culture/DST

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 100% Spec: 100% PPV: 100% NPV: 100%
	AMK-R	AMK-S	
	AMK-R	35	
AMK-S	0	65	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
AMK-R, amikacin-resistant; AM-S, amikacin-sensitive

Table 8: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of capreomycin (CAP) resistance from culture samples compared to culture/DST or to culture/DST and sequencing

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 84.6% Spec: 100% PPV: 100% NPV: 91.0%	Culture/DST + sequencing		Sens: 100% Spec: 100% PPV: 100% NPV: 100%
	CAP-R	CAP-S		CAP-R	CAP-S	
	CAP-R	33		0	CAP-R	
CAP-S	6	61	CAP-S	0	67	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
CAP-R, capreomycin-resistant; CAP-S, capreomycin-sensitive

Table 9: Performance characteristics of the **GenoType MTBDRsl** VER 2.0 for detection of kanamycin (KAN) resistance from culture samples compared to culture/DST or to culture/DST and sequencing.

GenoType MTBDRsl VER 2.0	Culture/DST		Sens: 94.4% Spec: 98.1% PPV: 97.1% NPV: 96.3%	Culture/DST + sequencing		Sens: 100% Spec: 100% PPV: 100% NPV: 100%
	KAN-R	KAN-S		KAN-R	KAN-S	
	KAN-R	34		1	KAN-R	
KAN-S	2	52	KAN-S	0	54	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;
KAN-R, kanamycin-resistant; KAN-S, kanamycin-sensitive

Further diagnostic performance characteristics of the **GenoType MTBDRsl** VER 2.0 have been published within the scope of an international multicenter study [21].

Analytical performance

Analytical specificity

The specificity of the **GenoType MTBDRsl** VER 2.0 is ensured by the accurate design of specific primers and probes which considers, among others, homology comparisons of the sequences published in gene databases, and by stringent reaction conditions.

The analytical specificity was determined with eight *M. tuberculosis* complex strains: *M. tuberculosis*, *M. africanum*, *M. bovis* BCG, *M. bovis* subsp. *bovis*, *M. bovis* subsp. *caprae*, *M. canettii*, *M. microti*, and *M. pinnipedii* (all sensitive to fluoroquinolones (FLQ) and to aminoglycosides/cyclic peptides (AG/CP)). The following 40 strains not detectable with the test system were also analyzed: *Bordetella pertussis*, *Corynebacterium spec.*, *Escherichia coli*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Mycobacterium abscessus*, *M. alvei*, *M. asiaticum*, *M. avium*, *M. celatum*, *M. chelonae*, *M. fortuitum*, *M. gastri*, *M. genavense*, *M. goodii*, *M. gordonae*, *M. haemophilum*, *M. immunogenum*, *M. interjectum*, *M. intermedium*, *M. intracellulare*, *M. kansasii*, *M. lentiflavum*, *M. mageritense*, *M. malmoense*, *M. marinum*, *M. mucogenicum*, *M. peregrinum*, *M. scrofulaceum*, *M. shimoidei*, *M. simiae*, *M. smegmatis*, *M. szulgai*, *M. triplex*, *M. ulcerans*, *M. xenopi*, *Nocardia spec.*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, and *S. pneumoniae*.

The eight *M. tuberculosis* complex isolates were correctly identified as FLQ- and AG/CP-sensitive MTBC strains. All other 40 isolates displayed invalid band patterns. Hence, an analytical specificity of 100% was achieved.

Analytical sensitivity (limit of detection, LOD)

To determine the LOD of the **GenoType MTBDRsl** for clinical samples, three BCG culture dilutions (FLQ- and AG/CP-sensitive, 1500, 150, and 15 CFU/ml) were prepared in triplicate. Including a negative control, DNA was extracted using the **GenoLyse**[®] kit and analyzed with the **GenoType MTBDRsl** applying the "MDR DIR" PCR protocol. An LOD of 150 CFU/ml was determined.

To determine the LOD of the **GenoType MTBDRsl** for culture samples, three BCG culture dilutions (FLQ- and AG/CP-sensitive, 1.65x 10⁶, 1.65x 10⁵, and 1.65x 10⁴ CFU/ml) were set up in triplicate. Including a negative control, DNA was extracted using the **GenoLyse**[®] kit and analyzed with the **GenoType MTBDRsl** applying the "MDR CUL" PCR protocol. An LOD of 1.65x 10⁵ CFU/ml was determined.

Reproducibility

In order to determine the intra-assay precision of the **GenoType MTBDRsl**, three BCG culture dilutions (FLQ- and AG/CP-sensitive; one above, one at, and one below the LOD) and one negative control were set up in four parallels and tested under identical conditions applying the "MDR DIR" PCR protocol. DNA extraction was performed using the **GenoLyse**[®] DNA extraction kit. All parallels showed identical and correct banding patterns and comparable signal strengths. Additionally, signal strengths between different sample dilutions were comparable. Hence, an intra-assay precision of 100% was achieved.

In order to determine the inter-assay precision of the **GenoType MTBDRsl**, three BCG culture dilutions (one above, one at, and one below the LOD) and a negative control were tested at three different points in time. Apart from the varied parameter, all other testing conditions were identical. DNA extraction was performed using the **GenoLyse**[®] DNA extraction kit and the isolates were analyzed with the **GenoType MTBDRsl** applying the "MDR DIR" PCR protocol. No deviations were detected between parallel samples, that is between runs banding patterns were identical and correct, and signal strengths were comparable. Hence, the inter-assay precision was 100%.

Interfering substances

There are substances that may inhibit PCR reactions. Such inhibitors may, for example, originate from the culture medium. In order to assess if the medium influences the **GenoType MTBDRsl**, 6 different *M. tuberculosis* complex samples (3x FLQ- and AG/CP-sensitive, 2x FLQ-sensitive and AG/CP-resistant, 1x FLQ- and AG/CP-resistant) were cultured in 4 different media (solid media: Loewenstein-Jensen, Stonebrink, and Middlebrook-7H10; liquid medium: MGIT (BD Diagnostics, Franklin Lakes, USA)). DNA extraction was performed using the **GenoLyse**[®] kit. Subsequently, the culture samples were tested with the **GenoType MTBDRsl** applying the "MDR CUL" PCR protocol.

All *M. tuberculosis* complex samples showed the same correct results. Hence, it can be excluded that the tested media import inhibitors into the **GenoType MTBDRsl**.

Interfering substances may also be carried over from the sample material. Hence, the substances indicated in table 10 were tested in order to assess a potential interference with the **GenoType MTBDRsl**. Defined BCG culture dilutions above, at, and below the detection limit of clinical samples were spiked with various amounts of the potential inhibitors. From all samples, DNA extraction was performed using the **GenoLyse**[®] kit. Then the culture dilutions were tested with the **GenoType MTBDRsl** applying the "MDR DIR" protocol for PCR.

Table 10: Tested potential interferents of the **GenoType MTBDRsl** VER 2.0

Substance/class	Description/active ingredient	Substance concentrations
Blood	Whole blood	2.5% v/v to 90% v/v
Blood	Hemoglobin	0.05% v/v to 13.5% v/v
Pus		0.5% v/v to 90% v/v

Interference of the **GenoType MTBDRsl** VER 2.0 (invalid test result) was observed in samples containing concentrations greater than 10% whole blood, 1% hemoglobin, and 2.5% pus.

Stability

Shelf life of the test kit when stored as recommended: see box label.

Stability is determined according to DIN EN ISO 23640.

References

- World Health Organization. Global tuberculosis report 2015. WHO/HTM/TB/2015.22. World Health Organization, Geneva, Switzerland 2015.
- Zhang Y, Yew WW. Mechanisms of drug resistance in *Mycobacterium tuberculosis*. *Int J Tuberc Lung Dis* 2009; 13: 1320-1330.
- Biosafety in microbiological and biomedical laboratories, 5th edition. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 2009.
- Protection of laboratory workers from occupationally acquired infections. Approved guideline. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards), USA, Document M29 (please refer to the latest version).
- Kent PT, Kubica GP. Public health mycobacteriology: a guide for the level III laboratory. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 1985.
- Isenberg HD. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C., USA 1992.
- Richter E, Beer J, Diel R, Hillemann D, Hoffmann H, Klotz M, Mauch H, Rüscher-Gerdes S. MiQ 5, Tuberkulose, Mykobakteriose. In: Podbielski A, Herrmann M, Kniehl E, Mauch H, Rüssmann H (eds): Mikrobiologisch-infektiologische Qualitätsstandards. Elsevier, Munich, Germany 2010.
- DIN, Deutsches Institut für Normung e.V. (ed). DIN 58943-4:2009-02: Medical microbiology - Diagnosis of tuberculosis - Part 4: Primary samples for the diagnosis of tuberculosis and mycobacteria – Qualitative and quantitative requirements, extraction, transport and storage. Beuth, Berlin, Germany 2009.
- Cheng AF, Yew WW, Chan EW, Chin ML, Hui MM, Chan RC. Multiplex PCR amplicon conformation analysis for rapid detection of *gyrA* mutations in fluoroquinolone-resistant *Mycobacterium tuberculosis* clinical isolates. *Antimicrob Agents Chemother* 2004; 48: 596-601.
- Aubry A, Veziris N, Cambau E, Truffot-Pernot C, Jarlier V, Fisher LM. Novel gyrase mutations in quinolone-resistant and -hypersusceptible clinical isolates of *Mycobacterium tuberculosis*: functional analysis of mutant enzymes. *Antimicrob Agents Chemother* 2006; 50: 104-112.
- Matrat S, Veziris N, Mayer C, Jarlier V, Truffot-Pernot C, Camuset J, Bouvet E, Cambau E, Aubry A. Functional analysis of DNA gyrase mutant enzymes carrying mutations at position 88 in the A subunit found in clinical strains of *Mycobacterium tuberculosis* resistant to fluoroquinolones. *Antimicrob Agents Chemother* 2006; 50: 4170-4173.
- Kocagöz T, Hackbarth CJ, Unsal I, Rosenberg EY, Nikaido H, Chambers HF. Gyrase mutations in laboratory-selected, fluoroquinolone-resistant mutants of *Mycobacterium tuberculosis* H37Ra. *Antimicrob Agents Chemother* 1996; 40: 1768-1774.
- Malik S, Willby M, Sikes D, Tsodikov OV, Posey JE. New insights into fluoroquinolone resistance in *Mycobacterium tuberculosis*: functional genetic analysis of *gyrA* and *gyrB* mutations. *PLoS One* 2012; 7: e39754. doi: 10.1371/journal.pone.0039754.
- Cole ST, Brosch R, Parkhill J, Garnier T, Churcher C, Harris D, Gordon SV, Eiglmeier K, Gas S, Barry CE 3rd, Tekaia F, Badcock K, Basham D, Brown D, Chillingworth T, Connor R, Davies R, Devlin K, Feltwell T, Gentles S, Hamlin N, Holroyd S, Hornsby T, Jagels K, Krogh A, McLean J, Moule S, Murphy L, Oliver K, Osborne J, Quail MA, Rajandream MA, Rogers J, Rutter S, Seeger K, Skelton J, Squares R, Squares S, Sulston JE, Taylor K, Whitehead S, Barrell BG. Deciphering the biology of *Mycobacterium tuberculosis* from the complete genome sequence. *Nature* 1998; 393: 537-44. Erratum in: *Nature* 1998; 396: 190.
- Maus CE, Plikaytis BB, Shinnick TM. Mutation of *tlyA* confers capreomycin resistance in *Mycobacterium tuberculosis*. *Antimicrob Agents Chemother* 2005; 49: 571-577.
- Maus CE, Plikaytis BB, Shinnick TM. Molecular analysis of cross-resistance to capreomycin, kanamycin, amikacin, and viomycin in *Mycobacterium tuberculosis*. *Antimicrob Agents Chemother* 2005; 49: 3192-3197.
- Georghiou SB, Magana M, Garfein RS, Catanzaro DG, Catanzaro A, Rodwell TC. Evaluation of genetic mutations associated with *Mycobacterium tuberculosis* resistance to amikacin, kanamycin and capreomycin: a systematic review. *PLoS One* 2012; 7: e33275. doi: 10.1371/journal.pone.0033275.
- Campbell PJ, Morlock GP, Sikes RD, Dalton TL, Metchock B, Starks AM, Hooks DP, Cowan LS, Plikaytis BB, Posey JE. Molecular detection of mutations associated with first- and second-line drug resistance compared with conventional drug susceptibility testing of *Mycobacterium tuberculosis*. *Antimicrob Agents Chemother* 2011; 55: 2032-2041.
- Chakravorty S, Lee JS, Cho EJ, Roh SS, Smith LE, Lee J, Kim CT, Via LE, Cho SN, Barry CE 3rd, Alland D. Genotypic susceptibility testing of *Mycobacterium tuberculosis* isolates for amikacin and kanamycin resistance by use of a rapid sloppy molecular beacon-based assay identifies more cases of low-level drug resistance than phenotypic Lowenstein-Jensen testing. *J Clin Microbiol* 2015; 53: 43-51.
- Zaunbrecher MA, Sikes RD Jr, Metchock B, Shinnick TM, Posey JE. Overexpression of the chromosomally encoded aminoglycoside acetyltransferase *eis* confers kanamycin resistance in *Mycobacterium tuberculosis*. *Proc Natl Acad Sci U S A* 2009; 106: 20004-20009.
- Tagliani E, Cabibbe AM, Miotto P, Borroni E, Toro JC, Mansjö M, Hoffner S, Hillemann D, Zalutskaya A, Skrahina A, Cirillo DM. Diagnostic performance of the new version (v2.0) of GenoType MTBDRsl assay for detection of resistance to fluoroquinolones and second-line injectable drugs: a multicenter study. *J Clin Microbiol* 2015; 53: 2961-2969.

Important Changes in IFU-317A-04

Chapter	Change
Reagents and Instruments	Three lot labels are included in the kit. New: "Conjugate Concentrate (CON-C) and Conjugate Buffer (CON-D) contain biological material. Hence, they must be considered as potentially infectious and must be handled accordingly."
Quality Control	New: "You can find the kit lot and the corresponding lots of the kit constituents on the lot labels included in the kit."
Limitations	New: "The test only works within the limits of the genomic regions the primers and probes were chosen from." New: "Please note that effects due to multiple mutations outside the investigated sequences cannot be detected by this test." New: "As any detection system based on hybridization, the test system on hand bears the possibility that sequence variations in the genomic regions the primers and probes were chosen from but the detection of which the test was not designed for may lead to false results. Due to the high variability of bacterial genomes, it is possible that certain subtypes might not be detected."

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Digitally signed by Gestion Documental Electronica
Date: 2022.02.11 14:53:52 -03:00

Digitally signed by Gestion Documental
Electronica
Date: 2022.02.11 14:53:52 -03:00

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NOMBRE Y DIRECCIÓN DEL FABRICANTE: Hain Lifescience GmbH, Hardwiesenstr 1. 72147 Nehren. Alemania.