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
ENGLISH

PALL Purecell PL HIGH EFFICIENCY LEUCOCYTE REMOVAL FILTER WITH ATTACHED SET FOR PLATELET TRANSFUSION

(Reorder No: PL3VAE)

INDICATIONS

The PALL Purecell PL High Efficiency* Leucocyte Removal Filter is indicated for the removal of leucocytes, microaggregates and the reduction of anaphylatoxin C3a† from 4 to 10 units of platelet concentrate or an equivalent single donor platelet collection.

Sterilised by irradiation.

STERILE and non-pyrogenic fluid pathway.

PRECAUTIONS

Follow instructions carefully.

For single use; do not resterilise or reuse. Do not use if packaging is damaged or end protectors are loose or displaced.

INSTRUCTIONS FOR USE

(Read and follow carefully)

1. PREPARATION OF SET AND PLATELET BAG

Agitate platelet bag and hang on drip stand.

Ensure cap (A) on vented spike is closed.

Ensure filter clamp (B) is closed and tethered protective cap (C) is hanging free.

Locate set clamp (D) directly downstream of drip chamber and close completely.

Insert filter set spike into platelet bag using a twisting motion.

2. TO PRIME FILTER AND DRIP CHAMBER

Ensure set is vertical and open filter clamp (B).

Allow filter and drip chamber to prime under gravity.

Flow will automatically stop when platelets reach air vent in drip chamber.

Once platelet flow into drip chamber has stopped completely, place tethered protective cap (C) on air vent and tighten.

3. TO TRANSFUSE platelets

Prime and connect tubing in usual manner.

Regulate flow with set clamp (D).

(20 Drops = Approximately 1 mL).

NOTE: At the end of the first unit, close filter clamp (B). Do not allow bag to empty and air to enter filter set.

4. WHEN PLATELET CONCENTRATES ARE GIVEN SEQUENTIALLY

Close set clamp (D). Agitate new platelet bag and hang on drip stand.

Remove spike from empty bag and insert into next unit using a twisting motion.**

Open filter clamp (B).

Grasp filter between thumb and forefinger of both hands and compress filter, hold for several seconds, then release. (This expels air, which may have been entrapped in tubing, into platelet bag).

Adjust flow with set clamp (D) to resume transfusion.

5. TO RECOVER PLATELETS

On completion of filtration open cap (A) on vented spike and allow upstream side of filter to drain. Filter is drained when flow into drip chamber stops.

NOTE:

1. The drip chamber should never be squeezed.

2. Saline priming is not required.

3. The filter should not be flushed with saline after filtration.

4. The cap on drip chamber must not be loosened or removed during transfusion.

5. This product is free of natural rubber latex.


* See product literature for performance summary.

** Always dispose of blood-contaminated products in a manner consistent with established biohazard procedures.

†C3a is a product of complement activation.

PALL, Pall and Purecell are trade marks of Pall Corporation.

Better Blood. Better Care. is a service mark of Pall Corporation


Filtration Separation Solutions
JORGE TROJAN
VICE-PRESIDENTE
PALLTECHNOLOGIES S.A.

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
ENGLISH

PALL Purecell PL HIGH EFFICIENCY LEUCOCYTE REMOVAL FILTER FOR PLATELETS (Re-order No: PL1BE)

INDICATIONS

The Pall Purecell PL High Efficiency* Leucocyte Removal Filter is indicated for the removal of leucocytes, microaggregates and the reduction of anaphylatoxin C3a† from a single unit of platelet concentrate prepared from one unit of whole blood or an equivalent volume of single donor platelet collection.

A 175 mL PVC administration bag with needleless access port is attached.

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PRECAUTIONS

Follow instructions carefully. For single use; do not resterilise or reuse. Do not use if packaging is damaged or end protectors are loose or displaced.

Filtered platelets should be transfused in accordance with standards and/or standard practices.

INSTRUCTIONS FOR USE

(Read and follow carefully)

Close clamp (A). Ensure cap of needleless access port (B) is secure.

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Remove protective cap from filter set spike and insert into platelet bag using a twisting motion.

Open clamp (A) to prime filter.

NOTE: Ensure filter hangs vertically below platelet bag during priming.

When filtration is complete, close clamp (A).

Hang administration bag and administer filtered platelets in usual manner.

NOTES:

1. To maintain traceability, do not remove platelet bag from filter set.
2. A needleless access port (B) is provided for syringe access, if required.
3. If an administration set is used for transfusion, any excess air should be removed from the bag by syringe or other appropriate means before transfusion.
4. This product is free of natural rubber latex.

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

† C3a is a product of Complement activation.

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Better Blood. Better Care. is a service mark of Pall Corporation.

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


Reorder No. PL6T

Purecell™ PL6T

High Efficiency Leukocyte Removal Filter for Platelets

ENGLISH - The Pall Purecell™ PL6T High Efficiency Leukocyte Removal Filter removes microaggregates, leukocytes and the anaphylatoxin C3a from platelet concentrates and single donor platelets. Designed for transfusion of 3 to 6 units of platelet concentrates or a single donor platelet collection. **STERILE**, non-pyrogenic, fluid pathway.

Indications: The Pall Purecell™ PL6T is a PL6T High Efficiency Leukocyte Removal Filter used to remove microaggregates, leukocytes and C3a from platelet concentrates.

Precautions: The Pall Purecell™ PL6T is a PL6T High Efficiency Leukocyte Removal Filter used to remove microaggregates, leukocytes and C3a from platelet concentrates. It is not intended for use with other blood components.

For more information, please contact Pall Corporation.

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Caution: Pall Corporation is not responsible for any adverse reactions to the use of this product.

Pall Corporation
2200 Northern Boulevard, East Hills, NY 11548 USA
USA 1-800-645-6578
<http://www.pall.com>

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
Reorder No. PL6KS

Purecell™ PL6KS

High Efficiency Leukocyte Removal Filter
with Attached Straight Administration Set
and Self-leveling Drip Chamber for Platelets

ENGLISH - The Pall Purecell™ PL6KS High Efficiency Leukocyte Removal Filter removes microaggregates, leucocytes and the anaphylatoxin C3a from platelet concentrates and single donor platelets. Designed for transfusion of 1 to 5 units of platelet concentrates or a single donor platelet collection. **STERILE** non-pyrogenic, latex free.

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Caution: Do not use this filter for the filtration of other types of blood components.

Pall Corporation
2200 Northern Boulevard, East Hills, NY 11548 USA
USA 1-800-645-6578
<http://www.pall.com>

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Medical

Purecell PL

MANUFACTURER
Pall Medical,
a division of Pall Europe Ltd.,
Europa House,
Havant Street,
Portsmouth,
PO1 3PD,
England

INTERNATIONAL OFFICES
Pall France S.A.
4 Rue des Chantiers,
BP 9290,
78170 St Germain Laval en France,
France

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CORPORATE HEADQUARTERS
Pall Corporation
2000 Northern Boulevard,
East Hill, NY 11742
USA
Visit us at the Web site
www.pall.com
14191 JVAE1

REF PL3VAE

Better Blood. Better Care.™

ENGLISH - PALL Purecell PL HIGH EFFICIENCY LEUCOCYTE REMOVAL FILTER WITH ATTACHED SET FOR PLATELET TRANSFUSION

STERILE and non-pyrogenic. Latex free. The product is not intended for distribution outside the USA.

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Pall Technologies S.A.
Nombre del fabricante: MEDSEP subsidiaria de PALL Corp.