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SD Biosensor, Inc.

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74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

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Digitally signed by GESTION DOCUMENTAL  
ELECTRONICA - GDE  
Date: 2023.11.08 22:36:35 -03:00

# STANDARD F TnI/CK-MB Combo FIA

STANDARD™ F TnI/CK-MB Combo FIA  
PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

## KIT CONTENTS

## MATERIALS REQUIRED BUT NOT PROVIDED


- STANDARD F Analyzer

## SPECIMEN COLLECTION AND STORAGE


- Whole blood (Venous whole blood)**
  - Collect the venous whole blood into the commercially available EDTA tube by venipuncture.
  - It is recommended that collected venous whole blood specimens are used immediately. If venous whole blood in an anticoagulant tube is stored at room temperature and refrigerated control (2-8°C/36-46°F), the specimen can be used for testing within 8 hours after collection.
  - Do not use hemolyzed blood specimens.
- Serum**
  - Collect the whole blood into the commercially available plain tube, NOT containing anti-coagulants such as heparin, EDTA or sodium citrate by venipuncture and leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
  - Serum specimen may be stored at room temperature and refrigerator conditions (2-8°C/36-46°F) for up to 24 hours prior to testing.
  - It should be brought to room temperature prior to use.

## PREPARATION AND TEST PROCEDURE

- Preparation**
  - Allow kit components and collected specimen to room temperature at least 30 minutes before starting the test.
  - Carefully read instructions for the STANDARD F TnI/CK-MB Combo FIA.
  - Check the expiry date at the back of the foil pouch. Use another lot if expiry date is expired.
  - Open the foil pouch and check the test device in the foil pouch.

The violet-line on the membrane of unused test device will disappear after use.

Do not write on the barcode or damage the barcode of the test device.

## Analysis of specimen

STANDARD TEST mode STANDARD F200 and F2400 Analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.
- |                         |   |
|-------------------------|---|
| STANDARD F2400 Analyzer | *Workplace* → *Run Test* → Scan or type patient ID and/or operator ID         |
| STANDARD F200 Analyzer  | *Standard Test* mode → Insert patient ID and / or operator ID on the analyzer |

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

- Using a disposable dropper, collect 100µl of the specimen (whole blood, serum) in the black line of the disposable dropper. Add the collected specimen into the extraction buffer and mix thoroughly.

- Apply 100µl of specimen mixture to the specimen well of the test device.

- After applying the specimen, immediately press the "TEST START" button.

- The analyzer will automatically display the test result within 10 minutes.

## PERFORMANCE CHARACTERISTICS

**1) Analytical performance**

**Accuracy (Method Comparison)**  
Results comparing the STANDARD F TnI/CK-MB Combo FIA with an automated chemiluminometric method are shown below.

	Regression Analysis			
	TnI	CK-MB	Whole Blood	CK-MB
Slope	0.9558	1.0016	1.0039	0.9914
Y-intercept	105.9	-0.0267	0.3879	0.2689
R	0.9954	0.9965	0.9961	0.9969
R <sup>2</sup>	0.9909	0.9931	0.9922	0.9939
n	210	210	210	210

**System Accuracy**

	Serum				Whole Blood			
	TnI	CK-MB	TnI	CK-MB	TnI	CK-MB	TnI	CK-MB
Below -25%	2 / 210 (1%)	0 / 210 (0%)	3 / 210 (1.4%)	0 / 210 (0%)				
Within ±25%	206 / 210 (98%)	210 / 210 (100%)	207 / 210 (98.6%)	210 / 210 (100%)				
Over +25%	2 / 210 (1%)	0 / 210 (0%)	0 / 210 (0%)	0 / 210 (0%)				

## INTERPRETATION OF TEST RESULTS

- Troponin I**
  - Measuring Range** 10 - 20,000 ng/L (0.01 - 20 ng/mL)
  - Out of measuring range**
    - If below 10 ng/L, "↓ 10 ng/L" message will be displayed.
    - If above 20,000 ng/L, "↑ 20,000" message will be displayed.
  - CV ≤ 10% at 99th Percentile URL** 70.0 ng/L (0.07 ng/mL)

Population	N	TnI Concentration (ng/L)		Maximum Observed
		99th Percentile URL	90% Confidence Interval	
732	732	70.0	54.0 - 81.0	100.0

Troponin I concentrations were determined using serum specimens obtained from 732 apparently healthy individuals. The upper reference limit (99th percentile) for TnI was 70 ng/L, and the lowest TnI concentration with CV 10% was found to be 70 ng/L.

Elevated troponin I is specific to myocardial injury and may be due to conditions other than MI. Acute MI can be established if there is a rise and/or fall in cardiac troponin with at least one value above the 99th percentile URL, together with a high clinical likelihood or ECG evidence of myocardial ischemia.

The default results are measured in ng/L. Should an alternate result unit be desired, the user can change unit through the "Setting" menu of the analyzer or manually calculate using the following conversion factors:

(Concentration in ng/L) x (0.001) = ng/mL = µg/L  
(Concentration in ng/L) x (1) = pg/mL

## CK-MB

- Measuring Range** 1-200 ng/mL
  - Out of measuring range**
    - If below 1 ng/mL, "1 ↓" message will be displayed.
    - If above 200 ng/mL, "200 ↑" message will be displayed.
  - Out-of-value for MI** <5.0 ng/mL
- People with skeletal muscle damage may have elevated CK-MB levels. To help differentiate heart attack from skeletal muscle damage, a CK index can be calculated using CK-MB and total CK (as follows: CK index = (CK-MB, ng/mL) x 100 / (total CK activity, IU/L). If CK-MB is elevated and the CK index is higher than 2.5 to 3.0, heart damage is likely.

The Troponin I and CK-MB reference ranges are provided for orientational purpose only. Clinicians should use the test results in conjunction with the patient's other diagnostic findings and clinical signs, and interpret the concrete values in the context of the patient's clinical situation.

## EXPLANATION AND SUMMARY

### Introduction

The Troponin (Tn), one of the subunits of the troponin regulatory complex, binds to actin and inhibits interactions between actin and myosin. Troponin complex plays an important role in the regulation of skeletal and cardiac muscle contraction. Troponin complex consists of three subunits: Troponin T (TnT), Troponin I (TnI) and Troponin C (TnC). The subunits are held together by non-covalent interactions. TnI is the most important subunit that binds to actin and inhibits interaction with myosin and the thin filament. The TnI subunit is responsible for inhibiting actomyosin formation of low intracellular Ca<sup>2+</sup> concentrations. The third subunit, TnC, binds Ca<sup>2+</sup> ions during the excitation of the muscle and changes the conformation of troponin complex, thus enabling the formation of actomyosin complex and the subsequent muscle contraction. Cardiac Troponin I (cTnI) has been found to be a reliable marker of cardiac cell death. It is now widely used for the diagnosis of acute myocardial infarction (AMI), postmyocardial trauma, chemotherapy cardiotoxicity, as well as some other diseases related to cardiac muscle injury. Being expressed cTnI only in cardiac tissue, cardiac TnI have been the preferred biomarkers for myocardial infarction for a long time. Creatine Kinase (CK), also known as Creatine Phosphokinase or Phospho-creatine Kinase is an enzyme expressed by various tissues and cell types. Disruption of cell membranes due to hypoxia or other injury releases CK from the cellular cytosol into the systemic circulation. CK enzyme is a dimer formed by the association of two subunits from muscle (Muscle type) or nerve cells (Brain type). These subunits associate to form three isoenzymic forms: CK-BB, CK-MM and CK-MB. These isoenzymes are expressed at different levels in various human tissues. Though CK-MM is the most abundant CK isoenzyme in the cardiac muscles, CK-MB constitutes about 20% of the total CK in the cardiac muscle tissue. Elevated levels of total CK is not specific to the myocardial tissue and may be observed in patients with skeletal muscle injury and certain other disorders but as CK-MB is more specific to cardiac muscle, CK-MB levels along with total CK can be considered as an important diagnostic indicator of myocardial infarction. The concentration of CK-MB in the healthy adult maintains low concentration but it shows great increases in several malignant diseases, mostly primary coronary syndrome, myocardial injury and infarction. CK-MB has been found to be more sensitive and early indicator of myocardial injury because it has a lower basal level and a much narrower normal range. Medical literature commonly reveals that following an acute myocardial infarction, CK-MB levels become elevated in 4 to 9 hours after the onset of chest pain, attain peak at 10 to 24 hours, and return to normal within 2 to 3 days. Use of CK-MB level as a percentage of total CK in the diagnosis of myocardial infarction is the most important of CK measurements in clinical chemistry. STANDARD F TnI/CK-MB Combo FIA, employing immunofluorescent detection system with STANDARD F Analyzer, provides significantly fast, easy and accurate system to identify cTnI and CK-MB levels in human serum and whole blood as an aid in the diagnosis of myocardial infarction and enables supportive treatment decisions.

### Intended use

STANDARD F TnI/CK-MB Combo FIA is a fluorescent immunoassay for the quantitative determination of cardiac troponin I and total creatine kinase isoenzyme-MB(CK-MB) levels in human serum and whole blood using STANDARD F Analyzers manufactured by SD BIOSENSORS. This test is an *in vitro* professional diagnostic use and intended for use as an aid in the screening and monitoring of myocardial infarction(MI).

### Test Principle

STANDARD F TnI/CK-MB Combo FIA uses antibodies that specifically bind to the TnI and CK-MB proteins. When a human specimen is applied into the specimen well of the test device, the specimen migrates through the membrane from the specimen well. If Troponin I (TnI) and/or CK-MB is present, it will be bound to the detector coupled to europium microparticle that migrates through the membrane. The europium microparticle complex will be captured by capture antibody on the test line where it is detected by STANDARD F Analyzer. If it is not present, the europium microparticle will not be trapped by the capture antibody nor detected by STANDARD F Analyzer. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer. STANDARD F Analyzer can analyze the presence of TnI and CK-MB simultaneously in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

### CK STORAGE AND STABILITY

For 18 months when stored at 2-30°C/36-86°F and 10-90% RH.

### WARNINGS AND PRECAUTIONS

- STANDARD F TnI/CK-MB Combo FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this Instructions for use before testing.
- STANDARD F TnI/CK-MB Combo FIA should be used with STANDARD F Analyzer.
- STANDARD F TnI/CK-MB Combo FIA should remain in its original sealed pouch until ready to use. Do not use the test device if pouch is damaged or the seal is broken.
- STANDARD F TnI/CK-MB Combo FIA is single use only. Do not re-use it.
- Do not use hemolyzed specimens or frozen specimens.
- Do not use any artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm and soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Desiccant in foil pouch should absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

- Check the expiration date printed at the pouch or package.
- Check the volume (100µl) of extraction buffer.
- Use the STANDARD F TnI/CK-MB Combo FIA at room temperature (15-32°C/59-89.6°F).
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

## LIMITATION OF TEST

- The test should be used for the detection of Troponin I and CK-MB levels in human whole blood and serum specimens.
- Failure to follow the test procedure and interpretation of test results may adversely affect the results.
- Invalid results may occur if a poor-quality specimen is obtained.
- The test result must always be evaluated with other clinical data available to the physician.

## QUALITY CONTROL

- Calibration**  
The calibration set test of STANDARD F Analyzer should be conducted according to the analyzer's manual.

### When to use calibration set

- Before using the analyzer for the first time.
- When you drop the analyzer.
- Whenever you do not agree with the final result.
- When you want to check the performance of the analyzer and test device.

### How to use calibration set

- Calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.
- Select the "Calibrator" in main menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 for white calibration, CAL-2 for UV LED calibration, and CAL-3 for RGB LED calibration in sequence.

STANDARD F Analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in STANDARD TEST™ mode. If "EE" message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the "EE" message still appears.

### Internal quality control

- The internal procedural control zone is on the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.
- The invalid result denotes that the fluorescence signal is not within the pre-set range. If the screen of STANDARD F Analyzer shows "Invalid", turn off and turn on the analyzer again and re-test with a new test device.

### External quality control

Quality control testing should be run to check the performance of STANDARD F TnI/CK-MB Combo FIA and STANDARD F Analyzer. STANDARD F TnI Control and STANDARD F CK-MB Control manufactured by SD BIOSENSORS can be used for the external quality control. Control test should be conducted in accordance with the Instructions for use of STANDARD F TnI Control and STANDARD F CK-MB Control. External quality control test should be run:

- once for each new lot.
- once for each trained operator.
- as required by test procedures in the Instructions for use of STANDARD F TnI Control and STANDARD F CK-MB Control in accordance with local, state and federal regulations or accreditation requirements.

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

## CONTENIDO DEL KIT

The violet-line on the membrane of unused test device will disappear after use.

Do not write on the barcode or damage the barcode of the test device.

## Analysis of specimen

STANDARD TEST mode STANDARD F200 and F2400 Analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.
- |                         |   |
|-------------------------|---|
| STANDARD F2400 Analyzer | *Workplace* → *Run Test* → Scan or type patient ID and/or operator ID         |
| STANDARD F200 Analyzer  | *Standard Test* mode → Insert patient ID and / or operator ID on the analyzer |

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer read the barcode data, and check the test device is valid.

- Using a disposable dropper, collect 100µl of the specimen (whole blood, serum) in the black line of the disposable dropper. Add the collected specimen into the extraction buffer and mix thoroughly.

- Apply 100µl of specimen mixture to the specimen well of the test device.

- After applying the specimen, immediately press the "TEST START" button.

- The analyzer will automatically display the test result within 10 minutes.

## PERFORMANCE CHARACTERISTICS

**1) Analytical performance**

**Accuracy (Method Comparison)**  
Results comparing the STANDARD F TnI/CK-MB Combo FIA with an automated chemiluminometric method are shown below.


	Regression Analysis			
	TnI	CK-MB	Whole Blood	CK-MB
Slope	0.9558	1.0016	1.0039	0.9914
Y-intercept	105.9	-0.0267	0.3879	0.2689
R	0.9954	0.9965	0.9961	0.9969
R <sup>2</sup>	0.9909	0.9931	0.9922	0.9939
n	210	210	210	210

**System Accuracy**

	Serum				Whole Blood			
	TnI	CK-MB	TnI	CK-MB	TnI	CK-MB	TnI	CK-MB
Below -25%	2 / 210 (1%)	0 / 210 (0%)	3 / 210 (1.4%)	0 / 210 (0%)				
Within ±25%	206 / 210 (98%)	210 / 210 (100%)	207 / 210 (98.6%)	210 / 210 (100%)				
Over +25%	2 / 210 (1%)	0 / 210 (0%)	0 / 210 (0%)	0 / 210 (0%)				

The violet-line on the membrane of the device will disappear after use.

Do not write on the barcode or damage the barcode of the test device.

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