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The CD4 easy count kit is a manual, two component, quantitative IVD test for subpopulation labelling of lymphocytes in adult venous EDTA whole blood, and subsequent enumeration with the Sysmex Partec CyFlow™ Counter IVD flow cytometer. The CD4 T-cell concentration of blood samples is a useful indicator for the initiation or follow-up of treatment for HIV positive patens in conjunction with other laboratory and clinical findings.

The test is intended to be performed by trained healthcare professionals.

### 2 Intended Use

**IND** For In Vitro Diagnostic Use. The CD4 easy count kit is a manual, two component, quantitative IVD test for subpopulation labelling of lymphocytes in adult venous EDTA whole blood, and subsequent enumeration with the Sysmex Partec CyFlow™ Counter IVD flow cytometer. The CD4 T-cell concentration of blood samples is a useful indicator for the initiation or follow-up of treatment for HIV positive patens in conjunction with other laboratory and clinical findings.

The test is intended to be performed by trained healthcare professionals.

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An aliquot of an EDTA whole-blood sample is mixed with the antibody (CD4) conjugated to the fluorochrome in a 1:1 ratio. After a fixed incubation time, the buffer is added and the sample is ready for analysis on the CyFlow™ Counter flow cytometer. The light source excites the fluorescent dye linked with the stained cell and the emitted light is detected while a certain volume of blood sample is running through the instrument. The integrated software calculates the concentration of the dedicated cell populations. For further information, please refer to the IFU of the CyFlow™ Counter.

### 4 Storage and shelf life after first opening

1. Storage Store the antibody and buffer reagents at 2 - 8 °C in the dark. Do not freeze or expose the reagents to elevated temperatures and keep them away from direct sunlight.

2. Shelf life after first opening: Always close the bottle after use and use new pipette tips each time the reagent is sampled to avoid contamination. Keep the antibody and buffer reagent at 2 - 8 °C in the dark. Do not freeze or expose the reagents to elevated temperatures and keep them away from direct sunlight.

Under these storage conditions, the CD4 easy count kit will be stable until the expiration date printed on the kit label.

### 5 Components

CD4 mAb PE is a murine monoclonal antibody supplied in PBS buffer with 0.2% BSA and 0.09% sodium azide. No lyse buffer is a PBS-based solution containing 0.09% sodium azide.

### 6 Evidence of deterioration

The antibody and buffer solutions are clear liquids. Do not use the reagents if there is the appearance of any kind of turbidity or contamination. For questions regarding the performance or quality of the product received, please contact your local Sysmex representative.

### 7 Precaution and warnings

Reagents contain 14 mM sodium azide as a preservative. The low concentration of sodium azide does not require hazard labelling, but the normal safety precautions for the handling of chemicals must be observed. Refer to the Safety Data Sheet (SDS) for a complete list of warnings and precautions.

### 8 Additional required equipment

Instrument required:	CyFlow™ Counter (REF No.: CY-S-3022 or CY-S-3023)
Required apparatus:	Sample Tubes 3.5 ml (REF No.: 04-2000) <p>A vented pipette 20 µl and pipette tips A vented pipette 100 – 1000 µl and pipette tips Venous blood collection system with EDTA as anticoagulant Stop watch</p>

### 9 Reagent preparation

CD4 mAb PE (REF No.: 05-8401-01), the antibody reagent is ready to use.

No lyse buffer (REF No.: 05-8401-02), the buffer solution is ready to use.

### 10 Primary sample collection, handling and storage

**WARNING** ⚠ *All biological specimens and material are considered as biohazards and should be handled as if capable of transmitting infection. Appropriate safety precautions and handling procedures must be applied in accordance with applicable laws and regulations.*

- The blood sample must be a venous, whole-blood sample, collected with a blood collection system which contains EDTA as anticoagulant.
- The blood sample must be transported in a dark container, protected from light and not exposed to elevated temperatures. A suitable transport container is, for example, a Styrofoam box with cooling material. Best transport conditions would be below 25 °C if transported and not take more than 6 hours after donation. For transport and storage time more than 6 hours, blood should be kept at 2 - 8 °C.
- Do not freeze and thaw the blood sample.
- For best conditions, the blood sample should be fresh, i.e. not more than six hours should have elapsed between sample collection and analysis.
- Blood samples can be used up to 24 hours after drawing if they are stored in the fridge at 2 – 8 °C.
- Before use, the sample tubes containing the blood must be inverted gently 8 to 10 times.
- Do not use clogged blood samples.

### 11 Examination procedure

**NOTICE** *Reverse pipetting is critical to accuracy especially when dispensing viscous and very small sample volumes. For whole blood pipetting, we recommend using a calibrated electronic pipette which is preprogrammed to operate in the reverse pipetting mode. If an electronic pipette is not available, follow these instructions for manual reverse pipetting:*

- Depress the operating button to the second stop. Let the operating button move up completely, excess sample is drawn up into the tip.*
- Depress the operating button to the first stop to expel a precise volume of blood, leaving excess sample in the tip.*
- Discard pipette tip with excess blood sample in the tip into medical waste container. Please note, you must not use reverse pipetting for antibody or buffer solutions.*

### 11.1 Staining

- Invert the blood sample in the EDTA blood collection tube gently 8 to 10 times.
- Pipette 20 µl EDTA whole blood sample to the bottom of a sample tube, avoid leaving a trail of blood on the inner tube wall from pipette tip. Discard the pipette tip.
- Add 20 µl CD4 mAb PE directly into the blood sample and mix it gently, avoid leaving a trail of blood on the inner tube wall from pipette tip. Discard the pipette tip.
- The sample should not be vortexed at this stage of preparation.

- Incubate the mixture for 15 minutes at room temperature in the dark.
- Add 800 µl no lyse buffer and vortex briefly, or tap tube gently in order to mix the sample. Discard the pipette tip.

- Vortex briefly, or tap tube gently before analysing with the CyFlow™ Counter.

**NOTICE** *The stained blood sample must be analysed within 2 hours after adding the no lyse buffer.*

### 11.2 Sample analysis

Please refer to the IFU of the CyFlow™ Counter for how to operate the instrument before analysing the samples. The start-up procedures and the internal quality control procedure must be successfully completed before sample analysis.

- Choose and load the configuration for CD4 measurement from the „menu bar of CyFlow™ Counter.“
- Insert the sample tube with the prepared blood sample into the sample port.
- Start measurement.
- After sample run, the instrument stops and cleans automatically.
- Remove and dispose of blood sample in accordance with local laboratory biohazard safety procedures.

Do not start a second measurement with the same sample tube. If a duplicate or repeat measurement is required, a second sample tube should be prepared as per 11.1.

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Please refer to the IFU of the CyFlow™ Counter for proper start-up procedure and maintenance of the device, if necessary.

The CD4 T-cell counting result will be displayed on the screen. If the histogram does not show a clear peak of the CD4 positive stained T-cells, the sample preparation and measurement must be repeated.

Results will be incorrect if gate positioning is not done precisely. Be careful not to gate monocytes as CD4 positive stained cells if CD4 concentration is extremely low in a patient blood sample.

### 12 Calculation of examination results

Calculation of the cell concentration is part of the software-based data analysis, which provides the result in CD4 cells per µl blood sample each of each measurement.

For further information, please refer to the IFU of the CyFlow™ Counter.

### 13 Interpretation of results

Human Immunodeficiency Virus (HIV) is one of the main reasons for CD4 cell depletion, which debilitates the immune system. Results should be interpreted by physicians, in conjunction with the locally applicable HIV treatment guidelines.

### 14 Control procedure

Appropriate quality control should be performed according to local and national regulations. A variety of standardized blood samples are available on the market, only a few are suitable when using certain flow cytometers. If the acceptance criterion of the individual control material is not met do not proceed with patient specimen testing and contact your local Sysmex representative.

For further information about the use of adequate control material, please contact your local Sysmex representative.

### 15 Performance characteristics

#### 15.1 Specificity

The mouse monoclonal antibody MEM-241 recognizes the human CD4 antigen, a transmembrane glycoprotein (55 kDa) of the immunoglobulin supergene family, present on a subset of T-lymphocytes (T-helper/T-H cells) and expressed at a lower level on monocytes, tissue macrophages and granulocytes. HCDM (former H.LDA VIII) Meeting, May 2006, Québec, Canada; WS Code M241 Detection capability. The evaluation for Limit of Blank (LoB) and Limit of Detection (LoD) and Limit of Quantitation (LoQ) was carried out according to the definition in the guideline. CLSI EP17-A2.

LoB:   ≤50% of LoD
LoD:   ≥20 CD4 cells/µl
LoQ:   ≥40 CD4 cells/µl

#### 15.2 Repeatability

Intra- and inter-assay variation:

CV values are less than 10% for CD4 values of 200 CD4 T-cells/µl or above.

CV values are less than 15% for CD4 absolute values below 200 CD4 T-cells/µl.

#### 15.3 Linearity

The measuring range of the assay and linearity of the measurement are valid over the range:
CD4 absolute: 40-2500 CD4 T-cells/µl

#### 15.4 Trueness

EDTA whole blood specimen were stained with CD4 easy count kit and analysed on the CyFlow™ Counter. CD4 absolute values were compared with results from a BD FACSCalibur using BD Tristest CD3/CD4/CD45 and BD Trucount tubes and showed a bias of -15,13% between both methods.

### 16 Limitations

**CAUTION** ⚠ *Sysmex Partec CD4 easy count kits (05-8401) have not been validated for use with pediatric and adolescent patients.*

**CAUTION** ⚠ *Patient measurements of CD4 and CD4% made using Sysmex Partec assays should not be used independently with other manufacturer's methods of determining CD4 and CD4%. Use the assays only with the Sysmex Partec instruments. Values obtained using other equipment or assays may not be interchangeable.*

Different endogenous and exogenous substances were tested regarding the CLSI EP07-A2 guideline for possible interfering effects which could have an influence to the blood sample analysis. Finally, ± 18% haemolysis (induced by using the protocol recommended in the CLSI EP07-A2, Appendix G, G1 – Osmotic Shock Procedure) was observed to show an interfering effect. For that reason haemolyzed samples should be rejected.

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### 19 Version number of IFU and date of issue

Revision: Rev-008\_05-04-2022, CN 2452  
Issued by: Sysmex Partec GmbH

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- 1) M. Fryland, P. Chaillet, R. Zachariah, A. Barnaba, L. Bonte, R. Andreassens, S. Charrondière, R. Teck, O. Didikus. (2006). The Partec CyFlow Counter™ could provide an option for CD4+ T-cell monitoring in the context of scaling-up antiretroviral treatment at the district level in Malawi. Transactions of the Royal Society of Tropical Medicine and Hygiene, 100(10):980-5
- 2) L. Zijenah, G. Kadzirange, S. Madzime, M. Borok, C. Mudwa, O. Tobaiwa, M. Mucheche, S. Rusakaniko, D. Katzenstein. (2006). Affordable flow cytometry for enumeration of absolute CD4+ T-lymphocytes to identify subtype C HIV-1 infected adults requiring antiretroviral therapy (ART) and monitoring response to ART in a resource-limited setting. Journal of Translational Medicine 4:33
- 3) H. Karcher, D. Bohning, R. Downing, S. Mashate, G. Harms. (2006). Comparison of two alternative methods for CD4+ T-cell determination (Coulter manual CD4 count and CyFlow) against standard dual platform flow cytometry in Uganda. Clinical Cytometry, 70B (3):163-69
- 4) S. Kohrsenudom, S. Dettraira. (2012). Evaluation of the Partec Single-Platform Volumetric CyFlow™ Counter System for Determining Percentage and Absolute Numbers of CD4 T Lymphocytes in HIV/AIDS Thai Patients. Thai AIDS J, 25:1-10
- 5) D. Wade, P.A. Diaw, G. Daneau, M. Camara, T.N. Dyeie, S. Mboup, L. Kestens. (2013). CD4 T-Cell Enumeration in a Field Setting: Evaluation of CyFlow Counter Using the CD4 Easy Count Kit-Dry and Pima CD4 Systems. PLoS ONE, 8(9): e75484

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## Product Labels

**Product:** CD4 easy count kit  
**Code No.:** 05-8401

Product Label	Container Type	Label Size
 <p>The label for the CD4 easy count kit includes the Sysmex logo, product name, and identification details. It specifies the reference number (REF 05-8401) and lot number (LOT). The description states it is for the identification and enumeration of the CD4+ helper/inducer T-lymphocyte subset in human blood samples. Contents include one vial of CD4 mAb PE (05-8401-01) and one lyse buffer bottle (05-8401-02). Storage instructions are 2°C to 8°C. The label also features CE, IVD, and I icons, the website www.sysmex-partec.com/services, and the manufacturer information: Sysmex Partec GmbH, Arndtstraße 11a-b, 02826 Görlitz, Germany. A barcode is present with the sample EAN (01)04250878904382(17)201113(10)700103 and Lab-0109-2.</p>	cardboard box (white)	70 x 50 mm
 <p>The label for the CD4 mAb PE component includes the Sysmex logo, product name, and identification details. It specifies the reference number (REF 05-8401-01) and lot number (LOT). The description states it is for the identification and enumeration of the CD4+ helper/inducer T-lymphocyte subset in human blood samples. Storage instructions are 2°C to 8°C. The label also features CE, IVD, and I icons, the website www.sysmex-partec.com/services, and the manufacturer information: Sysmex Partec GmbH, Arndtstraße 11a-b, 02826 Görlitz, Germany. Lab-0110-2 is also noted.</p>	3 ml glass vial (amber glass) with white PP cap	50 x 19 mm
 <p>The label for the no lyse buffer component includes the Sysmex logo, product name, and identification details. It specifies the reference number (REF 05-8401-02) and lot number (LOT). The description states it is for the identification and enumeration of the CD4+ helper/inducer T-lymphocyte subset in human blood samples. Storage instructions are 2°C to 8°C. The label also features CE, IVD, and I icons, the website www.sysmex-partec.com/services, and the manufacturer information: Sysmex Partec GmbH, Arndtstraße 11a-b, 02826 Görlitz, Germany. Lab-0111-2 is also noted.</p>	100 ml glass bottle (amber glass) with white PP cap	60 x 30 mm

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Identification of the IVD reagent

Product name: CD4% easy count kit

Content:	REF 05-8405
Vial containing CD4 mAb PE	REF 05-8405-01
Vial containing CD45 mAb PE-Cy5	REF 05-8405-02
Bottle containing Buffer 1	REF 05-8405-03
Bottle containing Buffer 2	REF 05-8405-04

1 Spezifikation	Human CD4	Human CD45
Spezifität	IgG1	IgG1
Isotyp	100 tests	100 tests
Contentu	PE	PE-Cy5
Fluorochrome	532/488	532/488
Excitation (nm)	578	670
Emission maximum (nm)		

2 Intended Use [VD] For In Vitro Diagnostic Use.

The CD4% easy count kit is a manual, four-component, quantitative IVD test for labelling of leukocytes and a subpopulation of lymphocytes in an adult, venous EDTA whole blood, which can then be enumerated with the Sysmex Partec CyFlow™ Counter IVD flow cytometer. The CD4 T-cell concentration and CD4% of lymphocytes are useful indicators for the initiation or follow-up of treatment for HIV positive patients in conjunction with other laboratory and clinical findings. The test is intended to be performed by trained healthcare professionals.

### 3 Principle of the examination method

An aliquot of an EDTA whole-blood sample is mixed with two antibodies (CD4 and CD45), each conjugated to a different fluorochrome for labelling dedicated cell populations. After a fixed incubation time, the two buffer solutions are added and the sample is ready for analysis on the CyFlow™ Counter flow cytometer. The light source excites the fluorescent dye linked to the stained cell and the emitted light is detected while a precise volume of blood sample is running through the instrument. The concentration of the dedicated cell populations is calculated by the integrated software.

4 Storage and shelf life after first opening

1. Storage: Store the antibody and buffer reagents at 2 - 8 °C in the dark. Do not freeze or expose the reagents to elevated temperatures and keep them away from direct sunlight.

2. Use: Always close the bottle after use and use new pipette tips each time the reagent is sampled to avoid contamination. Under these storage conditions, the CD4% easy count kit will be stable until the expiration date printed on the label.

5 Components

CD4 mAb PE and CD45 mAb PE-Cy5 are murine monoclonal antibodies supplied in PBS buffer with 0.2% BSA and 0.0% sodium azide. Buffer 1 and Buffer 2 are PBS-based solutions containing 0.09% sodium azide.

6 Evidence of deterioration

The antibody and buffer solutions are clear liquids. Do not use the reagents if there is the appearance of any kind of turbidity or contamination. For questions regarding the performance or quality of the product, please contact your local Sysmex representative.

7 Precaution and warnings

Reagents contained in this kit are used as a preservative. The low concentration of sodium azide does not require hazard labelling, but the normal safety precautions for the handling of chemicals must be observed. Refer to the Safety Data Sheet (SDS) for a complete list of warnings and precautions.

8 Additional required equipment

Instrument requirement: CyFlow™ Counter (REF No.: CY-S-3022 or CY-S-3023)

Required apparatus: Sample Tubes 3.5 ml (REF No.: 04-2000)

A verified pipette 10 µl and pipette tips

A verified pipette 20 µl and pipette tips

A verified pipette 100 - 1000 µl and pipette tips

Venous blood collection system with EDTA as anticoagulant

Stop watch

9 Reagent preparation

CD4 mAb PE (REF No.: 05-8405-01), the antibody reagent is ready to use.

CD45 mAb PE-Cy5 (REF No.: 05-8405-02), the antibody reagent is ready to use.

Buffer 1 (REF No.: 05-8405-03), the buffer solution is ready to use.

Buffer 2 (REF No.: 05-8405-04), the buffer solution is ready to use.

10 Primary sample collection, handling and storage

WARNING: All biological specimens and material are considered as biohazards and should be handled as if capable of transmitting infection. Appropriate safety precautions and handling procedures must be applied in accordance with applicable laws and regulations.

The blood sample must be a venous, whole-blood sample, collected with a blood collection system which contains EDTA as anticoagulant.

The blood sample must be transported in a dark container, protected from light and not exposed to elevated temperatures in a suitable transport container is, for example, a Styrofoam box with cooling material. Best transport conditions would be below 25 °C; if transport will not take more than 6 hours after donation.

For transport and storage time more than 6 hours, blood must be kept at 2 - 8 °C.

Do not freeze and do not thaw the blood sample.

For best conditions, the blood sample should be fresh, i.e. not more than six hours should have elapsed between sample collection and analysis.

Blood samples can be used up to 24 hours after drawing if they are stored in the fridge at 2 - 8 °C.

Do not use the sample tubes containing the blood must be inverted gently 8 to 10 times.

Do not use clogged blood samples.

11 Examination procedure

Reverse pipetting is critical to accuracy especially when dispensing viscous and very small sample volumes. For whole blood pipetting, we recommend using a calibrated electronic pipette which is pre-programmed to operate in the reverse pipetting mode.

If an electronic pipette is not available, follow these instructions for manual reverse pipetting:

Depress the operating button to the second stop. Let the operating button move up completely, excess sample is drawn up into the tip.

Depress the operating button to the first stop to expel a precise volume of blood, leaving excess sample in the tip.

Discard pipette tip with excess blood sample in the tip into medical waste container.

Please note, you must not use reverse pipetting for antibody or buffer solutions.

11.1 Staining

Invert the blood sample in the EDTA blood collection tube gently 8 to 10 times.

Pipette 20 µl EDTA whole blood sample to the bottom of a 3.5 ml sample tube, avoid leaving a trail of blood on the inner tube wall from pipette tip. Discard the pipette tip.

Add 10 µl CD4 mAb PE directly into the blood sample and mix it gently, avoid leaving a trail of blood on the inner tube wall from pipette tip. Discard the pipette tip.

Add 10 µl CD45 mAb PE-Cy5 directly into the blood sample and mix it gently, avoid leaving a trail of blood on the inner tube wall from pipette tip. Discard the pipette tip.

Incubate the mixture for 15 minutes at 15 - 30 °C in the dark.

Add 400 µl Buffer 1 and vortex briefly, or tap tube gently in order to mix the sample. Discard the pipette tip.

Add 400 µl Buffer 2 and vortex briefly, or tap tube gently in order to mix the sample. Discard the pipette tip.

Vortex briefly, or tap tube gently before analysing with the CyFlow™ Counter.

NOTICE: After addition of Buffer 2, samples must be stored for up to 2 - 8 °C in the dark. Do not use the sample later than that or if not prepared accordingly.

11.2 Sample analysis

Please refer to the IFU of the CyFlow™ Counter for how to operate the instrument before analysing the sample. The start-up procedure and the internal quality control procedure must be successfully completed before sample analysis.

Choose and load the configuration for CD4% measurement from the menu bar of CyFlow™ Counter.

Insert the sample tube with the prepared blood sample into the sample port.

Start measurement.

After sample run, the instrument stops and cleans automatically.

Remove and dispose of blood sample in accordance with local laboratory biohazard safety procedures.

Do not start a second measurement with the same sample tube. If a duplicate or repeat measurement is required, a second sample tube has to be prepared as per 11.1.

11.3 Data acquisition and analysis

Data analysis with the CyFlow™ Counter is only possible if the internal quality control of the instrument was successful.

Please refer to the IFU of the CyFlow™ Counter for start-up procedure and maintenance of the device, if necessary.

The CD4 T-cell counting result, CD4 percentage and total count of lymphocytes will be displayed on the screen. For results to be reliable there must be good separation of cell clusters as shown in Figure 1. If the clusters of stained cells towards the debris signals in the CD45 - SSC dot plot are not clearly separable (see Figure

not clearly separable, the sample preparation and measurement must be repeated.


Figure 1: Example of well separated clusters of granulocytes. Other CD45 positive stained leukocytes and debris signals.


Figure 2: Example of bad separation between debris and CD45. Other CD45 positive stained leukocytes. This sample is not analysable.

Results will be incorrect if gate positioning of CD45 positive stained cells (R3) is not done precisely.

12 Calculation of examination results

Calculation of the cell concentration is part of the software-based data analysis, which provides the result in CD4%, and CD4% cells per µl blood sample after each measurement. For further information, please refer to the IFU of the CyFlow™ Counter.

13 Interpretation of results

Human Immunodeficiency Virus (HIV) is one of the main reasons for CD4 cell depletion, which debilitates the immune system. Results should be interpreted by physicians, in conjunction with locally applicable HIV treatment guidelines.

14 Control procedure

Appropriate quality control should be performed according to local and national regulations. A variety of standardized blood samples are available on the market, only a few are suitable, when using certain flow cytometer. If the acceptance criterion of the individual control material is not met, do not proceed with patient specimen testing and contact your local Sysmex representative. For further information about the use of adequate control material, please contact your local Sysmex representative.

15 Performance characteristics

15.1 Specificity

The antibody MEM-241 recognizes CD4 antigen, a 55 kDa transmembrane glycoprotein expressed on a subset of T lymphocytes (‘helper’ T-cells) and on monocytes, tissue macrophages and granulocytes. HCM (former HLA VIII) Meeting, May 2006, Québec, Canada; W5 Code M241

The antibody MEM-28 reacts with all alternative forms of human CD45 antigen (Leukocyte Common Antigen), a 180-220 kDa single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except erythrocytes and platelets. HLA III; W5 Code NL 833a

15.2 Detection ability

The antibody MEM-28 reacts with all alternative forms of human CD45 antigen (Leukocyte Common Antigen), a 180-220 kDa single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except erythrocytes and platelets. HLA III; W5 Code NL 833a

15.3 Repeatability

Intra- and inter-assay variation: CV values are less than 10% for CD4 values and 20% CD4 T-cells/µl or above. CV values are less than 10% for CD4% values for 25% or above.

CV values are less than 15% for CD4 absolute values below 200 CD4 T-cells/µl. SD values are less or equal to 2.5% for CD4% values below 25%.

15.4 Linearity

The measuring range of the assay and linearity of the measurement are valid over the range: CD4 absolute: 40-2500 CD4 T-cells/µl

CD4%: 4 - 60%

15.5 Trueness

EDTA whole blood specimen were stained with CD4% easy count kit and analysed on the CyFlow™ Counter (CyView™ 2.11). CD4 absolute and CD4% values were compared with results from a BD FACSCalibur using BD Tristest CD3/CD4/CD45 and BD Trucount tubes and showed a bias of -22.56% (CD4 absolute) and +2.1% (CD4%) between both methods.

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CAUTION: Patient measurements of CD4 and CD4% made using Sysmex Partec assays should not be used alongside measurements obtained from other manufacturers' methods of determining CD4 and CD4%. Use the assays only with Sysmex Partec instruments. Values obtained using other equipment or assays are not interchangeable.

Different endogenous and exogenous substances were tested regarding the CLSI EP07-A2 guideline for possible interfering effects which could have an influence to the blood sample analysis. Finally, ≥ 18% haemolysis (induced by using the protocol recommended in the CLSI EP07-A2, Appendix G, G1 - Osmotic Shock Procedure) was observed to show an interfering effect. For that reason, haemolyzed samples should be rejected. For CD4 and CD4% determination in whole blood samples with CyFlow™ Counter assays, no other anticoagulant than EDTA should be used.

For further questions regarding the performance of the product, please contact your local Sysmex representative.

17 Literature references

1) M. Fryland, P. Chaillet, R. Zachariah, A. Barnaba, L. Bonte, R. Andreassen, S. Charrondière, R. Teck, O. Didakus. (2006). The Partec CyFlow Counter could provide an option for CD4+ T-cell monitoring in the context of scaling-up antiretroviral treatment at the district level in Malawi. Transactions of the Royal Society of Tropical Medicine and Hygiene, 100(10),980-5

2) J. Zienah, G. Kadzangra, S. Mazima, M. Borok, C. Mudwa, O. Tobaiva, M. Mucuche, S. Rusakaniko, D. Kateremba. (2006). Affordable flow cytometry for enumeration of CD4+ T lymphocytes to identify sub-type C-HIV infected adults requiring antiretroviral therapy (ART) and monitoring response to ART in a resource-limited setting. Journal of Translational Medicine 4:33

3) H. Karcher, D. Bohning, R. Downing, S. Mashate, G. Harms. (2006). Comparison of two alternative methods for CD4+ T-cell determination (Counter manual CD4 count and CyFlow) against standard dual platform flow cytometry in Uganda. Clinical Cytometry, 70B (3):163-69

4) S. Kohreandum, S. Detraira. (2012). Evaluation of the Partec Single-Platform Volumetric CyFlow® Counter System for Determining Percentage and Absolute Numbers of CD4 T Lymphocytes in HIV/AIDS Thai Patients. Thai AIDS J. 25:1-10

5) D. Wade, P.A. Diaw, G. Daneau, M. Camara, T.N. Dyeie, S. Mboup, L. Kestens. (2013). CD4 T-Cell Enumeration in a Field Setting: Evaluation of CyFlow Counter Using the CD4 Easy Count Kit-Dry and Pima CD4 Systems. PLoS ONE, 8(9): e75484

18 Contact of manufacturer

Manufacturer: Sysmex Partec GmbH

Amtdstraße 11 a-b, 02826 Görzitz, Germany

Tel +49 3581 8746 0, Fax +49 3581 8746 70, E-mail: info@sysmex-partec.com

19 Version number of IFU and date of issue

Revision: Rev-008\_05-04-2022\_CN 2452

Issued by: Sysmex Partec GmbH

Safety Data Sheet for this product is available at www.sysmex-partec.com/services.

Cy and CyDe are trademarks of GE Healthcare.

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## Product Labels

**Product:** CD4% easy count kit  
**Code No.:** 05-8405

Product Label	Container Type	Label Size
 <p>The label for the CD4% easy count kit kit includes the Sysmex logo, reference number REF 05-8405, lot number LOT, and a list of contents: 05-8405-01 Vial containing CD4 mAb PE, 05-8405-02 Vial containing CD45 mAb PE-Cy5, 05-8405-03 Bottle containing Buffer 1, and 05-8405-04 Bottle containing Buffer 2. It also features storage instructions for 2°C and 8°C, CE and IVD marks, the website www.sysmex-partec.com/services, and contact information for Sysmex Partec GmbH in Görlitz, Germany. A barcode with Sample EAN (01)04250878904382(17)201113(10)700103 and Lab-0112-2 is also present.</p>	cardboard box (white)	70 x 50 mm
 <p>The label for CD4 mAb PE includes the Sysmex logo, reference number REF 05-8405-01, lot number LOT 123456, and storage instructions for 2°C and 8°C. It also features CE and IVD marks, the website www.sysmex-partec.com/services, and contact information for Sysmex Partec GmbH in Görlitz, Germany. Lab-0113-2 is also indicated.</p>	3 ml glass vial (amber glass) with white PP cap	50 x 19 mm
 <p>The label for CD45 mAb PE-Cy5 includes the Sysmex logo, reference number REF 05-8405-02, lot number LOT 123456, and storage instructions for 2°C and 8°C. It also features CE and IVD marks, the website www.sysmex-partec.com/services, and contact information for Sysmex Partec GmbH in Görlitz, Germany. Lab-3503-1 is also indicated.</p>	3 ml glass vial (amber glass) with white PP cap	50 x 19 mm
 <p>The label for Buffer 1 includes the Sysmex logo, reference number REF 05-8405-03, lot number LOT, and storage instructions for 2°C and 8°C. It also features CE and IVD marks, the website www.sysmex-partec.com/services, and contact information for Sysmex Partec GmbH in Görlitz, Germany. Lab-0114-2 is also indicated.</p>	100 ml glass bottle (amber glass) with white PP cap	60 x 30 mm
 <p>The label for Buffer 2 includes the Sysmex logo, reference number REF 05-8405-04, lot number LOT, and storage instructions for 2°C and 8°C. It also features CE and IVD marks, the website www.sysmex-partec.com/services, and contact information for Sysmex Partec GmbH in Görlitz, Germany. Lab-0115-2 is also indicated.</p>	100 ml glass bottle (amber glass) with white PP cap	 <b>LILIANA E. PARODI</b> <b>FARMACEUTICA</b> <b>MAT. 9357</b>

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Digitally signed by GESTION DOCUMENTAL ELECTRONICA - GDE  
Date: 2023.10.31 08:46:05 -03:00

Digitally signed by GESTION DOCUMENTAL  
ELECTRONICA - GDE  
Date: 2023.10.31 08:46:06 -03:00

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