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Digitally signed by GESTION DOCUMENTAL
ELECTRONICA - GDE
Date: 2024.04.26 19:08:16 -03:00

ZenTech

NEONATAL PAP Screening ELISA

E-LC-192

English (en)

Español (es)

Français (fr)

ZenTech s.a.

Liège Science Park - Avenue Pré Aily, 10 - 4031 ANGLEUR, Belgium
Tel. : + 32-(0)4-361.42.32 - Fax : + 32-(0)4-367.00.63

info@zentech.be
www.zentech.be

ISO15223	MEDICAL DEVICES SYMBOL	SÍMBOLO DE DISPOSITIVOS MÉDICOS	SYMBOLE DES DISPOSITIFS MÉDICAUX
	STORAGE TEMPERATURE LIMITATION	LÍMITES DE TEMPERATURA DE ALMACENAMIENTO	LIMITES DE LA TEMPERATURE DE STOCKAGE
LOT	BATCH CODE	CÓDIGO DE LOTE	IDENTIFICATION DU LOT
	USE BY	UTILIZAR ANTES DE	DATE LIMITE D'UTILISATION
	CONSULT OPERATING INSTRUCTIONS	CONSULTAR LAS INSTRUCCIONES DE MANEJO O FUNCIONAMIENTO	CONSULTER LE MODE D'EMPLOI
IVD	IN VITRO DIAGNOSTIC DEVICE	DISPOSITIVO DE DIAGNÓSTICO IN VITRO	DISPOSITIF DE DIAGNOSTIC IN VITRO
	MANUFACTURED BY	FABRICADO POR	FABRIQUE PAR
REF	CATALOGUE NUMBER	NÚMERO DE CATÁLOGO	REFERENCE CATALOGUE

	SYMBOLS (EDMA recommendations)	SÍMBOLOS (recomendaciones de la EDMA)	SYMBOLES (recommandations de l'EDMA)
	Number of determinations	Número de determinaciones	Nombre de déterminations
CAL CONTROL	Calibrators & controls (dried blood spots)	Calibradores y controles, manchas de sangre	Calibrateurs et contrôles (gouttes de sang séché)
SORB MTP	Coated microtiter plate	Placa de microtitulación recubierta	Microplaque coatée
EXTR SOLN	Extraction solution	Solución de extracción	Solution d'extraction
Ab BIOT RCNS 300µL H2O	Lyophilized biotinylated anti-PAP (to be reconstituted in water)	Anticuerpo anti-PAP biotinilado liofilizado (que se deberá reconstituir en agua)	Anti-PAP biotinylé lyophilisé (à reconstituer avec de l'eau)
CONJ HRP	Enzyme conjugate	Conjugado enzimático	Conjugué enzymatique
SUBS TMB	Chromogenic substrate	Sustrato cromogénico	Substrat chromogène
H₂SO₄ 0.5 M	Blocking reagent	Reactivo de bloqueo	Solution d'arrêt
BUF WASH 10X	Washing solution (to be diluted 10-fold)	Solución de lavado (que se deberá diluir 10 veces)	Solution de lavage (à diluer 10 fois)

ENGLISH

Enzyme immunoassay for the titration of Pancreatitis-Associated Protein from dried blood spots of neonates

E-LC-192 (192 determinations: 2x (12 strips x 8 wells))

FOR IN VITRO DIAGNOSTIC (Professional use only)

1. INTENDED USE

The Neonatal PAP Screening ELISA from ZenTech is an enzyme immunoassay for the quantitative determination of pancreatitis-associated protein (PAP) in neonates using blood spot samples dried on 903® filter paper. This kit is used as an aid in identifying newborns that are at increased risk of having cystic fibrosis (CF).

2. CLINICAL APPLICATIONS

The pancreatitis-associated protein (PAP), also known as Reg3A, is a 175 amino acids secretory protein which is a member of the regeneration family of C-type lectins.^{1,2} It is absent from the healthy pancreas but synthesized in high amounts after a sustained pancreatic stress. Studies on acute pancreatitis led to consider using the PAP blood concentration for cystic fibrosis screening.³

Cystic fibrosis, the most common lethal genetic disorder affecting Caucasian population (1/2,500 births), is a multi-system illness, most frequently characterized by a childhood chronic obstructive pulmonary disease, a pancreatic exocrine insufficiency, and abnormal sweat electrolyte concentrations. The CF diagnostic is based on a combination of the above clinical findings and/or a positive family history of the illness in conjunction with an abnormal sweat test.

As PAP is already synthesized *in utero* in the CF pancreas, neonates with CF have raised PAP blood concentration in the first days of life. Similar to immunoreactive trypsinogen (IRT), PAP elevation is not strictly specific to CF.^{5,6,7} It was demonstrated that all newborns with CF show an elevation of both IRT and PAP, whereas those without CF show normal results or an elevation of IRT or PAP, but rarely both.^{4,6,7} Thus, a dried blood spot assay for PAP has the potential to screen CF in neonates, in addition to an IRT assay.

3. PRINCIPLE OF THE ASSAY

The Neonatal PAP Screening ELISA is an enzyme immunoassay to quantitate the pancreatitis-associated protein (PAP) using dried blood spots.

The Neonatal PAP Screening ELISA is a sandwich ELISA. The wells are coated with an anti-PAP antibody that captures the PAP present in the sample. After an incubation and a washing step to remove the unbound material, a biotinylated anti-PAP antibody is added. After another incubation and a second washing step, a streptavidin-horseradish peroxidase (HRP) conjugate is added to the wells and binds to the biotin with high affinity. After a third incubation and another washing step, the immuno-complex is detected by the reduction of 3,3',5,5'-tetramethylbenzidine (TMB) by the HRP. This reaction leads to a blue color development, which is directly proportional to the amount of antigen in the standard or in the sample. The enzymatic reaction is finally stopped by adding sulfuric acid and the absorbance at 450 nm is read using an ELISA microtiterplate reader.

4. REAGENTS PROVIDED WITH THE KIT

1. **Coated microtiterplates:** wells coated with polyclonal anti-human PAP antibodies.
2. **Calibrators and controls blood spots:** aluminum pouch with calibrators C0-C5 and internal controls L1-L2 (dried blood spots on 903® or 226 specimen collection paper using human matrix containing a specified amount of the human PAP).
3. **Extraction solution:** ready for use buffer with BSA and Tween® 20. Preservative: ProClin™ 300 0.05%.

4. **Biotinylated anti-PAP:** lyophilized biotinylated anti-human PAP antibody.
Rehydrate with 300 µL of fresh and germ-free distilled water (**shake min. 15 min. with vortex or rolling**). Dilute concentrated solution 1:24 with the extraction buffer, e.g. 300 µL + 6900 µL = 7200 µL. Preservative: ProClin™ 300 0.015%.
Stability after reconstitution: 12 days at 2-8°C.

The rinsing of the tube containing the antibody with the extraction buffer must be repeated several times to recover all material. This step is very important since poor rinsing will lead to poor results.

5. **Enzyme conjugate:** streptavidin-peroxidase complex in a stabilizing solution. Preservative: ProClin™ 300 0.05%. Ready for use. Conservation of the conjugate solution at low temperatures (2-8°C) may lead to the formation of a white precipitate. This precipitate dissolves under shaking at room temperature (RT).
6. **Washing solution:** 10x concentrated phosphate buffer with detergent. Preservative: Thimerosal 0.1%.
Dilute concentrated solution 1:10 with fresh and germ free distilled or deionized water, e.g. 55 mL + 495 mL = 550 mL. Conservation of the concentrated washing solution at low temperatures (2-8°C) may lead to crystals formation. These crystals dissolve upon heating (37°C) or dilution to the working concentration. Stability after dilution: 10 weeks at RT (20-25° C).
7. **Chromogenic substrate:** TMB (3,3'-5,5'-tetramethylbenzidine). Ready for use.
8. **Blocking reagent:** H₂SO₄, 0.5M. Ready for use.

Reagents	Quantity 192 tests	Physical state
Wells	2 x (12 x 8)	Ready for use
Calibrators 0-5 Controls 1-2	1 aluminum pouch	Dried blood spots
Extraction solution	1 x 55 mL	Ready for use
Biotinylated anti-PAP	4 vials	Lyophilized
Enzyme conjugate	1 x 28 mL	Ready for use
Washing solution	3 x 55 mL	10x concentrated
Chromogenic substrate	1 x 25 mL	Ready for use
Blocking reagent	1 x 28 mL	Ready for use

Additional chromogenic solution, blocking reagent and washing solution could be ordered for automated procedure using reference: E-LI-000: Newborn Screening Extra Reagents

5. STORAGE AND STABILITY OF THE KIT

- **Keep dried blood spots at -20°C in the original bag containing a desiccant.** Care should be taken to seal it tightly.
- Store the kit and all other reagents at 2-8°C in closed containers when not in use but allow all reagents to stabilize at room temperature before use.
- Microtiter wells must be stored at 2-8°C. Once the foil pouch has been opened, care should be taken to seal it tightly again.
- Unopened reagents will retain reactivity until expiration date shown on the label. Do not use reagents beyond this date.
- When opened and stored, reagents will retain reactivity two months.

6. MATERIALS REQUIRED BUT NOT PROVIDED

- Single or multichannel automatic pipettes to deliver volumes in the range of 10 to 1000 µL with an accuracy of ± 1.5% over the range 10-100 µL.
- Sealing tapes.
- For manual use, a microtiter plate reader capable of reading absorbance at 450 nm and 630 nm.
- A hole-puncher which produces 3.0 mm or 3.2 mm (1/8") discs

- A microtiterplate shaker (**900 rpm or 15 Hz – 300 rpm or 5Hz**).

$$\text{rpm} = \sqrt{\frac{\text{RCF}}{r \times 1.118}} \times 10^3$$

Where *r* is the orbital radius in mm (check the specifications of your device and, if necessary, divide the diameter by two).

RCF = 0,35g

- Blood collection cards. The minimum pre-printed information required is:

- *infant's name*
- *mother's name*
- *patient ID number*
- *date of birth*
- *sex*
- *specimen collection date*
- *submitter's identity and address*
- *physician's name and telephone number*
- *manufacturer and lot number of filter paper.*

Filter paper should be 903® or 226

7. WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic only.

In order to avoid personal and environmental contamination, the following precautions must be observed:

- Use disposable gloves while handling potentially infectious material and performing the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- All material of human origin used for the preparation of this kit is tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious; therefore, the assay waste must be decontaminated and disposed of, in accordance with established safety procedures. Disposable ignitable material must be incinerated; disposable non-ignitable material must be sterilized in autoclave for at least 1 hour at 121°C. Liquid waste must be added with sodium hypochlorite at a final concentration of 3%. The hypochlorite treatment should last for at least 30 minutes. Liquid waste containing acid must be neutralized with appropriate amounts of base before treating with sodium hypochlorite.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.
- Chromogenic substrate and blocking reagent solutions should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident, rinse thoroughly with running water.
- Washing solution contains thimerosal.
- Extraction solution, biotinylated anti-PAP and enzyme conjugate contain ProClin™ 300 (H315/H317/H318/H410). This product is highly toxic by inhalation, swallowing and contact with skin. Keep away from food and drink.
- Blocking reagent contains sulfuric acid (4.9% w/w) (H331/H314/H315)
- Wear protective clothes and gloves. In case of contact with skin or eyes, rinse thoroughly with water. In case of accident, consult immediately a physician and show him/her the product label.

Risk phrases:

- H314 causes severe skin burns and eye damage
- H315 causes skin irritation
- H317 may cause an allergic skin reaction
- H318 causes serious eye damages
- H331 toxic if inhaled
- H410 very toxic to aquatic life with long lasting effects

In order to obtain reproducible results, the following rules must be observed:

- **Using a multichannel automatic pipette is better to ensure good results**
- **Avoid any contact between the tip and the side of wells. To avoid any degradation of the active phase.**

- Do not mix reagents from different lot numbers or from other manufacturers.
- Do not freeze kits.
- Strict compliance to the specific time and temperature of incubations is recommended for accurate results.
- Allow the microtiter plate strips and reagents to stabilize at room temperature before opening and using them.
- Do not use reagents after their expiration dates.
- Incomplete or inefficient washing will cause poor precision and high background.
- Use thoroughly clean glassware, free from contamination of metal ions or oxidizing substances.
- Use distilled water, stored in clean containers.
- Microbial contaminated serum or specimens containing heavy, visible particles should not be used.
- Cross contaminations of reagents or samples could cause false results. Use a clear, fresh, disposable pipette tip for each reagent of specimen manipulation.
- Do not expose the chromogenic substrate to light during storage or incubation.
- Follow exact incubation times. Dispense chromogenic substrate and blocking reagent solutions in no more than 3-4 minutes; dispense the two reagents in the same sequence.
- Traces of hypochlorite of soda can destroy many reagents biological activity.

A variety of factors influence the assay performances. These include the accuracy and reproducibility of pipetting technique, the photometer used or the timing bias during the assay.

8. SPECIMEN COLLECTION

Blood should ideally be collected between the second and the third day of life (48 to 72 hours after birth) as there is an age-related increase in PAP levels. A specimen collection after 168 hours could lead to false-positive results.⁸

Collect blood from the infant's heel only using the medial (closest to body center line) or lateral (furthest from body center line) portion of the plantar (walking) surface. Blood collection from another area of the infant foot, e.g. arch, may result in nerve, tendon or cartilage injury.

Fill out the required information on the blood collection form or card (see section 6).

To increase local blood flow, warm the skin puncture site using a warm moist towel (<42°C) applied to the site for 2-3 minutes. Similarly, holding the infant's limb below the level of the heart will increase the venous flow.

Clean the skin with 70% (v/v) isopropyl alcohol. Wipe dry with a sterile gauze and allow the skin to air dry.

Puncture the heel with a sterile lancet (2.4-2.5 mm long) or with an automatic lancet device. Wipe away the first drop of blood with a sterile gauze.

Gently touch one side of the filter paper against a large drop of blood and, in one step, allow a sufficient quantity of blood to soak into the paper to completely fill the preprinted circle. Ensure that the blood has penetrated and saturated the paper by examining both sides of the paper. Repeat the procedure to fill the required number of preprinted circles on the specimen collection card.

Do not milk or squeeze the puncture site during collection as it may cause hemolysis or dilution of the blood with tissue fluids.

Do not apply successive drops of blood to the same preprinted circle. If the blood flow diminishes before the complete fill of the preprinted repeat the sampling in a new circle.

Do not touch or smear the blood spots.

Cord blood shall not be used as a sample.

The sample is then air-dried for 4 hours at room temperature and stored in sealed paper envelopes or containers that will provide protection from moisture, light, heat and contact with other materials. Arrange transport of the collection card to the screening laboratory within 24 hours of collection.

The receiving laboratory should store the samples at 2-8°C for short-term storage in a moisture-proof environment protected from direct light. For long-term storage, samples should be stored at below -20°C in a moisture-proof environment.

The sample discs should be punched from similar areas on each individual blood spot. Do not punch sample discs from areas that include printed marks or that are near the edges of the blood spot.

9. ASSAY PROCEDURE

The results reliability depends on strict adherence to the described procedure.

Suggested assay plate layout (included only as a guide).

	1	2	3	4	5	6	7	8	9	10	11	12
A	C0	C0	Smp									
B	C1	C1										
C	C2	C2										
D	C3	C3										
E	C4	C4										
F	C5	C5										
G	L1	L1										
H	L2	L2										

Manual procedure:

Sample preparation

- Bring all reagents and samples to room temperature before use. Gently mix the reagents for at least 15 minutes to homogenize them.
- Select the required number of microtiter strips and place them in the strip holder, taking account of the number of calibrators C0-C5, controls L1-L2 and samples to be assayed. It is recommended to analyze the calibrators and controls in duplicate.
- Following the suggested plate layout, punch 3.0 mm or 1/8" diameter blood spots from the calibrators (C0 to C5) and controls (L1 and L2) into their respective wells. Punch patient specimens into the remaining individual wells, noting the position of each patient sample on the plate. Spots should be punched in the center of the discs.

Do not autoclave calibrators, internal controls or patients' samples prior to punching.

PAP extraction

- Add 100 µL of extraction solution to each well.

Each well should be visually examined to ensure the blood spots are completely immersed in the solution.

Cover wells with a sealing tape to prevent evaporation.

- Incubate the plate for 1h15 at 37°C on a microtiter plate shaker at **900 rpm**. (see section 6) Ensure by periodic visual inspection that each disc is fully immersed in the extraction solution.
- During this incubation, prepare the working washing solution and the biotinylated anti-PAP solution (See section 4).

- After incubation, empty all wells completely by plate inversion. Ensure that all blood spots are discarded.**

A video of an appropriate elimination of biological samples is available:

https://youtu.be/ozS_ADEVnJE.

- Wash each well 3 times with 350 µL of diluted washing solution. **A soft flushing rate is required for good results (dispensing rate of 100 µL/s is recommended).**

After each washing step, ensure that all the wells are empty (remaining liquid < 15 µL). After the last washing, carefully remove remaining fluid by tapping the strips on tissue paper prior to next step. **Be sure to set up the washer. Please note that incomplete or inefficient washing will cause poor precision and high background.**

Immuno-complex formation

- Add 100 µL of the diluted biotinylated anti-PAP solution to each well.
- Cover wells with a sealing tape.
- Incubate for 2h at 37°C on a microtiter plate shaker (the recommended speed is **300 rpm**).
- After the incubation, repeat washing as described above.

Amplification reaction

- Add 100 µL of the ready for use enzyme conjugate to each well.
- Cover wells with a sealing tape.
- Incubate for 20 min at room temperature (15-30°C) on a microtiter plate shaker (the recommended speed is **300 rpm**).
- After the incubation, repeat washing as described above.

Color generation

- Dispense 100 µL of the chromogenic substrate to each well.
- Incubate for 15 minutes at room temperature (15-30°C) keeping the plate away from light, without shaking.
- After the incubation, stop the enzymatic reaction by adding 100 µL of the blocking reagent to each well.

Plate reading

- Measure the absorbance at 450 nm. Dual wavelength reading using 630 nm as reference wavelength is recommended but not compulsory.

Absorbance should be read within 20 minutes after completing the assay.

Wells Reagents	Calibrator (C0-C5)	Control (L1-L2)	Samples
Calibrators	1 x 3.0mm	----	----
Controls	----	1 x 3.0mm	----
Samples	----	----	1 x 3.0mm
Extraction solution	100 µL	100 µL	100 µL
- Cover the plate and elute for 1h15 at 900 rpm and 37°C - Aspirate and wash 3 x 350 µL			
Biotinylated anti-PAP	100 µL	100 µL	100 µL
- Cover the plate and incubate for 2h at 300 rpm and 37°C - Aspirate and wash 3 x 350 µL			
Enzyme conjugate	100 µL	100 µL	100 µL
- Cover the plate and incubate for 20 min at 300 rpm and RT - Aspirate and wash 3 x 350 µL			
Chrom. substrate	100 µL	100 µL	100 µL
- Cover the plate and incubate for 15 min at RT away from light			
Blocking reagent	100 µL	100 µL	100 µL
- Read: 450 nm and 630 nm as reference			

Automate procedure: Timing on automate: 4h30

Prepare the samples following the instructions given in the manual procedure (See section sample preparation).

Prepare the working washing solution and biotinylated anti-PAP (See section 4).

The automate procedure follows the protocol provided with the machine.

Be careful that the punches will have to be manually removed from the microtiter plate after the first incubation.

Additional chromogenic substrate, blocking reagent and washing solution could be ordered for automated procedure using reference:

E-LI-000: Newborn Screening Extra Reagents.

10. CALCULATION OF RESULTS

Titration curve

CHECK THE ACCURATE VALUES OF THE CALIBRATORS AND CONTROLS ON THE INSERT IN THE BOX OF EACH LOT.

The PAP amount in the patient's sample is calculated from a standard curve prepared from a known amount of human PAP calibrated by an in-house method. Currently, no standardized reference material is available for PAP.

The standard curve may change between lots. **The label on the current lot of calibrators and controls should be consulted for calibration values to be used in calculations.**

PAP concentration values are expressed in µg/L.


Use a quadratic regression curve to calculate the results after averaging the absorbance duplicates (if applicable).

- X-axis: use concentration on linear way
- Y-axis: use optical density on linear way

Results of a typical assay are shown below. **This curve is for illustration only and shall not be used to calculate patients' results.**

The maximum concentration that can be reported is limited by the linear performance characteristics of the photometer used. If the OD value of the highest calibrator is above the range of the photometer, then this calibrator must be omitted from the plot of the calibrator curve. Similarly, any sample measuring above the range of the microplate reader should be simply noted as greater than the highest acceptable calibrator or must be dilute.

Calibrator	[PAP] µg/L	OD 450nm
Cal 0	0	0.101
Cal 1	2.2	0.346
Cal 2	4.3	0.559
Cal 3	5.8	0.712
Cal 4	16.4	1.662
Cal 5	30.7	2.649

11. QUALITY CONTROL

Internal controls (L1-L2) included in the kit should be routinely monitored to check that measured concentrations are within the stated values. These controls provide valuable information regarding the validity of the test according to manufacturer specifications.

12. SUGGESTED NORMAL RANGE

In order to estimate normal range and cut-off values, 1,383 dried blood samples were tested. These samples came from the neonates' screening program in Liège Province (Belgium). These samples are all negative for CF.

An initial approach using the 99th percentile allow to determine a provisional cut-off at 4.06 µg/L.

The suggested normal range and cut-off are guidelines only. Each laboratory should establish specific thresholds and range based on performance of the assay in their laboratory and with their constituent demographic population.

N	Median (µg/L)	Percentile values			
		95th	97.5th	99th	99.5th
1,383	1.24	2.73	3.14	4.06	5.69

13. PERFORMANCES OF THE ASSAY

Performances are established following the automate procedure.

a. Detection limits

The limit of detection (LoD) for PAP is 0.559 µg/L, determined according to the guidelines in CLSI document EP17-A2 and with proportions of α and β error risks less than 5% based on 120 determinations, with 60 blanks and 60 low-level samples: LoB = 0.304 µg/L.

b. Precision

The Neonatal PAP Screening ELISA was evaluated for precision according to CLSI document EP05-A3.

It was achieved with internal controls, consisting of blood enriched with PAP at three different levels. The precision definition was based on 80 determinations (20 days, 2 runs per day, 2 replicates per run) performed with three lots of reagents.

Sample	Mean (µg/L)	Repeatability		Between-Run		Between-Day		Within-Laboratory	
		SD	CV	SD	CV	SD	CV	SD	CV
1	3.07	0.47	15.2	0.28	9.3	0.22	7.2	0.52	17.0
2	5.32	0.57	10.7	0.40	7.5	0.27	5.2	0.80	15.1
3	7.16	0.61	8.6	0.58	8.1	0.46	6.4	0.97	13.5

c. Interferences

The Neonatal PAP Screening ELISA was evaluated for interference according to CLSI document EP7-A2.

The following substances, when added to blood enriched with 1.5 µg/L of PAP, were found not to interfere at the concentration indicated. These include neonates' medication, sample additives, abnormal biochemical metabolites and one PAP structural analog.

A bias less than ± 0.2 µg/L is not considered as a significant interference.

Substance tested	Test concentration	Comments
Amoxicillin	95 mg/L	3x median serum concentration ⁹
Cefotaxime	306 mg/L	CLSI EP7-A2, Appendix C
Acyclovir	300 mg/L	3x median serum concentration ¹⁰
EDTA	1.3 mg/L	CLSI EP7-A2, Appendix C
Heparin	3,000 U/L	CLSI EP7-A2, Appendix C
Bilirubin	150 mg/L	CLSI EP7-A2, Appendix C
Vitamin K ₁	5 mg/L	3x neonates' plasma concentration ¹¹
Reg1A	1.3 µg/L	3x median serum concentration ¹²

The following substance, when added to blood enriched with 1.5 µg/L of PAP, was found to interfere above the indicated concentration.

Substance tested	Interfering concentration	Comments
Sodium citrate	8.08 g/L	Short draw*

*if the collection tube is correctly filled, the blood sodium citrate concentration is between 3 and 5.5 g/L (1:5 or 1:10 total volume).

14. LIMITATIONS OF THE PROCEDURE

- The Neonatal PAP Screening ELISA is a screening method for measuring the PAP concentration in newborn blood spot specimens. This assay is designed to be used with dried blood specimens that are collected on 903® and 2226 filter paper.
- Elevated results are not a diagnostic per se of cystic fibrosis but indicate the need for further study of the newborn from which a presumptive positive specimen was received.
- To ensure accurate and reliable results, make sure that all blood spot disks are within the extraction solution during the incubation period.
- Strict compliance to the protocol is advised to obtain reliable results. Any modification or change made to the kit or the assay procedure is under the responsibility of the user.
- Do not use cord blood.
- Sex, birth weight and gestational age can lead to small differences in PAP concentrations but do not significantly affect the screening protocol.⁸
- A specimen collected between 48 and 72 hours of age is optimal. There is an age-related increase in PAP levels and a specimen collection after 168 hours could lead to false-positive results.⁸
- Blood transfusion leads to clinically significant higher PAP levels and to a higher risk on a false-positive screening test result. It is probably due to the much higher levels of PAP found in adult blood.⁸
- Sick infants may have a false-positive screen due to increased stress on the body.⁸
- Premature infants or newborns with mild CF phenotype may exhibit false negative result due to limited pancreatic damage or immaturity of the pancreas.^{4,8}
- Newborns with meconium ileus may exhibit false negative result.³

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15. BIBLIOGRAPHY-BIBLIOGRAFÍA-BIBLIOGRAPHIE

1. ORELLE, Beatrice, KEIM, Volker, MASCIOTRA, Luis, *et al.* Human pancreatitis-associated protein. Messenger RNA cloning and expression in pancreatic diseases. *The Journal of clinical investigation*, 1992, vol. 90, no 6, p. 2284-2291.
2. PARIKH, Abhirath, STEPHAN, Anne-Fleur, et TZANAKAKIS, Emmanuel S. Regenerating proteins and their expression, regulation, and signaling. *Biomolecular concepts*, 2012, vol. 3, no 1, p. 57-70.
3. SARLES, Jacques, BARTHELLEMY, Sandrine, FÉREC, Claude, *et al.* Blood concentrations of pancreatitis associated protein in neonates: relevance to neonatal screening for cystic fibrosis. *Archives of Disease in Childhood-Fetal and Neonatal Edition*, 1999, vol. 80, no 2, p. F118-F122.
4. SOMMERBURG, Olaf, LINDNER, Martin, MUCKENTHALER, Martina, *et al.* Initial evaluation of a biochemical cystic fibrosis newborn screening by sequential analysis of immunoreactive trypsinogen and pancreatitis-associated protein (IRT/PAP) as a strategy that does not involve DNA testing in a Northern European population. *Journal of inherited metabolic disease*, 2010, vol. 33, no 2, p. 263-271.
5. SOMMERBURG, Olaf, HAMMERMANN, Jutta, LINDNER, Martin, *et al.* Five years of experience with biochemical cystic fibrosis newborn screening based on IRT/PAP in Germany. *Pediatric pulmonology*, 2015, vol. 50, no 7, p. 655-664.
6. FARRELL, Philip M., WHITE, Terry B., HOWENSTINE, Michelle S., *et al.* Diagnosis of cystic fibrosis in screened populations. *The Journal of pediatrics*, 2017, vol. 181, p. S33-S44. e2.
7. SARLES, Jacques, BERTHÉZÈNE, Patrice, LE LOUARN, Christian, *et al.* Combining immunoreactive trypsinogen and pancreatitis-associated protein assays, a method of newborn screening for cystic fibrosis that avoids DNA analysis. *The Journal of pediatrics*, 2005, vol. 147, no 3, p. 302-305.
8. VERNOOIJ-VAN LANGEN, Annette MM, LOEBER, J. Gerard, ELVERS, Bert, *et al.* The influence of sex, gestational age, birth weight, blood transfusion, and timing of the heel prick on the pancreatitis-associated protein concentration in newborn screening for cystic fibrosis. *Journal of inherited metabolic disease*, 2013, vol. 36, no 1, p. 147-154.
9. GRAS-LE GUEN, Christele, BOSCHER, Cecile, GODON, Nathalie, *et al.* Therapeutic amoxicillin levels achieved with oral administration in term neonates. *European journal of clinical pharmacology*, 2007, vol. 63, no 7, p. 657-662.
10. TOD, M., LOKIEC, F., BIDAULT, R., *et al.* Pharmacokinetics of oral acyclovir in neonates and in infants: a population analysis. *Antimicrobial agents and chemotherapy*, 2001, vol. 45, no 1, p. 150-157.
11. LIPPI, Giuseppe et FRANCHINI, Massimo. Vitamin K in neonates: facts and myths. *Blood transfusion*, 2011, vol. 9, no 1, p. 4.
12. BACON, Siobhan, KYITHAR, Ma Peyh, SCHMID, Jasmin, *et al.* Serum levels of pancreatic stone protein (PSP)/reg1A as an indicator [...]. *BMC endocrine disorders*, 2012, vol. 12, no 1, p. 13.
13. Joanne V Mei, Sherri D Zobel, Elisabeth M Hall, Victor R De Jesus, Barbara W Adam & W Harry Hannon. Performance properties of filter paper devices for whole blood collection. *Bioanalysis* (2010) 2 (8) 1397-1403

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



Digitally signed by GESTION DOCUMENTAL ELECTRONICA - GDE
Date: 2024.04.09 11:42:38 -03:00

Digitally signed by GESTION DOCUMENTAL
ELECTRONICA - GDE
Date: 2024.04.09 11:42:39 -03:00

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