



República Argentina - Poder Ejecutivo Nacional
2018 - Año del Centenario de la Reforma Universitaria

Disposición

Número: DI-2018-3638-APN-ANMAT#MS

CIUDAD DE BUENOS AIRES
Martes 17 de Abril de 2018

Referencia: 1-47-3110-5247/16-9

VISTO el expediente N° 1-47-3110-5247/16-9 del Registro de la Administración Nacional de Medicamentos Alimentos y Tecnología Médica y,

CONSIDERANDO:

Que por los presentes actuados la firma TECNOLAB S.A. solicita la modificación del registro del Producto para diagnóstico de uso "in vitro" denominado: 1) ANTICUERPOS MONOCLONALES Y POLICLONALES DISEÑADOS PARA SER UTILIZADOS EN TÉCNICAS DE INMUNOHISTOQUÍMICA; 2) CONTROLES POSITIVOS; 3) REACTIVOS AUXILIARES.

Que lo solicitado se encuadra dentro de los alcances de la Disposición ANMAT N° 2674/99 y la documentación aportada ha satisfecho los requisitos de la normativa aplicable.

Que en el expediente de referencia consta el informe técnico producido por el Servicio de Productos para Diagnóstico que establece que se autoriza la modificación solicitada.

Que la Dirección Nacional de Productos Médicos ha tomado la intervención de su competencia.

Que la presente se dicta en virtud de las facultades conferidas por los Decretos N° 1490/92 el por el Decreto N° 101 de fecha 16 de diciembre de 2015.

Por ello;

EL ADMINISTRADOR NACIONAL DE LA ADMINISTRACIÓN NACIONAL DE
MEDICAMENTOS, ALIMENTOS Y TECNOLOGÍA MÉDICA

DISPONE:

ARTÍCULO 1°.- Autorízase la modificación del Certificado N° 8342 del producto para diagnóstico de uso in vitro denominado: 1) ANTICUERPOS MONOCLONALES Y POLICLONALES DISEÑADOS PARA SER UTILIZADOS EN TÉCNICAS DE INMUNOHISTOQUÍMICA; 2) CONTROLES POSITIVOS; 3) REACTIVOS AUXILIARES.

ARTICULO 2º.- Acéptase la incorporación de una nueva familia de productos que se detallan en el listado que consta al pie de la presente Disposición y forma parte integrante de la misma.

ARTÍCULO 3º.- Dése de baja al registro otorgado mediante Certificado N° 007847 y al PM 1252-77 pertenecientes a la firma TECNOLAB S.A.

ARTICULO 4º.- Practíquese la atestación de la presente disposición al Certificado de Inscripción N° 008342.

ARTÍCULO 5º.- En los nuevos proyectos de Manual de Instrucciones deberán constar las modificaciones descritas en el artículo 2º precedente

ARTÍCULO 6º.- Regístrese. Inscríbase en el Registro Nacional de Productores y Productos de Tecnología Médica al nuevo producto. Por el Departamento de Mesa de Entrada, notifíquese al interesado, haciéndole entrega de la presente Disposición. Gírese a la Dirección de Gestión de Información Técnica a los fines de confeccionar el legajo correspondiente. Cumplido, archívese.

LISTADO DE PRODUCTOS

1. ANTICUERPOS

Nº Catalog.	Descripción
223A-1 series	A-1-Antitrypsin Rabbit Polyclonal Antibody
206A-7 series	ACTH Rabbit Polyclonal Antibody
201M-9 series	Actin, Muscle Specific (HHF35) Mouse Monoclonal Antibody
202M-9 series	Actin, Smooth Muscle (1A4) Mouse Monoclonal Antibody
212M-1 series	Adenovirus (20/11 & 2/6) Mouse Monoclonal Antibody
393A-1 series	Adipophilin Rabbit Polyclonal Antibody
203A-1 series	Alpha-Fetoprotein Rabbit Polyclonal Antibody
203R-1 series	Alpha-Fetoprotein (EP209) Rabbit Monoclonal Primary Antibody
400M-1 series	ALDH1A1 (44) Mouse Monoclonal Antibody
204M-1 series	ALK Protein (ALK-1) Mouse Monoclonal Antibody
200R-2 series	Androgen Receptor (EP120) Rabbit Monoclonal Primary Antibody
200M-1 series	Androgen Receptor (AR441) Mouse Monoclonal Antibody
200R-1 series	Androgen Receptor (SP107) Rabbit Monoclonal Antibody
221M-1 series	Annexin A1 (MRQ-3) Mouse Monoclonal Antibody
380R-1 series	Arginase-1 (SP156) Rabbit Monoclonal Antibody
380R-2 series	Arginase-1 (EP261) Rabbit Monoclonal Primary Antibody
225M-1 series	BCA-225 (Cu-18) Mouse Monoclonal Antibody
226M-9 series	BCL2 (124) Mouse Monoclonal Antibody
226R-1 series	BCL2 (E17) Rabbit Monoclonal Antibody
226R-2 series	BCL2 (SP66) Rabbit Monoclonal Antibody
227M-9 series	BCL6 (GI191E/A8) Mouse Monoclonal Antibody
227R-2 series	BCL6 (EP278) Rabbit Monoclonal Primary Antibody
224M-1 series	Beta-Catenin (14) Mouse Monoclonal Antibody

224A-1 series Beta-Catenin Rabbit Polyclonal Antibody
228M-1 series BG8, LewisY (F3) Mouse Monoclonal Antibody
294M-1 series BOB.1 (MRQ-35) Mouse Monoclonal Antibody
249R-1 series BOB.1 (SP92) Rabbit Monoclonal Antibody
395R-1 series c-Myc (EP121) Rabbit Monoclonal Antibody
403A-7 series C3d Rabbit Polyclonal Antibody
404A-1 series C4d Rabbit Polyclonal Antibody
404R-1 series C4d (SP91) Rabbit Monoclonal Antibody
325M-1 series CA-125 (OC125) Mouse Monoclonal Antibody
399M-1 series CA19-9 (121SLE) Mouse Monoclonal Antibody
378R-1 series Cadherin-17 (SP183) Rabbit Monoclonal Antibody
229A-1 series Calcitonin Rabbit Polyclonal Antibody
229R-1 series Calcitonin (SP17) Rabbit Monoclonal Antibody
230R-1 series Caldesmon (E89) Rabbit Monoclonal Antibody
231M-1 series Calponin (CALP) Mouse Monoclonal Antibody
231R-1 series Calponin-1 (EP798Y) Rabbit Monoclonal Antibody
232A-7 series Calretinin Rabbit Polyclonal Antibody
232R-1 series Calretinin (SP13) Rabbit Monoclonal Antibody
379M-1 series Carbonic Anhydrase IX (CA IX) (MRQ-54) Mouse Monoclonal Antibody
379R-1 series Carbonic Anhydrase IX (CA IX) (EP161) Rabbit Monoclonal Primary
Antibody
402M-1 series Cathepsin K (3F9) Mouse Monoclonal Antibody
412M-1 series Caveolin-1 (2297) Mouse Monoclonal Antibody
101R-1 series CD1a (EP3622) Rabbit Monoclonal Antibody
101M-1 series CD1a (MTB1) Mouse Monoclonal Antibody
102R-1 series CD2 (EP222) Rabbit Monoclonal Primary Antibody
102M-1 series CD2 (MRQ-11) Mouse Monoclonal Antibody
103R-9 series CD3 (MRQ-39) Rabbit Monoclonal Antibody
103A-7 series CD3 Rabbit Polyclonal Antibody
104R-2 series CD4 (EP204) Rabbit Monoclonal Antibody
104R-1 series CD4 (SP35) Rabbit Monoclonal Antibody
205M-1 series CD5 (4C7) Mouse Monoclonal Antibody
205R-2 series CD5 (EP77) Rabbit Monoclonal Antibody
205R-1 series CD5 (SP19) Rabbit Monoclonal Antibody
107M-2 series CD7 (MRQ-56) Mouse Monoclonal Antibody
107R-1 series CD7 (EP132) Rabbit Monoclonal Primary Antibody
108M-9 series CD8 (C8/144B) Mouse Monoclonal Antibody
108R-1 series CD8 (SP16) Rabbit Monoclonal Antibody
110M-1 series CD10 (56C6) Mouse Monoclonal Antibody
111M-1 series CD11c (5D11) Mouse Monoclonal Antibody
111R-1 series CD11c (EP157) Rabbit Monoclonal Primary Antibody
113R-1 series CD13 (SP187) Rabbit Monoclonal Antibody
114R-1 series CD14 (EPR3653) Rabbit Monoclonal Antibody
115M-1 series CD15 (MMA) Mouse Monoclonal Antibody
116R-1 series CD16 (SP175) Rabbit Monoclonal Antibody
119M-1 series CD19 (MRQ-36) Mouse Monoclonal Antibody
119R-1 series CD19 (EP169) Rabbit Monoclonal Primary Antibody
120M-8 series CD20 (L26) Mouse Monoclonal Antibody
120R-1 series CD20 (SP32) Rabbit Monoclonal Antibody

121M-1 series CD21 (2G9) Mouse Monoclonal Antibody
121R-1 series CD21 (EP3093) Rabbit Monoclonal Antibody
123R-2 series CD23 (EP75) Rabbit Monoclonal Primary Antibody
123M-2 series CD23 (MRQ-57) Mouse Monoclonal Antibody
123R-1 series CD23 (SP23) Rabbit Monoclonal Antibody
125M-1 series CD25 (4C9) Mouse Monoclonal Antibody
130M-9 series CD30 (Ber-H2) Mouse Monoclonal Antibody
131M-1 series CD31 (1A10) Mouse Monoclonal Antibody
131R-2 series CD31 (EP78) Rabbit Monoclonal Primary Antibody
131M-9 series CD31 (JC70) Mouse Monoclonal Antibody
133M-1 series CD33 (PWS44) Mouse Monoclonal Antibody
134M-1 series CD34 (QBEnd/10) Mouse Monoclonal Antibody
134R-1 series CD34 (EP88) Rabbit Monoclonal Primary Antibody
135M-1 series CD35 (RLB25) Mouse Monoclonal Antibody
118R-1 series CD38 (SP149) Rabbit Monoclonal Antibody
143M-1 series CD43 (MT1) Mouse Monoclonal Antibody
143R-1 series CD43 (SP55) Rabbit Monoclonal Antibody
144M-9 series CD44 (MRQ-13) Mouse Monoclonal Antibody
145M-9 series CD45 (LCA) (2B11 & PD7/26) Mouse Monoclonal Antibody
146M-1 series CD45R (MB1) Mouse Monoclonal Antibody
147M-9 series CD45RO (UCHL-1) Mouse Monoclonal Antibody
156M-8 series CD56 (123C3.D5) Mouse Monoclonal Antibody
156R-9 series CD56 (MRQ-42) Rabbit Monoclonal Antibody
157M-9 series CD57 (NK-1) Mouse Monoclonal Antibody
161M-1 series CD61 (2f2) Mouse Monoclonal Antibody
263M-1 series CD63 (NK1/C3) Mouse Monoclonal Antibody
168M-9 series CD68 (Kp-1) Mouse Monoclonal Antibody
171M-9 series CD71 (MRQ-48) Mouse Monoclonal Antibody
174M-1 series CD74 (LN2) Mouse Monoclonal Antibody
174R-1 series CD74 (EP167) Rabbit Monoclonal Primary Antibody
179R-2 series CD79a (EP82) Rabbit Monoclonal Primary Antibody
179M-9 series CD79a (JCB117) Mouse Monoclonal Antibody
179R-1 series CD79a (SP18) Rabbit Monoclonal Antibody
199R-1 series CD99 (EPR3097Y) Rabbit Monoclonal Antibody
199M-1 series CD99 (H036-1.1) Mouse Monoclonal Antibody
437R-1 series CD103 (EP206) Rabbit Monoclonal Primary Antibody
117R-1 series CD117, c-kit (YR145) Rabbit Monoclonal Antibody
198M-1 series CD123 (6H6) Mouse Monoclonal Antibody
138R-2 series CD138 (EP201) Rabbit Monoclonal Primary Antibody
138M-1 series CD138/syndecan-1 (B-A38) Mouse Monoclonal Antibody
163M-1 series CD163 (MRQ-26) Mouse Monoclonal Antibody
235R-1 series CDX-2 (EPR2764Y) Rabbit Monoclonal Antibody
236M-9 series CEA (CEA31) Mouse Monoclonal Antibody
236A-1 series CEA Rabbit Polyclonal Antibody
238M-9 series Chromogranin A (LK2H10) Mouse Monoclonal Antibody
424M-1 series CITED1 (5H6) Mouse Monoclonal Antibody
359A-1 series Claudin 1 Rabbit Polyclonal Antibody
418M-1 series Claudin 7 (5D10F3) Mouse Monoclonal Antibody

213M-1 series CMV Mouse Monoclonal Antibody
239M-1 series Collagen Type IV (CIV22) Mouse Monoclonal Antibody
240R-1 series COX-2 (SP21) Rabbit Monoclonal Antibody
241R-4 series Cyclin D1 (EP12) Rabbit Monoclonal Antibody
241R-1 series Cyclin D1 (SP4) Rabbit Monoclonal Antibody
334M-8 series Cytokeratin (34betaE12) Mouse Monoclonal Antibody
335M-9 series Cytokeratin (35betaH11) Mouse Monoclonal Antibody
452M-1 series Cytokeratin (CAM 5.2) Mouse Monoclonal Antibody
300M-1 series Cytokeratin (OSCAR) Mouse Monoclonal Antibody
305R-1 series Cytokeratin 5 (EP1601Y) Rabbit Monoclonal Antibody
356M-1 series Cytokeratin 5 & 6 (D5 & 16B4) Mouse Monoclonal Antibody
307M-9 series Cytokeratin 7 (OV-TL 12/30) Mouse Monoclonal Antibody
818M-9 series Cytokeratin 8 & 18 Mouse Monoclonal Antibody
410R-1 series Cytokeratin 10 (EP97) Rabbit Monoclonal Antibody
314M-1 series Cytokeratin 14 (LL002) Mouse Monoclonal Antibody
314R-1 series Cytokeratin 14 (SP53) Rabbit Monoclonal Antibody
317M-1 series Cytokeratin 17 (Ks 17.E3) Mouse Monoclonal Antibody
317R-1 series Cytokeratin 17 (EP98) Rabbit Monoclonal Primary Antibody
319M-1 series Cytokeratin 19 (A53-B/A2.26) Mouse Monoclonal Antibody
320M-1 series Cytokeratin 20 (Ks20.8) Mouse Monoclonal Antibody
313M-1 series Cytokeratin Cocktail Mouse Monoclonal Antibody
303M-1 series Cytokeratin, HMW (AE3) Mouse Monoclonal Antibody
301M-1 series Cytokeratin, LMW (AE1) Mouse Monoclonal Antibody
243M-1 series Desmin (D33) Mouse Monoclonal Antibody
243R-1 series Desmin (EP15) Rabbit Monoclonal Primary Antibody
436R-1 series Desmoglein 3 (EP306) Rabbit Monoclonal Primary Antibody
244R-1 series DOG1 (SP31) Rabbit Monoclonal Antibody
246R-1 series E-cadherin (EP700Y) Rabbit Monoclonal Antibody
414R-2 series EGFR (EP22) Rabbit Monoclonal Primary Antibody
414R-1 series EGFR (SP84) Rabbit Polyclonal Antibody
247M-9 series EMA (E29) Mouse Monoclonal Antibody
248M-1 series Ep-CAM/Epithelial Specific Antigen (MOC-31) Mouse Monoclonal Antibody
248M-9 series Ep-CAM/Epithelial Specific Antigen (Ber-EP4) Mouse Monoclonal Antibody
245M-1 series Epstein-Barr Virus (CS1-4) Mouse Monoclonal Antibody
245R-1 series Epstein-Barr Virus (MRQ-47) Rabbit Monoclonal Antibody
434R-1 series ERG (EP111) Rabbit Monoclonal Primary Antibody
249R-2 series Estrogen Receptor (EP1) Rabbit Monoclonal Antibody
249R-1 series Estrogen Receptor (SP1) Rabbit Monoclonal Antibody
415M-1 series EZH2 (11) Mouse Monoclonal Antibody
250A-1 series Factor VIII-R Ag. Rabbit Polyclonal Antibody
251M-1 series Factor XIIIa (AC-1A1) Mouse Monoclonal Antibody
251R-1 series Factor XIIIa (EP3372) Rabbit Monoclonal Antibody
252M-1 series Fascin (55k-2) Mouse Monoclonal Antibody
254M-1 series FLI-1 (MRQ-1) Mouse Monoclonal Antibody
405M-1 series FOXA1 (2F83) Mouse Monoclonal Antibody
250R-2 series FOXP1 (EP137) Rabbit Monoclonal Primary Antibody
350R-1 series FOXP1 (SP133) Rabbit Monoclonal Antibody

207A-7 series FSH Rabbit Polyclonal Antibody
207R-1 series FSH (EP257) Rabbit Monoclonal Primary Antibody
255M-1 series Galectin-3 (9C4) Mouse Monoclonal Antibody
256A-1 series Gastrin Rabbit Polyclonal Antibody
417M-1 series GATA1 (4F5) Mouse Monoclonal Antibody
390M-1 series GATA3 (L50-823) Mouse Monoclonal Antibody
257M-1 series GCDFP-15 (23A3) Mouse Monoclonal Antibody
257R-1 series GCDFP-15 (EP1582Y) Rabbit Monoclonal Antibody
906H-0 series GCDFP-15 (23A3) + Mammaglobin Cocktail (304-1A5 & 31A5)
208R-1 series GH (EP267) Rabbit Monoclonal Primary Antibody
208A-7 series GH Rabbit Polyclonal Antibody
258M-1 series Glial Fibrillary Acidic Protein (G-A-5) Mouse Monoclonal Antibody
258R-1 series Glial Fibrillary Acidic Protein (EP672Y) Rabbit Monoclonal Antibody
258R-2 series Glial Fibrillary Acidic Protein (SP78) Rabbit Monoclonal Antibody
259R-1 series Glucagon (EP74) Rabbit Monoclonal Primary Antibody
259A-1 series Glucagon Rabbit Polyclonal Antibody
355A-1 series GLUT1 Rabbit Polyclonal Antibody
413A-1 series GLUT3 Rabbit Polyclonal Antibody
389M-1 series Glutamine Synthetase (GS-6) Mouse Monoclonal Antibody
260M-1 series Glycophorin A (GA-R2) Mouse Monoclonal Antibody
261M-9 series Glypican-3 (1G12) Mouse Monoclonal Antibody
262A-1 series Granzyme B Rabbit Polyclonal Antibody
262R-1 series Granzyme B (EP230) Rabbit Monoclonal Primary Antibody
234A-1 series hCG Rabbit Polyclonal Antibody
398M-1 series Heat Shock Protein 27 (G3.1) Mouse Monoclonal Antibody
215A-7 series Helicobacter pylori Rabbit Polyclonal Antibody
360R-1 series Hemoglobin A (EPR3608) Rabbit Monoclonal Antibody
283M-1 series HBME-1 (HBME-1) Mouse Monoclonal Antibody
216A-1 series Hepatitis B Virus Core Antigen Rabbit Polyclonal Antibody
217M-2 series Hepatitis B Virus Surface Antigen (A10F1) Mouse Monoclonal Antibody
217M-1 series Hepatitis B Virus Surface Antigen (S1-210) Mouse Monoclonal Antibody
264M-9 series Hepatocyte Specific Antigen (Hep-Par1) (OCH1E5) Mouse Monoclonal Antibody
264R-1 series Hepatocyte Specific Antigen (Hep-Par1) (EP265) Rabbit Monoclonal
237M-1 series Her2/Neu (c-erbB-2) (CB-11) Mouse Monoclonal Antibody
237R-1 series Her2/Neu (SP3) Rabbit Monoclonal Antibody
237R-2 series Her2/Neu (EP3) Rabbit Monoclonal Antibody
361A-1 series Herpes Simplex Virus I Rabbit Polyclonal Antibody
361M-1 series Herpes Simplex Virus I (10A3) Mouse Monoclonal Antibody
362A-1 series Herpes Simplex Virus II Rabbit Polyclonal Antibody
375M-9 series HGAL (MRQ-49) Mouse Monoclonal Antibody
265M-1 series HHV-8 (13B10) Mouse Monoclonal Antibody
282M-9 series HMB-45 (HMB-45) Mouse Monoclonal Antibody
904H-0 series HMB-45 + MART-1 (Melan A) (A103) + Tyrosinase (T311)
266A-1 series Human Placental Lactogen (hPL) Rabbit Polyclonal Antibody
267R-1 series IgA (EP170) Rabbit Monoclonal Primary Antibody
267A-1 series IgA Rabbit Polyclonal Antibody
268R-1 series IgD (EP173) Rabbit Monoclonal Primary Antibody
268A-1 series IgD Rabbit Polyclonal Antibody

269A-1 series IgG Rabbit Polyclonal Antibody
367R-1 series IgG4 (EP138) Rabbit Monoclonal Primary Antibody
367M-1 series IgG4 (MRQ-44) Mouse Monoclonal Antibody
270A-1 series IgM Rabbit Polyclonal Antibody
433R-1 series IMP3 (EP286) Rabbit Monoclonal Primary Antibody
271M-1 series Inhibin, alpha (R1) Mouse Monoclonal Antibody
271R-2 series Inhibin, alpha (MRQ-63) Rabbit Monoclonal Antibody
272M-1 series INI-1 (MRQ-27) Mouse Monoclonal Antibody
273A-1 series Insulin Guinea Pig Polyclonal Antibody
431R-1 series Islet-1 (EP283) Rabbit Monoclonal Primary Antibody
274M-9 series Kappa (L1C1) Mouse Monoclonal Antibody
274R-1 series Kappa (EP171) Rabbit Monoclonal Primary Antibody
366M-9 series KBA.62 (KBA.62) Mouse Monoclonal Antibody
275R-1 series Ki-67 (SP6) Rabbit Monoclonal Antibody
276M-9 series Ksp-Cadherin (MRQ-33) Mouse Monoclonal Antibody
277M-9 series Lambda (Lamb14) Mouse Monoclonal Antibody
277R-1 series Lambda (EP172) Rabbit Monoclonal Primary Antibody
392M-1 series Langerin (12D6) Mouse Monoclonal Antibody
442R-1 series LEF1 (EP310) Rabbit Monoclonal Primary Antibody
209A-1 series LH Rabbit Polyclonal Antibody
464R-1 series LIN28 (EP150) Rabbit Monoclonal Primary Antibody
370R-1 series LMO2 (SP51) Rabbit Monoclonal Antibody
278R-1 series Lysozyme (EP134) Rabbit Monoclonal Primary Antibody
278A-1 series Lysozyme Rabbit Polyclonal Antibody
279M-1 series Macrophage (HAM-56) Mouse Monoclonal Antibody
280R-1 series Mammaglobin (31A5) Rabbit Monoclonal Antibody
280R-2 series Mammaglobin (EP249) Rabbit Monoclonal Primary Antibody
280C-1 series Mammaglobin Cocktail (304-1A5 & 31A5) Mouse/Rabbit Monoclonal Antibodies
281M-8 series MART-1 (Melan A) (A103) Mouse Monoclonal Antibody
281M-9 series MART-1 (Melan A) (M2-7C10) Mouse Monoclonal Antibody
903H-0 series MART-1 (M2-7C10) + Tyrosinase (T311) Mouse Monoclonal Antibody
435R-1 series MCM3 (EP202) Rabbit Monoclonal Primary Antibody
439R-1 series Mesothelin (EP140) Rabbit Monoclonal Primary Antibody
284M-9 series Microphthalmia Transcription Factor (MiTF) (C5/D5) Mouse Monoclonal
285M-1 series MLH1 (G168-728) Mouse Monoclonal Antibody
286M-1 series MSH2 (G219-1129) Mouse Monoclonal Antibody
287M-1 series MSH6 (44) Mouse Monoclonal Antibody
287R-2 series MSH6 (SP93) Rabbit Monoclonal Antibody
290M-1 series MUC1 (MRQ-17) Mouse Monoclonal Antibody
290R-1 series MUC1 (EP85) Rabbit Monoclonal Primary Antibody
291M-1 series MUC2 (MRQ-18) Mouse Monoclonal Antibody
406M-1 series MUC4 (8G7) Mouse Monoclonal Antibody
292M-9 series MUC5AC (MRQ-19) Mouse Monoclonal Antibody
293M-9 series MUC6 (MRQ-20) Mouse Monoclonal Antibody
358M-1 series MUM1 (MRQ-8) Mouse Monoclonal Antibody
358R-7 series MUM1 (MRQ-43) Rabbit Monoclonal Antibody
358R-1 series MUM1 (EP190) Rabbit Monoclonal Primary Antibody

295A-1 series Myelin Basic Protein Rabbit Polyclonal Antibody
289A-1 series Myeloperoxidase Rabbit Polyclonal Antibody
289R-2 series Myeloperoxidase (EP151) Rabbit Monoclonal Primary Antibody
289R-1 series Myeloperoxidase (SP72) Rabbit Monoclonal Antibody
386R-1 series MyoD1 (EP212) Rabbit Monoclonal Antibody
296M-1 series Myogenin (F5D) Mouse Monoclonal Antibody
296R-1 series Myogenin (EP162) Rabbit Monoclonal Primary Antibody
297A-7 series Myoglobin Rabbit Polyclonal Antibody
298M-1 series Myosin, Smooth Muscle (SMMS-1) Mouse Monoclonal Antibody
298R-1 series Myosin Smooth Muscle (EP166) Rabbit Monoclonal Primary Antibody
443R-1 series Nanog (EP225) Rabbit Monoclonal Primary Antibody
352A-7 series Napsin A Rabbit Polyclonal Antibody
352M-9 series Napsin A (MRQ-60) Mouse Monoclonal Antibody
352R-1 series Napsin A (EP205) Rabbit Monoclonal Primary Antibody
304M-1 series Nerve Growth Factor Receptor (NGFR) (MRQ-21) Mouse Monoclonal Antibody
388M-1 series Nestin (10C2) Mouse Monoclonal Antibody
302M-1 series Neurofilament (2F11) Mouse Monoclonal Antibody
302R-1 series Neurofilament (EP79) Rabbit Monoclonal Primary Antibody
445R-1 series NKX2.2 (EP336) Rabbit Monoclonal Primary Antibody
441R-1 series NKX3.1 (EP356) Rabbit Monoclonal Primary Antibody
306M-2 series NSE (MRQ-55) Mouse Monoclonal Antibody
308M-1 series Oct-2 (MRQ-2) Mouse Monoclonal Antibody
309M-1 series Oct-4 (MRQ-10) Mouse Monoclonal Antibody
387M-1 series Olig2 (211F1.1) Mouse Monoclonal Antibody
387R-1 series Olig2 (EP112) Rabbit Monoclonal Primary Antibody
421M-1 series p21/WAF1 (DCS-60.2) Mouse Monoclonal Antibody
427M-9 series p27/Kip1 (SX53G8) Mouse Monoclonal Antibody
453M-9 series p53 (D07) Mouse Monoclonal Antibody
453R-1 series P53 (SP5) Rabbit Monoclonal Antibody
457M-9 series p57/Kip2 (Kp10) Mouse Monoclonal Antibody
420M-1 series p120 Catenin (MRQ-5) Mouse Monoclonal Antibody
504R-1 series P504s Rabbit Monoclonal Antibody
310M-2 series Parathyroid Hormone (PTH) (MRQ-31) Mouse Monoclonal Antibody
396M-1 series Parvalbumin (2E11) Mouse Monoclonal Antibody
396R-1 series Parvalbumin (EP300) Rabbit Monoclonal Primary Antibody
218M-1 series Parvovirus B19 (R92F6) Mouse Monoclonal Antibody
311R-1 series PAX-2 (EP235) Rabbit Monoclonal Primary Antibody
311A-1 series PAX-2 Rabbit Polyclonal Antibody
312M-1 series PAX-5 (24) Mouse Monoclonal Antibody
312R-2 series PAX-5 (EP156) Rabbit Monoclonal Antibody
312R-1 series PAX-5 (SP34) Rabbit Monoclonal Antibody
363A-1 series PAX-8 Rabbit Polyclonal Antibody
363M-1 series PAX-8 (MRQ-50) Mouse Monoclonal Antibody
315M-9 series PD-1 (NAT105) Mouse Monoclonal Antibody
315R-1 series PD-1 (MRQ-22) Mouse Monoclonal Antibody
438R-1 series PD-L1 (28-8) Rabbit Monoclonal Primary Antibody
438R-2 series PD-L1 (ZR3) Rabbit Monoclonal Antibody
316M-1 series Perforin (MRQ-23) Mouse Monoclonal Antibody

318A-1 series PGP 9.5 Rabbit Polyclonal Antibody
369A-1 series Phosphohistone H3 (PHH3) Rabbit Polyclonal Antibody
422R-1 series PHOX2B (EP312) Rabbit Monoclonal Primary Antibody
321M-1 series PLAP (NB10) Mouse Monoclonal Antibody
321R-2 series PLAP (EP194) Rabbit Monoclonal Primary Antibody
321R-1 series PLAP (SP15) Rabbit Monoclonal Antibody
288R-1 series PMS2 (EPR3947) Rabbit Monoclonal Antibody
288M-1 series PMS2 (MRQ-28) Mouse Monoclonal Antibody
365M-9 series PNL2 (PNL2) Mouse Monoclonal Antibody
219M-1 series Pneumocystis jiroveci (carinii) (3F6) Mouse Monoclonal Antibody
322M-1 series Podoplanin (D2-40) Mouse Monoclonal Antibody
323R-1 series Progesterone Receptor (Y85) Rabbit Monoclonal Antibody
323R-3 series Progesterone Receptor (SP42) Rabbit Monoclonal Antibody
210R-1 series Prolactin (EP193) Rabbit Monoclonal Primary Antibody
210A-1 series Prolactin Rabbit Polyclonal Antibody
324M-1 series PSA (ER-PR8) Mouse Monoclonal Antibody
326M-1 series PSAP (PASE/4LJ) Mouse Monoclonal Antibody
328R-1 series PU.1 (EPR3158Y) Rabbit Monoclonal Antibody
329M-9 series Renal Cell Carcinoma (PN-15) Mouse Monoclonal Antibody
330M-1 series S-100 (4C4.9) Mouse Monoclonal Antibody
449R-1 S100 beta(EP32) Rabbit Monoclonal Primary Antibody
408R-1 series S100A1 (EP184) Rabbit Monoclonal Antibody
376M-9 series S100P (16/f5) Mouse Monoclonal Antibody
385M-1 series SALL4 (6E3) Mouse Monoclonal Antibody
385M-1 series SALL4 (EP299) Rabbit Monoclonal Primary Antibody
384R-1 series SATB2 (EP281) Rabbit Monoclonal Primary Antibody
377M-1 series Smoothelin (R4A) Mouse Monoclonal Antibody
332R-1 series Somatostatin (EP130) Rabbit Monoclonal Primary Antibody
332A-1 series Somatostatin Rabbit Polyclonal Antibody
371R-1 series SOX-2 (SP76) Rabbit Monoclonal Antibody
371R-2 series SOX-2 (EP103) Rabbit Monoclonal Primary Antibody
383A-7 series SOX-10 Rabbit Polyclonal Antibody
383R-1 series SOX-10 (EP268) Rabbit Monoclonal Primary Antibody
382M-1 series SOX-11 (MRQ-58) Mouse Monoclonal Antibody
333R-1 series Spectrin (EP251) Rabbit Monoclonal Antibody
333M-1 series Spectrin (RBC2/3D5) Mouse Monoclonal Antibody
426R-1 series STAT6 (EP325) Rabbit Monoclonal Primary Antibody
394R-1 series Stathmin (SP49) Rabbit Monoclonal Antibody
446R-1 series Steroidogenic Acute Regulatory Protein (STAR) (EP226) Rabbit
351M-1 series SV40 (MRQ-40) Mouse Monoclonal Antibody
336R-1 series Synaptophysin (EP158) Rabbit Monoclonal Primary Antibody
336R-9 series Synaptophysin (MRQ-40) Rabbit Monoclonal Antibody
336A-7 series Synaptophysin Rabbit Polyclonal Antibody
368R-7 series T-bet (MRQ-46) Rabbit Monoclonal Antibody
337M-8 series TAG-72 (B72.3) Mouse Monoclonal Antibody
357M-1 series TCL1 (MRQ-7) Mouse Monoclonal Antibody
338A-7 series TdT Rabbit Polyclonal Antibody
338R-2 series TdT (EP266) Rabbit Monoclonal Primary Antibody

354R-1 series TFE3 (MRQ-37) Rabbit Monoclonal Antibody
 339M-1 series Thrombomodulin (1009) Mouse Monoclonal Antibody
 340M-1 series Thyroglobulin (2H11+6E1) Mouse Monoclonal Antibody
 340M-9 series Thyroglobulin (MRQ-41) Mouse Monoclonal Antibody
 409R-1 series Thyroid Peroxidase (EP159) Rabbit Monoclonal Antibody
 381R-1 series TIA-1 (EP243) Rabbit Monoclonal Antibody
 401M-1 series TLE1 (1f5) Mouse Monoclonal Antibody
 220A-1 series Toxoplasma Gondii Rabbit Polyclonal Antibody
 341M-9 series TRAcP (9C5) Mouse Monoclonal Antibody
 423M-1 series Transgelin (2A10C2) Mouse Monoclonal Antibody
 397A-1 series Treponema Pallidum Rabbit Polyclonal Antibody
 342M-1 series Tryptase (G3) Mouse Monoclonal Antibody
 342R-1 series Tryptase (EP259) Rabbit Monoclonal Primary Antibody
 211R-1 series TSH (EP254) Rabbit Monoclonal Primary Antibody
 211A-1 series TSH Rabbit Polyclonal Antibody
 343M-9 series TTF-1 (8G7G3/1) Mouse Monoclonal Antibody
 343R-1 series TTF-1 (EP229) Rabbit Monoclonal Primary Antibody
 344M-9 series Tyrosinase (T311) Mouse Monoclonal Antibody
 345M-1 series Uroplakin III (AU-1) Mouse Monoclonal Antibody
 345R-1 series Uroplakin III (SP73) Rabbit Monoclonal Antibody
 364M-1 series Varicella Zoster Virus (SG1-1, SG1-SG4, NCP-1 & IE-62) Mouse
 346M-1 series Villin (CWWB1) Mouse Monoclonal Antibody
 346R-1 series Villin (EP163) Rabbit Monoclonal Primary Antibody
 347M-1 series Vimentin (V9) Mouse Monoclonal Antibody
 347R-2 series Vimentin (EP21) Rabbit Monoclonal Primary Antibody
 347R-1 series Vimentin (SP20) Rabbit Monoclonal Antibody
 348M-9 series WT1 (6F-H2) Mouse Monoclonal Antibody
 349M-9 series ZAP-70 (2F3.2) Mouse Monoclonal Antibody
 905H-0 series Cytokeratin 5 (EP1601Y) + Cytokeratin 14 (LL002) Mouse/Rabbit
 902H-0 series Cytokeratin 5 & 6 (D5/16B4) + TTF-1 (8G7G3/1) Mouse Monoclonal
 Antibody

1. POSITIVE CONTROLS SLIDES (Portaobjetos control positivo)

N° Catalog.	Descripción
222S	A-1-Antichymotrypsin Positive Controls Slides
223S	A-1-Antitrypsin Positive Controls Slides
206S	ACTH Positive Controls Slides
201S	Actin, Muscle Specific Positive Controls Slides
202S	Actin, Smooth Muscle Positive Controls Slides
212S	Adenovirus Positive Controls Slides
393S	Adipophilin Positive Control Slides
400S	ALDH1A1 Positive Control Slides
204S	ALK Protein Positive Controls Slides
203S	Alpha-Fetoprotein Positive Controls Slides
200S	Androgen Receptor Positive Controls Slides
221S	Annexin A1 Positive Controls Slides

225S BCA-225 Positive Controls Slides
226S bcl-2 Positive Controls Slides
227S bcl-6 Positive Controls Slides
224S Beta-Catenin Positive Controls Slides
228S BG8, Lewis Y Positive Controls Slides
294S BOB.1 Positive Controls Slides
395S c-Myc Positive Control Slides
403S C3d Positive Controls Slides
404S C4d Positive Controls Slides
325S CA-125 Positive Controls Slides
399S CA19-9 Positive Controls Slides
378S Cadherin-17 Positive Controls Slides
229S Calcitonin Positive Controls Slides
230S Caldesmon Positive Controls Slides
231S Calponin Positive Controls Slides
231S Calponin-1 Positive Controls Slides
232S Calretinin Positive Controls Slides
402S Cathespin K Positive Control Slides
412S Caveolin-1 (2297) Positive Controls Slides
101S CD1a Positive Controls Slides
102S CD2 Positive Controls Slides
103S CD3 Positive Controls Slides
104S CD4 Positive Controls Slides
205S CD5 Positive Controls Slides
107S CD7 Positive Controls Slides
108S CD8 Positive Controls Slides
110S CD10 Positive Controls Slides
113S CD13 Positive Controls Slides
114S CD14 Positive Controls Slides
115S CD15 Positive Controls Slides
116S CD16 Positive Controls Slides
119S CD19 Positive Controls Slides
120S CD20 Positive Controls Slides
121S CD21 Positive Controls Slides
123S CD23 Positive Controls Slides
125S CD25 Positive Controls Slides
130S CD30 Positive Controls Slides
131S CD31 Positive Controls Slides
134S CD34 Positive Controls Slides
135S CD35 Positive Controls Slides
143S CD43 Positive Controls Slides
144S CD44 Positive Controls Slides
145S CD45 (LCA) Positive Controls Slides
146S CD45R Positive Controls Slides
147S CD45RO Positive Controls Slides
156S CD56 Positive Controls Slides
157S CD57 Positive Controls Slides
161S CD61 Positive Controls Slides

263S CD63 Positive Controls Slides
168S CD68 Positive Controls Slides
171S CD71 Positive Controls Slides
174S CD74 Positive Controls Slides
179S CD79a Positive Controls Slides
199S CD99 Positive Controls Slides
437S CD103 Positive Control Slides
105S CD105 (Endoglin) Positive Controls Slides
117S CD117, c-kit Positive Controls Slides
198S CD123 Positive control slides
138S CD138/syndecan-1 Positive Controls Slides
163S CD163 Positive Controls Slides
235S CDX-2 Positive Controls Slides
236S CEA Positive Controls Slides
238S Chromogranin A Positive Controls Slides
424S CITED1 Positive Control Slides
359S Claudin 1 Positive Controls Slides
418S Claudin 7 Positive Control Slides
213S CMV Positive Controls Slides
239S Collagen Type IV Positive Controls Slides
240S COX-2 Positive Controls Slides
241S Cyclin D1 Positive Controls Slides
242S Cystatin A Positive Controls Slides
334S Cytokeratin 34betaE12 Positive Controls Slides
335S Cytokeratin 35betaH11 Positive Controls Slides
452S Cytokeratin CAM 5.2 Positive Controls Slides
300S Cytokeratin OSCAR Positive Controls Slides
305S Cytokeratin 5 Positive Controls Slides
356S Cytokeratin 5&6 Positive Controls Slides
307S Cytokeratin 7 Positive Controls Slides
818S Cytokeratin 8&18 Positive Controls Slides
410S Cytokeratin 10 EP97 Positive Controls Slides
314S Cytokeratin 14 Positive Controls Slides
317S Cytokeratin 17 Positive Controls Slides
319S Cytokeratin 19 Positive Controls Slides
320S Cytokeratin 20 Positive Controls Slides
313S Cytokeratin Cocktail Positive Controls Slides
303S Cytokeratin, HMW Positive Controls Slides
301S Cytokeratin, LMW Positive Controls Slides
243S Desmin Positive Controls Slides
436S Desmoglein 3 Positive Control Slides
244S DOG1 Positive Controls Slides
246S E-Cadherin Positive Controls Slides
414S EGFR Positive Control Slides
247S EMA Positive Controls Slides
248S EP-CAM/Epithelial Specific Antigen Positive Controls Slides
245S Epstein Barr Virus Positive Controls Slides
434S ERG Positive Control Slides

249S Estrogen Receptor Positive Controls Slides
415S EZH2 Positive Control Slides
250S Factor VIII-R Ag. Positive Controls Slides
251S Factor XIIIa Positive Controls Slides
252S Fascin Positive Controls Slides
254S FLI-1 Positive Controls Slides
253S Follicular Dendritic Cell Positive Controls Slides
207S FSH Positive Control Slides
413S GLUT3 Positive Control Slides
405S MUC4 Positive Control Slides
350S FOXP1 Positive Controls Slides
207S FSH Positive Controls Slides
255S Galectin-3 Positive Controls Slides
256S Gastrin Positive Controls Slides
417S GATA1 Positive Control Slides
390S GATA3 Positive Control Slides
257S GCDFP-15 Positive Controls Slides
208S GH Positive Controls Slides
258S Glial Fibrillary Acidic Protein Positive Controls Slides
259S Glucagon Positive Controls Slides
355S GLUT1 Positive Controls Slides
389S Glutamine Synthetase Positive Control Slides
260S Glycophorin A Positive Controls Slides
261S Glypican-3 Positive Controls Slides
262S Granzyme B Positive Controls Slides
234S hCG Positive Controls Slides
398S Heat Shock Protein 27 G3.1 Positive Controls Slides
215S Helicobacter pylori Positive Controls Slides
360S Hemoglobin A Positive Controls Slides
216S Hepatitis B Virus Core Antigen Positive Controls Slides
217S Hepatitis B Virus Surface Antigen Positive Controls Slides
264S Hepatocyte Specific Antigen Positive Controls Slides
237S Her2/Neu (c-erbB-2) Positive Controls Slides
214S Herpes Simplex Virus I and II Positive Controls Slides
375S HGAL Positive Controls Slides
265S HHV-8 Positive Controls Slides
266S Human Placental Lactogen (hPL) Positive Controls Slides
267S IgA Positive Controls Slides
268S IgD Positive Controls Slides
269S IgG Positive Controls Slides
367S IgG4 Positive Controls Slides
270S IgM Positive Controls Slides
433S IMP3 Positive Control Slides
271S Inhibin, alpha Positive Controls Slides
272S INI-1 Positive Controls Slides
273S Insulin Positive Controls Slides
431S Islet-1 Positive Control Slides
274S Kappa Positive Controls Slides

275S Ki-67 Positive Controls Slides
276S Ksp-Cadherin Positive Controls Slides
277S Lambda Positive Controls Slides
392S Langerin Positive Control Slides
442S LEF1 Positive Control Slides
209S LH Positive Controls Slides
464S LIN28 Positive Control Slides
370S LMO2 Positive Controls Slides
278S Lysozyme Positive Controls Slides
279S Macrophage Positive Controls Slides
280S Mammaglobin Positive Controls Slides
281S MART-1 Positive Controls Slides
435S MCM Positive Control Slides
282S Melanoma Positive Controls Slides
366S Melanoma Associated Antigen KBA.62 Positive Controls Slides
365S Melanoma Associated Antigen PNL2 Positive Controls Slides
283S Mesothelial Cell Positive Controls Slides
439S Mesothelin Positive Control Slides
284S Microphthalmia Transcription Factor Positive Controls Slides
285S MLH1 Positive Controls Slides
286S MSH2 Positive Controls Slides
287S MSH6 Positive Controls Slides
290S MUC1 Positive Controls Slides
291S MUC2 Positive Controls Slides
292S MUC5AC Positive Controls Slides
293S MUC6 Positive Controls Slides
358S MUM1 Positive Controls Slides
295S Myelin Basic Protein Positive Controls Slides
289S Myeloperoxidase Positive Controls Slides
386S MyoD1 Positive Control Slides
296S Myogenin Positive Controls Slides
297S Myoglobin Positive Controls Slides
298S Myosin, Smooth Muscle Positive Controls Slides
443S Nanog Positive Control Slides
352S Napsin A Positive Controls Slides
302S Neurofilament Positive Controls Slides
304S Nerve Growth Factor Receptor (NGFR) Positive Controls Slides
388S Nestin Positive Controls Slides
445S NKX2.2 Positive Control Slides
306S NSE Positive Controls Slides
308S Oct-2 Positive Controls Slides
309S Oct-4 Positive Controls Slides
387S Olig2 Positive Controls Slides
416S p16 Positive Controls Slides
421S p21/WAF1 Positive Controls Slides
427S p27/Kip1 Positive Controls Slides
453S p53 Positive Controls Slides
457S p57/Kip2 Positive Controls Slides

463S p63 Positive Controls Slides
420S p120 Catenin Positive Controls Slides
504S P504s Positive Controls Slides
310S Parathyroid Hormone (PTH) Positive Controls Slides
396S Parvalbumin Positive Control Slides
218S Parvovirus B19 Positive Controls Slides
311S PAX-2 Positive Controls Slides
312S PAX-5 Positive Controls Slides
363S PAX-8 Positive Controls Slides
315S PD-1 Positive Controls Slides
316S Perforin Positive Controls Slides
318S PGP 9.5 Positive Controls Slides
369S Phosphohistone H3 Positive Controls Slides
422S PHOX2B Positive Control Slides
321S PLAP Positive Controls Slides
288S PMS2 Positive Controls Slides
219S Pneumocystis jiroveci (carinii) Positive Controls Slides
322S Podoplanin Positive Controls Slides
323S Progesterone Receptor Positive Controls Slides
210S Prolactin Positive Controls Slides
324S PSA Positive Controls Slides
326S PSAP Positive Controls Slides
328S PU.1 Positive Controls Slides
329S Renal Cell Carcinoma Positive Controls Slides
330S S-100 Protein Positive Controls Slides
449S S100 beta Positive Control Slides
408S S100A1 EP184 Positive Controls Slides
376S S100P Positive Controls Slides
385S SALL4 Positive Controls Slides
384S SATB2 Positive Control Slides
377S Smoothelin Positive Controls Slides
332S Somatostatin Positive Controls Slides
371S SOX-2 Positive Controls Slides
383S SOX-10 Positive Controls Slides
333S Spectrin Positive Controls Slides
426S STAT6 Positive Control Slides
394S Stathmin Positive Control Slides
446S Steroidogenic Acute Regulatory Protein (STAR) Positive Control Slides
351S SV40 Positive Controls Slides
336S Synaptophysin Positive Controls Slides
368S T-bet Positive Controls Slides
337S TAG-72 Positive Controls Slides
357S TCL1 Positive Controls Slides
338S TdT Positive Controls Slides
354S TFE3 Positive Controls Slides
339S Thrombomodulin Positive Controls Slides
340S Thyroglobulin Positive Controls Slides
409S Thyroid Peroxidase EP159 Positive Controls Slides

- 381S TIA-1 Positive Control Slides
- 220S Toxoplasma Gondii Positive Controls Slides
- 341S TRAcP Positive Controls Slides
- 423S Transgelin Positive Control Slides
- 397S Treponema Pallidum Positive Controls Slides
- 342S Tryptase Positive Controls Slides
- 211S TSH Positive Controls Slides
- 343S TTF-1 Positive Controls Slides
- 344S Tyrosinase Positive Controls Slides
- 345S Uroplakin III Positive Controls Slides
- 364S Varicella Zoster Virus Positive Controls Slides
- 346S Villin Positive Controls Slides
- 347S Vimentin Positive Controls Slides
- 348S WT1 Positive Controls Slides
- 349S ZAP-70 Positive Controls Slides

1. GENERAL PURPOSE REAGENTS (Reactivos auxiliares)

Nº Catalog.	Periodo de vida útil y condiciones de conservación:
920P series Trilogy™	36 meses conservado entre 20 – 26 °C
921P series Declere®	36 meses conservado entre 20 – 26 °C
925B series Peroxide Block	36 meses conservado entre 20 – 26 °C
926B series PeroxFree™ Block	36 meses conservado entre 20 – 26 °C
927B series Background Block	36 meses conservado entre 2 – 8 °C
928B series Avidin-Biotin Blocking Reagents	36 meses conservado entre 2 – 8 °C
929B series IHC Wash Buffer	36 meses conservado entre 20 – 26 °C
930B series Hematoxylin	24 meses conservado entre 20 – 26 °C
931B series Permanent Acqueous Mounting Medium	36 meses conservado entre 20 – 26 °C
932B series Mouse Negative Control Serum,	36 meses conservado entre 2 – 8 °C
933B series Rabbit Negative Control Serum	36 meses conservado entre 2 – 8 °C
934B series 20X PBS IHC Wash Buffer + Tween® 20	36 meses conservado entre 2 – 8 °C
935B series 20X TBS IHC Wash Buffer + Tween®20	36 meses conservado entre 2 – 8 °C
936B series Antibody Diluent	36 meses conservado entre 2 – 8 °C
937B series Protease	24 meses conservado entre 2 – 8 °C
938B series Diamond: Antibody Diluent	36 meses conservado entre 2 – 8 °C
951D series CytoScan™ HRP Detection System	24 meses conservado entre 2 – 8 °C
952D series CytoScan™ Alk Phos Detection System	24 meses conservado entre 2 – 8 °C
954D series HiDef Detection™ HRP Polymer System	24 meses conservado entre 2 – 8 °C
955D series CytoScan™ Alk Phos Detection System w/ Perm Magenta	24 meses conservado entre 2 – 8 °C
956D series PolyScan™ HRP / DAB Detection	

System	24 meses conservado entre 2 – 8 °C
957D series DAB Substrate Kit	24 meses conservado entre 2 – 8 °C
958D series AEC	24 meses conservado entre 2 – 8 °C
959D series Permanent Magenta Chromogen Kit	24 meses conservado entre 2 – 8 °C
960D series Permanent Red Chromogen Kit	24 meses conservado entre 2 – 8 °C
962D series HiDef Detection™ Alk Phos Polymer	24 meses conservado entre 2 – 8 °C
System	
963D series Firenze Detection HRP Polymer	24 meses conservado entre 2 – 8 °C
System	

Expediente N° 1-47-3110-5247/16-9

Digitally signed by LEDE Roberto Luis
 Date: 2018.04.17 09:24:52 ART
 Location: Ciudad Autónoma de Buenos Aires

Roberto Luis Ledo
 SubAdministrador
 Administración Nacional de Medicamentos, Alimentos y Tecnología
 Médica