



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A. N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

DISPOSICIÓN N° 7817

BUENOS AIRES,

03 DIC 2010

VISTO el expediente N° 1-47-22430/08-6 del Registro de la Administración Nacional de Medicamentos Alimentos y Tecnología Médica y,

CONSIDERANDO:

Que por los presentes actuados la firma INVITROGEN ARGENTINA S.A. solicita autorización para la venta a laboratorios de análisis clínicos del Reactivo de diagnóstico de uso "in vitro" denominado Anticuerpos Monoclonales / PARA TINTAR EL ANTIGENO EN SECCIONES DE TEJIDO CONGELADAS E INFILTRADAS EN PARAFINA, que se detallan en el Anexo I.

Que a fojas 298 a 300 consta el informe técnico producido por el Servicio de Productos para Diagnóstico que establecen que los productos reúnen las condiciones de aptitud requeridas para su autorización.

Que la Dirección de Tecnología Médica ha tomado la intervención de su competencia.

Que se ha dado cumplimiento a los términos que establece la Ley N° 16.463, Resolución Ministerial N° 145/98 y Disposición ANMAT N° 2674/99.

Que la presente se dicta en virtud de las facultades conferidas por el Artículo 8° inciso 11) del Decreto 1490/92 y Decreto N° 425/10.

Por ello;

EL INTERVENTOR DE LA ADMINISTRACIÓN NACIONAL DE
MEDICAMENTOS, ALIMENTOS Y TECNOLOGÍA MÉDICA

DISPONE:



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N.M.A.T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

DISPOSICIÓN Nº **7817**

ARTÍCULO 1º.- Autorízase la venta a laboratorios de análisis clínicos del producto de diagnóstico para uso in Vitro denominado Anticuerpos Monoclonales / PARA TINTAR EL ANTIGENO EN SECCIONES DE TEJIDO CONGELADAS E INFILTRADAS EN PARAFINA, el que será elaborado por Invitrogen Corporation, 542 Flynn Road, Camarillo, CA 93012 (U.S.A.) e importado terminado por la firma INVITROGEN ARGENTINA S.A. en envases que se detallan en el Anexo I, con una vida útil de DIECIOCHO (18) MESES, conservado entre 2-8°C, y que la composición se detalla a fojas 69 a 73.

ARTICULO 2º.- Acéptense los proyectos de rótulos y Manual de Instrucciones a fojas 220 a 228, 232 a 234 y 238 a 240 debiendo constar en los mismos que la fecha de vencimiento es la declarada por el elaborador impreso en los rótulos de cada partida.

ARTÍCULO 3º.- Extiéndase el Certificado correspondiente.

ARTICULO 4º.- LA ADMINISTRACIÓN NACIONAL DE MEDICAMENTOS, ALIMENTOS Y TECNOLOGÍA MEDICA se reserva el derecho de reexaminar los métodos de control, estabilidad y elaboración cuando las circunstancias así lo determinen.

ARTÍCULO 5º.- Anótese, gírese al Departamento de Registro a sus efectos, por Mesa de Entradas notifíquese al interesado y hágasele entrega de la copia de la presente Disposición junto con la copia de los proyectos de rótulos, manual de instrucciones y el certificado correspondiente. Remítase una copia de la presente Disposición y Certificado a la Dirección de Tecnología Médica. Cumplido, Archívese PERMANENTE.-

EXPEDIENTE Nº 1-47-22430/08-6

DISPOSICIÓN Nº: **7817**

Dr. OTTO A. ORSINGER
SUB-INTERVENTOR
A.N.M.A.T.



7817

ANEXO I

Expediente N° 1-47-22430/08-6

PRODUCTO/USO: Anticuerpos Monoclonales / PARA TINTAR EL ANTIGENO EN SECCIONES DE TEJIDO CONGELADAS E INFILTRADAS EN PARAFINA.

N° de catalogo	Description - Volumen
180270	MOUSE ANTI- EPITHELIAL GLYCOPR 1 ML
180133	MOUSE ANTI- METALLOTHIONEIN (C 1 ML
187328	MOUSE ANTI- VEGF (CLONE VG1) C 0.5 ML
81268	MOUSE ANTI-ACTIN (CLONE HHF35) 6 ML
80053	MOUSE ANTI-ACTIN (CLONE ZCA34) 6 ML
81053	MOUSE ANTI-ACTIN (CLONE ZCA34) 6 ML
180177Z	MOUSE ANTI-ACTIN (SARCOMERIC) 1 ML
081177	MOUSE ANTI-ACTIN (SARCOMERIC) 6 ML
180106	MOUSE ANTI-ACTIN (SMOOTH MUSCL 1 ML
84106	MOUSE ANTI-ACTIN (SMOOTH MUSCL 10 ML
180003	MOUSE ANTI-ALPHA FETOPROTEIN (1 ML
080003	MOUSE ANTI-ALPHA FETOPROTEIN (6 ML
081003	MOUSE ANTI-ALPHA FETOPROTEIN (6 ML
180225	MOUSE ANTI-ALPHA-CATENIN (CLON 1 ML
181215	MOUSE ANTI-ALPHA-SYNUCLEIN (CL 1 ML
081215	MOUSE ANTI-ALPHA-SYNUCLEIN (CL 6 ML
180092	MOUSE ANTI-ALPHA-TUBULIN (CLON 1 ML
81292	MOUSE ANTI-ANDROGEN RECEPTOR (6 ML
180088	MOUSE ANTI-B CELL (CD20) (CLON 1 ML
080088	MOUSE ANTI-B CELL (CD20) (CLON 6 ML
081088	MOUSE ANTI-B CELL (CD20) (CLON 6 ML
180218	MOUSE ANTI-BAX (CLONE 2D2) CON 1 ML
180275	MOUSE ANTI-BCL-10 (CLONE 151) 1 ML
081275	MOUSE ANTI-BCL-10 (CLONE 151) 6 ML
084193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 10 ML
080193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 6 ML
081193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 6 ML
81426	MOUSE ANTI-BCL-6 (CLONE PIF6) 6 ML
180217	MOUSE ANTI-BCL-XL (CLONE 2H12) 1 ML
180226	MOUSE ANTI-BETA-CATENIN (CLONE 1 ML
180093	MOUSE ANTI-BETA-TUBULIN (CLONE 1 ML
187242	MOUSE ANTI-CA 125 ANTIGEN (CLO 0.5 ML
187260	MOUSE ANTI-CA 125 ANTIGEN (CLO 0.5 ML
081260	MOUSE ANTI-CA 125 ANTIGEN (CLO 6 ML
187239	MOUSE ANTI-CA 15-3 ANTIGEN (CL 0.5 ML
187240	MOUSE ANTI-CA 19-9 ANTIGEN (CL 0.5 ML
187265	MOUSE ANTI-CA 19-9 ANTIGEN (CL 0.5 ML



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N.M.A.T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

7817

081265 MOUSE ANTI-CA 19-9 ANTIGEN (CL 6 ML
180291 MOUSE ANTI-CALRETININ (CLONE Z 1 ML
081291 MOUSE ANTI-CALRETININ (CLONE Z 6 ML
081287 MOUSE ANTI-CD10 (CLONE 56C6) 2 6 ML
182287 MOUSE ANTI-CD10 (CLONE 56C6) C 0.2 ML
182409 MOUSE ANTI-CD13 (CLONE 38C12) 0.2 ML
081409 MOUSE ANTI-CD13 (CLONE 38C12) 6 ML
180121 MOUSE ANTI-CD14 (CLONE RPA-M1) 1 ML
081122 MOUSE ANTI-CD15 (CLONE MY1) 2N 6 ML
180122 MOUSE ANTI-CD15 (CLONE MY1) CO 1 ML
080122 MOUSE ANTI-CD15 (CLONE MY1) PR 6 ML
180380 MOUSE ANTI-CD18 (CLONE MEM-48) 1 ML
81269 MOUSE ANTI-CD3 (T CELL) (CLONE 6 ML
187155 MOUSE ANTI-CD30 (KI-1) (CLONE 0.5 ML
080155 MOUSE ANTI-CD30 (KI-1) (CLONE 6 ML
081155 MOUSE ANTI-CD30 (KI-1) (CLONE 6 ML
180155 MOUSE ANTI-CD30 1 ML
081425 MOUSE ANTI-CD31 (CLONE 1A10) 2 6 ML
182425 MOUSE ANTI-CD31 (CLONE 1A10) C 0.2 ML
180375 MOUSE ANTI-CD31 (CLONE MEM-05) 1 ML
180227 MOUSE ANTI-CD34 (CLONE BI-3C5) 1 ML
084227 MOUSE ANTI-CD34 (CLONE BI-3C5) 10 ML
081282 MOUSE ANTI-CD4 (CLONE 1F6) 2ND 6 ML
182282 MOUSE ANTI-CD4 (CLONE 1F6) CON 0.2 ML
187250 MOUSE ANTI-CD43 (CLONE DF-T1) 0.5 ML
080186 MOUSE ANTI-CD43 (CLONE DF-T1) 6 ML
080184 MOUSE ANTI-CD44 (CLONE DF1485) 6 ML
180367 MOUSE ANTI-CD45 (LCA) (CLONES 1 ML
081367 MOUSE ANTI-CD45 (LCA) (CLONES 6 ML
180365 MOUSE ANTI-CD45RO (CLONE UCHL- 1 ML
081365 MOUSE ANTI-CD45RO (CLONE UCHL- 6 ML
081283 MOUSE ANTI-CD5 (CLONE 4C7) 2ND 6 ML
180173Z MOUSE ANTI-CD54 (I-CAM) (CLONE 1 ML
180152 MOUSE ANTI-CD56 (SMALL CELL LU 1 ML
080152 MOUSE ANTI-CD56 (SMALL CELL LU 6 ML
081152 MOUSE ANTI-CD56 (SMALL CELL LU 6 ML
180167Z MOUSE ANTI-CD57 (CLONE NK-1) C 1 ML
080167 MOUSE ANTI-CD57 (CLONE NK-1) P 6 ML
180124 MOUSE ANTI-CD61 (CLONE Y2/51) 1 ML
187300 MOUSE ANTI-CD63 (CLONE FC-5.01 0.5 ML
080125 MOUSE ANTI-CD68 (CLONE KP1) PR 6 ML
081289 MOUSE ANTI-CD8 (CLONE 1A5) 2ND 6 ML
182289 MOUSE ANTI-CD8 (CLONE 1A5) CON 0.2 ML
180235 MOUSE ANTI-CD99 (EWING'S SARCO 1 ML
081235 MOUSE ANTI-CD99 (EWING'S SARCO 6 ML
180164 MOUSE ANTI-CDW75 (B CELL) (CLO 1 ML



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

7817

180057	MOUSE ANTI-CEA (CLONE COL-1) C 1 ML
081013	MOUSE ANTI-CEA (CLONE ZC23) 2N 6 ML
180013	MOUSE ANTI-CEA (CLONE ZC23) CO 1 ML
81384	MOUSE ANTI-C-KIT (CLONE 2E4) 2 6 ML
180384	MOUSE ANTI-C-KIT (CLONE 2E4) C 1 ML
187363	MOUSE ANTI-CLAUDIN 2 (CLONE 12 0.5 ML
187364	MOUSE ANTI-CLAUDIN 5 (CLONE 4C 0.5 ML
187391	MOUSE ANTI-CLAUDIN 7 (CLONE 5D 0.5 ML
187341	MOUSE ANTI-CLAUDIN-4 (CLONE 3E 0.5 ML
180176Z	MOUSE ANTI-C-MYC (CLONE 9E10) 1 ML
180128	MOUSE ANTI-COLLAGEN IV (CLONE 1 ML
187381	MOUSE ANTI-COX-1 (CLONE COX111 0.5 ML
187379	MOUSE ANTI-COX-2 (COX-229) CON 0.5 ML
180220	MOUSE ANTI-CYCLIN D1 (CLONE AM 1 ML
081220	MOUSE ANTI-CYCLIN D1 (CLONE AM 6 ML
180404	MOUSE ANTI-CYCLIN E (CLONE HE1 1 ML
180153	MOUSE ANTI-CYTOKERATIN (ACIDIC 1 ML
080153	MOUSE ANTI-CYTOKERATIN (ACIDIC 6 ML
180154	MOUSE ANTI-CYTOKERATIN (BASIC) 1 ML
081212	MOUSE ANTI-CYTOKERATIN (HMW) (6 ML
180132	MOUSE ANTI-CYTOKERATIN (PAN) (1 ML
084132	MOUSE ANTI-CYTOKERATIN (PAN) (10 ML
080132	MOUSE ANTI-CYTOKERATIN (PAN) (6 ML
081132	MOUSE ANTI-CYTOKERATIN (PAN) (6 ML
180158Z	MOUSE ANTI-CYTOKERATIN 18 (CLO 1 ML
080190	MOUSE ANTI-CYTOKERATIN 19 (CLO 6 ML
182200	MOUSE ANTI-CYTOKERATIN 20 (CLO 0.2 ML
081200	MOUSE ANTI-CYTOKERATIN 20 (CLO 6 ML
180267	MOUSE ANTI-CYTOKERATIN 5/6 (CL 1 ML
081267	MOUSE ANTI-CYTOKERATIN 5/6 (CL 6 ML
180234	MOUSE ANTI-CYTOKERATIN 7 (CLON 1 ML
081234	MOUSE ANTI-CYTOKERATIN 7 (CLON 6 ML
180185Z	MOUSE ANTI-CYTOKERATIN 8 (CLON 1 ML
180213	MOUSE ANTI-CYTOKERATIN 8/18 (C 1 ML
081213	MOUSE ANTI-CYTOKERATIN 8/18 (C 6 ML
280001	MOUSE ANTI-CYTOKERATIN PAN (CL 2.5 ML
182410	MOUSE ANTI-D2-40 CONCENTRATE 0.2 ML
180016	MOUSE ANTI-DESMIN (CLONE ZC18) 1 ML
080016	MOUSE ANTI-DESMIN (CLONE ZC18) 6 ML
081016	MOUSE ANTI-DESMIN (CLONE ZC18) 6 ML
180223	MOUSE ANTI-E-CADHERIN (CLONE 4 1 ML
081223	MOUSE ANTI-E-CADHERIN (CLONE 4 6 ML
81222	MOUSE ANTI-E-CADHERIN (CLONE H 6 ML
084205	MOUSE ANTI-EGFR (CLONE 31G7) V 10 ML
180017	MOUSE ANTI-EMA (CLONE ZCE113) 1 ML
080017	MOUSE ANTI-EMA (CLONE ZCE113) 6 ML

07



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A. N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

7817

81017	MOUSE ANTI-EMA (CLONE ZCE113) 6 ML
182394	MOUSE ANTI-EPHB4 RECEPTOR (CLO 0.2 ML
182357	MOUSE ANTI-EZRIN (CLONE 3C12) 0.2 ML
80061	MOUSE ANTI-FACTOR VIII RELATED 6 ML
81061	MOUSE ANTI-FACTOR VIII RELATED 6 ML
180475	MOUSE ANTI-GALECTIN-3 (ZG001) CONC 1 ML
081475	MOUSE ANTI-GALECTIN-3 (ZG001) PRED 6 ML
187284	MOUSE ANTI-GCDFP-15 (CLONE D6) 0.5 ML
081021	MOUSE ANTI-GFAP (CLONE ZCG29) 6 ML
182306	MOUSE ANTI-GMP-1 (SENTRIN-1- S 0.2 ML
187326	MOUSE ANTI-GPCR-K2 (CLONE 5D5) 0.5 ML
180005	MOUSE ANTI-HCG- BETA CHAIN SPE 1 ML
080005	MOUSE ANTI-HCG- BETA CHAIN SPE 6 ML
187329	MOUSE ANTI-HER4 (CLONE HFR1) C 0.5 ML
180168	MOUSE ANTI-HLA-DR (CLONE LN-3) 1 ML
180031	MOUSE ANTI-HUMAN KAPPA LIGHT C 1 ML
081031	MOUSE ANTI-HUMAN KAPPA LIGHT C 6 ML
180032	MOUSE ANTI-HUMAN LAMBDA LIGHT 1 ML
081032	MOUSE ANTI-HUMAN LAMBDA LIGHT 6 ML
180066	MOUSE ANTI-INSULIN (CLONE Z006 1 ML
080066	MOUSE ANTI-INSULIN (CLONE Z006 6 ML
180192Z	MOUSE ANTI-KI-67 (CLONE 7B11) 1 ML
084192	MOUSE ANTI-KI-67 (CLONE 7B11) 10 ML
081156	MOUSE ANTI-KI-67 2NDGEN PREDIL 6 ML
080156	MOUSE ANTI-KI-67 PREDILUTE 6 ML
180473	MOUSE ANTI-LAT 1 ML
180037	MOUSE ANTI-LUTEINIZING HORMONE 1 ML
081410	MOUSE ANTI-LYMPHATIC ENDOTHELI 6 ML
180165	MOUSE ANTI-MACROPHAGES (CLONE 1 ML
180376	MOUSE ANTI-MAP KINASE (ERK1+ER 1 ML
187259	MOUSE ANTI-MART-1 (MELAN-A) (C 0.5 ML
187263	MOUSE ANTI-MART-1 (MELAN-A) CO 0.5 ML
182403	MOUSE ANTI-MDM2 (CLONE IF2) CO 0.2 ML
180485	MOUSE ANTI-MELANOMA PNL2 1 ML
182050	MOUSE ANTI-MELANOSOME (CLONE H 0.2 ML
084050	MOUSE ANTI-MELANOSOME (CLONE H 10 ML
080050	MOUSE ANTI-MELANOSOME (CLONE H 6 ML
081050	MOUSE ANTI-MELANOSOME (CLONE H 6 ML
187337	MOUSE ANTI-MGMT (CLONE MT23.2) 0.5 ML
180369	MOUSE ANTI-MITF (CLONE C5+D5) 1 ML
081369	MOUSE ANTI-MITF (CLONE C5+D5) 6 ML
081342	MOUSE ANTI-MLH1 (CLONE 14) 2ND 6 ML
182342	MOUSE ANTI-MLH1 (CLONE 14) CON 0.2 ML
081445	MOUSE ANTI-MLH1 (ZY001) 2ND GEN PRED 6 ML
180445	MOUSE ANTI-MLH1 (ZY001) CONCENTRATE 1 ML
081343	MOUSE ANTI-MSH2 (CLONE FE11) 2 6 ML

5

A



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

7817

182343	MOUSE ANTI-MSH2 (CLONE FE11) C 0.2 ML
081443	MOUSE ANTI-MSH6 (2D4B5) 2ND GEN 6 ML
180443	MOUSE ANTI-MSH6 (2D4B5) CONCENTRATE 1 ML
182298	MOUSE ANTI-MUC1 (CLONE VU-4-H5) 0.2 ML
182299	MOUSE ANTI-MUC2 (CLONE CCP58) 0.2 ML
182322	MOUSE ANTI-MUC4 (CLONE 1G8) CO 0.2 ML
182261	MOUSE ANTI-MUC5AC (CLONE 45M1) 0.2 ML
080160	MOUSE ANTI-MULTI-DRUG RESEISTA 6 ML
187243	MOUSE ANTI-MULTI-DRUG RESISTAN 0.5 ML
187244	MOUSE ANTI-MULTI-DRUG RESISTAN 0.5 ML
180160Z	MOUSE ANTI-MULTI-DRUG RESISTAN 1 ML
180073	MOUSE ANTI-MYOGLOBIN (CLONE Z0) 1 ML
180105	MOUSE ANTI-MYOSIN (SKELETAL MU) 1 ML
180224	MOUSE ANTI-N-CADHERIN (CLONE 3) 1 ML
080139	MOUSE ANTI-NEUROFILAMENT (NF-M) 6 ML
180171Z	MOUSE ANTI-NEUROFILAMENT (PAN) 1 ML
187397	MOUSE ANTI-NNOS (CLONE 3G6B10) 0.5 ML
180196	MOUSE ANTI-NSE (CLONE NSE-1G4) 1 ML
080196	MOUSE ANTI-NSE (CLONE NSE-1G4) 6 ML
081196	MOUSE ANTI-NSE (CLONE NSE-1G4) 6 ML
187288	MOUSE ANTI-NUCLEOPHOSMIN (NPM) 0.5 ML
182359	MOUSE ANTI-NY-ESO-1 (CLONE E97) 0.2 ML
180471	MOUSE ANTI-OCT2 (ZO0001) CONC 1 ML
081471	MOUSE ANTI-OCT2 (ZO0001) PREDILUTE 6 ML
81401	MOUSE ANTI-P21 (CLONE EA10) 2N 6 ML
180401	MOUSE ANTI-P21 (CLONE EA10) CO 1 ML
182370	MOUSE ANTI-P27 (CLONE 57) CONC 0.2 ML
180437	MOUSE ANTI-P34 CONCENTRATE 1 ML
081437	MOUSE ANTI-P34 PREDILUTE 6 ML
081129	MOUSE ANTI-P53 (CLONE BP53.12) 6 ML
187251	MOUSE ANTI-P53 (CLONE DO7) CON 0.5 ML
180414	MOUSE ANTI-P63 (4B1E2) CONCENTRATE 1 ML
081414	MOUSE ANTI-P63 (4B1E2) PREDILUTE 6 ML
180406	MOUSE ANTI-P73 (CLONE P73-9F7) 1 ML
180472	MOUSE ANTI-PAX5 (ZP007) CONCENTRATE 1 ML
081472	MOUSE ANTI-PAX5(ZP007) PREDILUTE 6 ML
180407	MOUSE ANTI-PAXILLIN (CLONE ZO3) 1 ML
081110	MOUSE ANTI-PCNA (CLONE PC10) 2 6 ML
180110	MOUSE ANTI-PCNA (CLONE PC10) C 1 ML
182400	MOUSE ANTI-PHOSPHO-STAT-1 (TYR) 0.2 ML
182262	MOUSE ANTI-PNRA (PROXIMAL NEPH) 0.2 ML
180162Z	MOUSE ANTI-PS2 (CLONE BC04) CO 1 ML
180159Z	MOUSE ANTI-RETINOBLASTOMA GENE 1 ML
180307	MOUSE ANTI-SKP2 (CLONE 1G12E9) 1 ML
187248	MOUSE ANTI-TAG-72 (CA 72.4 ANT) 0.5 ML
080148	MOUSE ANTI-TAU (CLONE T14) PRE 6 ML

57



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N. M. A. T.

"2010 - Año del Bicentenario de la Revolución de Mayo"

187319	MOUSE ANTI-THYMIDYLATE SYNTHAS 0.5 ML
180399	MOUSE ANTI-THYMIDYLATE SYNTHAS 1 ML
180405	MOUSE ANTI-THYMIDYLATE SYNTHAS 1 ML
081399	MOUSE ANTI-THYMIDYLATE SYNTHAS 6 ML
081405	MOUSE ANTI-THYMIDYLATE SYNTHAS 6 ML
180199	MOUSE ANTI-TRACP (HAIRY CELLS) 1 ML
081199	MOUSE ANTI-TRACP (HAIRY CELLS) 6 ML
180221	MOUSE ANTI-TTF-1 (CLONE 8G7G3/ 1 ML
081221	MOUSE ANTI-TTF-1 (CLONE 8G7G3/ 6 ML
187352	MOUSE ANTI-TYROSINASE (CLONE T 0.5 ML
080147	MOUSE ANTI-UBIQUITIN (CLONE UB 6 ML
180052	MOUSE ANTI-VIMENTIN (CLONE V9) 1 ML
080052	MOUSE ANTI-VIMENTIN (CLONE V9) 6 ML
081052	MOUSE ANTI-VIMENTIN (CLONE V9) 6 ML
182416	MOUSE ANTI-ZAP-70 (CLONE 2F3.2) 0.2 ML
180416	MOUSE ANTI-ZAP-70 (CLONE 2F3.2) 1 ML
081416	MOUSE ANTI-ZAP-70 (CLONE2F3.2) 6 ML

7817

Dr. OTTO A. ORSINGER
SUB-INTERVENTOR
A.N.M.A.T.



CERTIFICADO DE AUTORIZACIÓN DE VENTA DE
REACTIVOS DE DIAGNOSTICO

Expediente nº 1-47-22430/08-6

Se autoriza a la firma INVITROGEN ARGENTINA S.A. a comercializar el Reactivo de diagnóstico de uso in vitro denominado:

Anticuerpos Monoclonales / PARA TINTAR EL ANTIGENO EN SECCIONES DE TEJIDO CONGELADAS E INFILTRADAS EN PARAFINA.

Nº de catalogo	Description - Volumen
180270	MOUSE ANTI- EPITHELIAL GLYCOPR 1 ML
180133	MOUSE ANTI- METALLOTHIONEIN (C 1 ML
187328	MOUSE ANTI- VEGF (CLONE VG1) C 0.5 ML
81268	MOUSE ANTI-ACTIN (CLONE HHF35) 6 ML
80053	MOUSE ANTI-ACTIN (CLONE ZCA34) 6 ML
81053	MOUSE ANTI-ACTIN (CLONE ZCA34) 6 ML
180177Z	MOUSE ANTI-ACTIN (SARCOMERIC) 1 ML
081177	MOUSE ANTI-ACTIN (SARCOMERIC) 6 ML
180106	MOUSE ANTI-ACTIN (SMOOTH MUSCL 1 ML
84106	MOUSE ANTI-ACTIN (SMOOTH MUSCL 10 ML
180003	MOUSE ANTI-ALPHA FETOPROTEIN (1 ML
080003	MOUSE ANTI-ALPHA FETOPROTEIN (6 ML
081003	MOUSE ANTI-ALPHA FETOPROTEIN (6 ML
180225	MOUSE ANTI-ALPHA-CATENIN (CLON 1 ML
181215	MOUSE ANTI-ALPHA-SYNUCLEIN (CL 1 ML
081215	MOUSE ANTI-ALPHA-SYNUCLEIN (CL 6 ML
180092	MOUSE ANTI-ALPHA-TUBULIN (CLON 1 ML
81292	MOUSE ANTI-ANDROGEN RECEPTOR (6 ML
180088	MOUSE ANTI-B CELL (CD20) (CLON 1 ML
080088	MOUSE ANTI-B CELL (CD20) (CLON 6 ML
081088	MOUSE ANTI-B CELL (CD20) (CLON 6 ML
180218	MOUSE ANTI-BAX (CLONE 2D2) CON 1 ML
180275	MOUSE ANTI-BCL-10 (CLONE 151) 1 ML
081275	MOUSE ANTI-BCL-10 (CLONE 151) 6 ML
084193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 10 ML
080193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 6 ML
081193	MOUSE ANTI-BCL-2 (CLONE BCL-2- 6 ML
81426	MOUSE ANTI-BCL-6 (CLONE P1F6) 6 ML
180217	MOUSE ANTI-BCL-XL (CLONE 2H12) 1 ML
180226	MOUSE ANTI-BETA-CATENIN (CLONE 1 ML

180093 MOUSE ANTI-BETA-TUBULIN (CLONE 1 ML
 187242 MOUSE ANTI-CA 125 ANTIGEN (CLO 0.5 ML
 187260 MOUSE ANTI-CA 125 ANTIGEN (CLO 0.5 ML
 081260 MOUSE ANTI-CA 125 ANTIGEN (CLO 6 ML
 187239 MOUSE ANTI-CA 15-3 ANTIGEN (CL 0.5 ML
 187240 MOUSE ANTI-CA 19-9 ANTIGEN (CL 0.5 ML
 187265 MOUSE ANTI-CA 19-9 ANTIGEN (CL 0.5 ML
 081265 MOUSE ANTI-CA 19-9 ANTIGEN (CL 6 ML
 180291 MOUSE ANTI-CALRETININ (CLONE Z 1 ML
 081291 MOUSE ANTI-CALRETININ (CLONE Z 6 ML
 081287 MOUSE ANTI-CD10 (CLONE 56C6) 2 6 ML
 182287 MOUSE ANTI-CD10 (CLONE 56C6) C 0.2 ML
 182409 MOUSE ANTI-CD13 (CLONE 38C12) 0.2 ML
 081409 MOUSE ANTI-CD13 (CLONE 38C12) 6 ML
 180121 MOUSE ANTI-CD14 (CLONE RPA-M1) 1 ML
 081122 MOUSE ANTI-CD15 (CLONE MY1) 2N 6 ML
 180122 MOUSE ANTI-CD15 (CLONE MY1) CO 1 ML
 080122 MOUSE ANTI-CD15 (CLONE MY1) PR 6 ML
 180380 MOUSE ANTI-CD18 (CLONE MEM-48) 1 ML
 81269 MOUSE ANTI-CD3 (T CELL) (CLONE 6 ML
 187155 MOUSE ANTI-CD30 (KI-1) (CLONE 0.5 ML
 080155 MOUSE ANTI-CD30 (KI-1) (CLONE 6 ML
 081155 MOUSE ANTI-CD30 (KI-1) (CLONE 6 ML
 180155 MOUSE ANTI-CD30 1 ML
 081425 MOUSE ANTI-CD31 (CLONE 1A10) 2 6 ML
 182425 MOUSE ANTI-CD31 (CLONE 1A10) C 0.2 ML
 180375 MOUSE ANTI-CD31 (CLONE MEM-05) 1 ML
 180227 MOUSE ANTI-CD34 (CLONE BI-3C5) 1 ML
 084227 MOUSE ANTI-CD34 (CLONE BI-3C5) 10 ML
 081282 MOUSE ANTI-CD4 (CLONE 1F6) 2ND 6 ML
 182282 MOUSE ANTI-CD4 (CLONE 1F6) CON 0.2 ML
 187250 MOUSE ANTI-CD43 (CLONE DF-T1) 0.5 ML
 080186 MOUSE ANTI-CD43 (CLONE DF-T1) 6 ML
 080184 MOUSE ANTI-CD44 (CLONE DF1485) 6 ML
 180367 MOUSE ANTI-CD45 (LCA) (CLONES 1 ML
 081367 MOUSE ANTI-CD45 (LCA) (CLONES 6 ML
 180365 MOUSE ANTI-CD45RO (CLONE UCHL- 1 ML
 081365 MOUSE ANTI-CD45RO (CLONE UCHL- 6 ML
 081283 MOUSE ANTI-CD5 (CLONE 4C7) 2ND 6 ML
 180173Z MOUSE ANTI-CD54 (I-CAM) (CLONE 1 ML
 180152 MOUSE ANTI-CD56 (SMALL CELL LU 1 ML
 080152 MOUSE ANTI-CD56 (SMALL CELL LU 6 ML
 081152 MOUSE ANTI-CD56 (SMALL CELL LU 6 ML
 180167Z MOUSE ANTI-CD57 (CLONE NK-1) C 1 ML



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A. N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

080167	MOUSE ANTI-CD57 (CLONE NK-1) P 6 ML
180124	MOUSE ANTI-CD61 (CLONE Y2/51) 1 ML
187300	MOUSE ANTI-CD63 (CLONE FC-5.01 0.5 ML
080125	MOUSE ANTI-CD68 (CLONE KP1) PR 6 ML
081289	MOUSE ANTI-CD8 (CLONE 1A5) 2ND 6 ML
182289	MOUSE ANTI-CD8 (CLONE 1A5) CON 0.2 ML
180235	MOUSE ANTI-CD99 (EWING'S SARCO 1 ML
081235	MOUSE ANTI-CD99 (EWING'S SARCO 6 ML
180164	MOUSE ANTI-CDW75 (B CELL) (CLO 1 ML
180057	MOUSE ANTI-CEA (CLONE COL-1) C 1 ML
081013	MOUSE ANTI-CEA (CLONE ZC23) 2N 6 ML
180013	MOUSE ANTI-CEA (CLONE ZC23) CO 1 ML
81384	MOUSE ANTI-C-KIT (CLONE 2E4) 2 6 ML
180384	MOUSE ANTI-C-KIT (CLONE 2E4) C 1 ML
187363	MOUSE ANTI-CLAUDIN 2 (CLONE 12 0.5 ML
187364	MOUSE ANTI-CLAUDIN 5 (CLONE 4C 0.5 ML
187391	MOUSE ANTI-CLAUDIN 7 (CLONE 5D 0.5 ML
187341	MOUSE ANTI-CLAUDIN-4 (CLONE 3E 0.5 ML
180176Z	MOUSE ANTI-C-MYC (CLONE 9E10) 1 ML
180128	MOUSE ANTI-COLLAGEN IV (CLONE 1 ML
187381	MOUSE ANTI-COX-1 (CLONE COX111 0.5 ML
187379	MOUSE ANTI-COX-2 (COX-229) CON 0.5 ML
180220	MOUSE ANTI-CYCLIN D1 (CLONE AM 1 ML
081220	MOUSE ANTI-CYCLIN D1 (CLONE AM 6 ML
180404	MOUSE ANTI-CYCLIN E (CLONE HE1 1 ML
08153	MOUSE ANTI-CYTOKERATIN (ACIDIC 1 ML
080153	MOUSE ANTI-CYTOKERATIN (ACIDIC 6 ML
180154	MOUSE ANTI-CYTOKERATIN (BASIC) 1 ML
081212	MOUSE ANTI-CYTOKERATIN (HMW) (6 ML
180132	MOUSE ANTI-CYTOKERATIN (PAN) (1 ML
084132	MOUSE ANTI-CYTOKERATIN (PAN) (10 ML
080132	MOUSE ANTI-CYTOKERATIN (PAN) (6 ML
081132	MOUSE ANTI-CYTOKERATIN (PAN) (6 ML
180158Z	MOUSE ANTI-CYTOKERATIN 18 (CLO 1 ML
080190	MOUSE ANTI-CYTOKERATIN 19 (CLO 6 ML
182200	MOUSE ANTI-CYTOKERATIN 20 (CLO 0.2 ML
081200	MOUSE ANTI-CYTOKERATIN 20 (CLO 6 ML
180267	MOUSE ANTI-CYTOKERATIN 5/6 (CL 1 ML
081267	MOUSE ANTI-CYTOKERATIN 5/6 (CL 6 ML
180234	MOUSE ANTI-CYTOKERATIN 7 (CLON 1 ML
081234	MOUSE ANTI-CYTOKERATIN 7 (CLON 6 ML
180185Z	MOUSE ANTI-CYTOKERATIN 8 (CLON 1 ML
180213	MOUSE ANTI-CYTOKERATIN 8/18 (C 1 ML
081213	MOUSE ANTI-CYTOKERATIN 8/18 (C 6 ML

280001 MOUSE ANTI-CYTOKERATIN PAN (CL 2.5 ML
 182410 MOUSE ANTI-D2-40 CONCENTRATE 0.2 ML
 180016 MOUSE ANTI-DESMIN (CLONE ZC18) 1 ML
 080016 MOUSE ANTI-DESMIN (CLONE ZC18) 6 ML
 081016 MOUSE ANTI-DESMIN (CLONE ZC18) 6 ML
 180223 MOUSE ANTI-E-CADHERIN (CLONE 4 1 ML
 081223 MOUSE ANTI-E-CADHERIN (CLONE 4 6 ML
 81222 MOUSE ANTI-E-CADHERIN (CLONE H 6 ML
 084205 MOUSE ANTI-EGFR (CLONE 31G7) V 10 ML
 180017 MOUSE ANTI-EMA (CLONE ZCE113) 1 ML
 080017 MOUSE ANTI-EMA (CLONE ZCE113) 6 ML
 81017 MOUSE ANTI-EMA (CLONE ZCE113) 6 ML
 182394 MOUSE ANTI-EPHB4 RECEPTOR (CLO 0.2 ML
 182357 MOUSE ANTI-EZRIN (CLONE 3C12) 0.2 ML
 80061 MOUSE ANTI-FACTOR VIII RELATED 6 ML
 81061 MOUSE ANTI-FACTOR VIII RELATED 6 ML
 180475 MOUSE ANTI-GALECTIN-3 (ZG001) CONC 1 ML
 081475 MOUSE ANTI-GALECTIN-3 (ZG001) PRED 6 ML
 187284 MOUSE ANTI-GCDFP-15 (CLONE D6) 0.5 ML
 081021 MOUSE ANTI-GFAP (CLONE ZCG29) 6 ML
 182306 MOUSE ANTI-GMP-1 (SENTRIN-1- S 0.2 ML
 187326 MOUSE ANTI-GPCR-K2 (CLONE 5D5) 0.5 ML
 180005 MOUSE ANTI-HCG- BETA CHAIN SPE 1 ML
 080005 MOUSE ANTI-HCG- BETA CHAIN SPE 6 ML
 187329 MOUSE ANTI-HER4 (CLONE HFR1) C 0.5 ML
 180168 MOUSE ANTI-HLA-DR (CLONE LN-3) 1 ML
 180031 MOUSE ANTI-HUMAN KAPPA LIGHT C 1 ML
 081031 MOUSE ANTI-HUMAN KAPPA LIGHT C 6 ML
 180032 MOUSE ANTI-HUMAN LAMBDA LIGHT 1 ML
 081032 MOUSE ANTI-HUMAN LAMBDA LIGHT 6 ML
 180066 MOUSE ANTI-INSULIN (CLONE Z006 1 ML
 080066 MOUSE ANTI-INSULIN (CLONE Z006 6 ML
 180192Z MOUSE ANTI-KI-67 (CLONE 7B11) 1 ML
 084192 MOUSE ANTI-KI-67 (CLONE 7B11) 10 ML
 081156 MOUSE ANTI-KI-67 2NDGEN PREDIL 6 ML
 080156 MOUSE ANTI-KI-67 PREDILUTE 6 ML
 180473 MOUSE ANTI-LAT 1 ML
 180037 MOUSE ANTI-LUTEINIZING HORMONE 1 ML
 081410 MOUSE ANTI-LYMPHATIC ENDOTHEL 6 ML
 180165 MOUSE ANTI-MACROPHAGES (CLONE 1 ML
 180376 MOUSE ANTI-MAP KINASE (ERK1+ER 1 ML
 187259 MOUSE ANTI-MART-1 (MELAN-A) (C 0.5 ML
 187263 MOUSE ANTI-MART-1 (MELAN-A) CO 0.5 ML
 182403 MOUSE ANTI-MDM2 (CLONE IF2) CO 0.2 ML

07

[Handwritten signature]



Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N. M. A. T.

"2010 – Año del Bicentenario de la Revolución de Mayo"

180485	MOUSE ANTI-MELANOMA PNL2 1 ML
182050	MOUSE ANTI-MELANOSOME (CLONE H 0.2 ML
084050	MOUSE ANTI-MELANOSOME (CLONE H 10 ML
080050	MOUSE ANTI-MELANOSOME (CLONE H 6 ML
081050	MOUSE ANTI-MELANOSOME (CLONE H 6 ML
187337	MOUSE ANTI-MGMT (CLONE MT23.2) 0.5 ML
180369	MOUSE ANTI-MITF (CLONE C5+D5) 1 ML
081369	MOUSE ANTI-MITF (CLONE C5+D5) 6 ML
081342	MOUSE ANTI-MLH1 (CLONE 14) 2ND 6 ML
182342	MOUSE ANTI-MLH1 (CLONE 14) CON 0.2 ML
081445	MOUSE ANTI-MLH1 (ZY001) 2ND GEN PRED 6 ML
180445	MOUSE ANTI-MLH1 (ZY001) CONCENTRATE 1 ML
081343	MOUSE ANTI-MSH2 (CLONE FE11) 2 6 ML
182343	MOUSE ANTI-MSH2 (CLONE FE11) C 0.2 ML
081443	MOUSE ANTI-MSH6 (2D4B5) 2ND GEN 6 ML
180443	MOUSE ANTI-MSH6 (2D4B5) CONCENTRATE 1 ML
182298	MOUSE ANTI-MUC1 (CLONE VU-4-H5 0.2 ML
182299	MOUSE ANTI-MUC2 (CLONE CCP58) 0.2 ML
182322	MOUSE ANTI-MUC4 (CLONE 1G8) CO 0.2 ML
182261	MOUSE ANTI-MUC5AC (CLONE 45M1) 0.2 ML
080160	MOUSE ANTI-MULTI-DRUG RESEISTA 6 ML
187243	MOUSE ANTI-MULTI-DRUG RESISTAN 0.5 ML
187244	MOUSE ANTI-MULTI-DRUG RESISTAN 0.5 ML
180160Z	MOUSE ANTI-MULTI-DRUG RESISTAN 1 ML
180073	MOUSE ANTI-MYOGLOBIN (CLONE Z0 1 ML
180105	MOUSE ANTI-MYOSIN (SKELETAL MU 1 ML
180224	MOUSE ANTI-N-CADHERIN (CLONE 3 1 ML
080139	MOUSE ANTI-NEUROFILAMENT (NF-M 6 ML
180171Z	MOUSE ANTI-NEUROFILAMENT (PAN) 1 ML
187397	MOUSE ANTI-NNOS (CLONE 3G6B10) 0.5 ML
180196	MOUSE ANTI-NSE (CLONE NSE-1G4) 1 ML
080196	MOUSE ANTI-NSE (CLONE NSE-1G4) 6 ML
081196	MOUSE ANTI-NSE (CLONE NSE-1G4) 6 ML
187288	MOUSE ANTI-NUCLEOPHOSMIN (NPM/ 0.5 ML
182359	MOUSE ANTI-NY-ESO-1 (CLONE E97 0.2 ML
180471	MOUSE ANTI-OCT2 (ZO0001) CONC 1 ML
081471	MOUSE ANTI-OCT2 (ZO0001) PREDILUTE 6 ML
81401	MOUSE ANTI-P21 (CLONE EA10) 2N 6 ML
180401	MOUSE ANTI-P21 (CLONE EA10) CO 1 ML
182370	MOUSE ANTI-P27 (CLONE 57) CONC 0.2 ML
180437	MOUSE ANTI-P34 CONCENTRATE 1 ML
081437	MOUSE ANTI-P34 PREDILUTE 6 ML
081129	MOUSE ANTI-P53 (CLONE BP53.12) 6 ML
187251	MOUSE ANTI-P53 (CLONE DO7) CON 0.5 ML

180414	MOUSE ANTI-P63 (4B1E2) CONCENTRATE 1 ML
081414	MOUSE ANTI-P63 (4B1E2) PREDILUTE 6 ML
180406	MOUSE ANTI-P73 (CLONE P73-9F7) 1 ML
180472	MOUSE ANTI-PAX5 (ZP007) CONCENTRATE 1 ML
081472	MOUSE ANTI-PAX5(ZP007) PREDILUTE 6 ML
180407	MOUSE ANTI-PAXILLIN (CLONE ZO3 1 ML
081110	MOUSE ANTI-PCNA (CLONE PC10) 2 6 ML
180110	MOUSE ANTI-PCNA (CLONE PC10) C 1 ML
182400	MOUSE ANTI-PHOSPHO-STAT-1 (TYR 0.2 ML
182262	MOUSE ANTI-PNRA (PROXIMAL NEPH 0.2 ML
180162Z	MOUSE ANTI-PS2 (CLONE BC04) CO 1 ML
180159Z	MOUSE ANTI-RETINOBLASTOMA GENE 1 ML
180307	MOUSE ANTI-SKP2 (CLONE 1G12E9) 1 ML
187248	MOUSE ANTI-TAG-72 (CA 72.4 ANT 0.5 ML
080148	MOUSE ANTI-TAU (CLONE T14) PRE 6 ML
187319	MOUSE ANTI-THYMIDYLATE SYNTHAS 0.5 ML
180399	MOUSE ANTI-THYMIDYLATE SYNTHAS 1 ML
180405	MOUSE ANTI-THYMIDYLATE SYNTHAS 1 ML
081399	MOUSE ANTI-THYMIDYLATE SYNTHAS 6 ML
081405	MOUSE ANTI-THYMIDYLATE SYNTHAS 6 ML
180199	MOUSE ANTI-TRACP (HAIRY CELLS) 1 ML
081199	MOUSE ANTI-TRACP (HAIRY CELLS) 6 ML
180221	MOUSE ANTI-TTF-1 (CLONE 8G7G3/ 1 ML
081221	MOUSE ANTI-TTF-1 (CLONE 8G7G3/ 6 ML
187352	MOUSE ANTI-TYROSINASE (CLONE T 0.5 ML
080147	MOUSE ANTI-UBIQUITIN (CLONE UB 6 ML
180052	MOUSE ANTI-VIMENTIN (CLONE V9) 1 ML
080052	MOUSE ANTI-VIMENTIN (CLONE V9) 6 ML
081052	MOUSE ANTI-VIMENTIN (CLONE V9) 6 ML
182416	MOUSE ANTI-ZAP-70 (CLONE 2F3.2) 0.2 ML
180416	MOUSE ANTI-ZAP-70 (CLONE 2F3.2) 1 ML
081416	MOUSE ANTI-ZAP-70 (CLONE2F3.2) 6 ML

Vida útil: DIECIOCHO (18) MESES, conservado entre 2-8°C Se le asigna la categoría: venta a Laboratorios de análisis clínicos por hallarse en la condiciones establecidas en la Ley N° 16.463 y Resolución Ministerial N° 145/98. Lugar de elaboración: Invitrogen Corporation, 542 Flynn Road, Camarillo, CA 93012 (U.S.A.). En las etiquetas de los envases, anuncios y prospectos deberá constar PRODUCTO PARA DIAGNOSTICO





Ministerio de Salud
Secretaría de Políticas, Regulación
e Institutos
A.N. M. A. T.

"2010 - Año del Bicentenario de la Revolución de Mayo"

USO IN VITRO AUTORIZADO POR LA ADMINISTRACIÓN NACIONAL DE
MEDICAMENTOS, ALIMENTOS Y TECNOLOGIA MEDICA. Certificado n°

006631

ADMINISTRACIÓN NACIONAL DE MEDICAMENTOS, ALIMENTOS Y
TECNOLOGIA MEDICA

Buenos Aires,

03 DIC 2010

Firma y sello

Dr. OTTO A. ORSINGER
SUB-INTERVENTOR
A.N.M.A.T.